

QTRACE and DigitalTRACE Analysis System Operator's Manual


Table of contents

Welcome	3
Product Use Limitations	3
Limited License Agreement	4
Disclaimer	4
Copyright and Trademarks	5
Introduction	6
QTRACE and DigitalTRACE System Workflow	9
Materials	11
Key to Symbols	11
Materials Provided	12
Materials Sold Separately	16
DNA Sample Requirements	18
Warnings	19
Shipping Conditions	20
Storage and Handling Requirements	20
Technical Support	21
Genotyping Test	23
Protocol	25
Genotyping Data Analysis and Report	37
Multiple Donor Analysis Using QTRACE Plates	42
Custom Genotyping Panels	46
Virtual Typing	50
Monitoring Test	51
Assigning Informative Markers	53
Protocol	55
Monitoring Data Analysis and Report	61
Assay Filtering from Monitoring Reports	64
External Reference Functions	66
Simultaneous Genotyping and Monitoring	68
Software Overview	70
Preferences	70
Software Access And User Profile Management	78
File Management	80
Data Exports - csv and xlsx	82
HPRIM Data Export	84
Software Buttons	87
Drag and Drop Plate Layouts	90
Anonymized Reporting	91
Modification of Existing Data	92
Data Analysis Algorithms	94
Bibliography	100
Glossary	101
End User Software License Agreement	103

Welcome

Welcome to the QTRACE® and DigitalTRACE™ Analysis System Operator's Manual and Help System. This file serves as both the QTRACE® System Operator's Manual (PN 3 31047, JETA Molecular), the DigitalTRACE™ Operator's Manual (PN 331307) and the help system found within the TRACE Analysis™ Software package.

Navigate to your subject of interest and find the solution to your question. If your topic of interest is not listed or you would like to receive additional information, do not hesitate to contact us. The content of this help system is regularly updated. We encourage you to inform us on inaccuracies or suggestions. We do our utmost to implement your suggestions swiftly, such that you and other QTRACE® and DigitalTRACE™ System users may benefit from it.



Product Use Limitations

This version of the TRACE Analysis™ Software is for Research Use Only. It is not intended for use in diagnostic procedures.

No claim or representation is intended to provide information for the diagnosis, prevention, or treatment of a disease.

Limited License Agreement

Use of this product signifies the agreement of any purchaser or user of the QTRACE® and DigitalTRACE™ Analysis System kits or components with the following terms:

The QTRACE® and DigitalTRACE™ Kits may be used solely in accordance with the QTRACE® and DigitalTRACE™ Analysis System Operator's Manual and for use with components contained in the kits only. JETA Molecular grants no license under any of its intellectual property to use or incorporate the enclosed components of these kits with any components not included within these kits except as described in the QTRACE® and DigitalTRACE™ Analysis System Operator's Manual and additional protocols available at www.jetamolecular.com.

Other than expressly stated licenses, JETA Molecular makes no warranty that these kits and/or their use(s) do not infringe the rights of third-parties.

The QTRACE® and DigitalTRACE™ Kits and their components are licensed for one-time use and may not be re-used, re-furbished, re-sold or reverse engineered.

JETA Molecular specifically disclaims any other licenses, expressed or implied other than those expressly stated.

The purchaser and user of the QTRACE® and DigitalTRACE™ Kits agree not to take or permit anyone else to take any steps that could lead to or facilitate any acts prohibited above. JETA Molecular may enforce the prohibitions of this Limited License Agreement in any Court, and shall recover all its investigative and Court costs, including attorney fees, in any action to enforce this Limited License Agreement or any of its intellectual property rights relating to the kit and/or its components.

Disclaimer

JETA Molecular makes every effort to ensure that this manual is accurate. JETA Molecular disclaims liability for any inaccuracies or omissions that may have occurred. Information in this manual is subject to change without notice.

JETA Molecular assumes no responsibility for any inaccuracies that may be contained in this manual.

JETA Molecular reserves the right to make improvements to this manual and/or to the products described in this manual, at any time without notice.

If you find information in this manual that is incorrect, misleading, or incomplete, we would appreciate your comments and suggestions. Please send them to info@jetabv.com.

Copyright and Trademarks

This publication, including all photographs, illustrations, is protected under international copyright laws, with all rights reserved. Neither this manual, nor any of the material contained herein, may be reproduced without written consent of JETA Molecular.

JETA Molecular and its logo design are registered trademarks of JETA Molecular BV in the U.S. and/or certain other countries.

QTRACE® and DigitalTRACE™ are registered trademarks of ElsworthMolecular Holding BV in the U.S. and/or certain other countries. QTRACE® Products and Software are licensed exclusively to JETA Molecular BV from ElsworthMolecular Holding BV.

MicroAmp® is a registered trademark of Thermo Fisher.

Moq

Copyright © 2007. Clarius Consulting, Manas Technology Solutions, InSTEDD
<http://www.moqthis.com/>
All rights reserved.

NUnit

Portions Copyright © 2002-2013 Charlie Poole or Copyright © 2002-2004 James W. Newkirk, Michael C. Two, Alexei A. Vorontsov or Copyright © 2000-2002 Philip A. Craig

PDFsharp

Copyright © 2005-2007 empira Software GmbH, Cologne (Germany)

SharpZipLib

WPF Toolkit

All other trademarks, copyrights, patents, service marks, logos and trade names are the sole property of their respective owners.

Introduction

Background Information

The QTRACE®, DigitalTRACE™ and MultiTRACE™ Genotyping Plates, QTRACE® and DigitalTRACE™ Monitoring Assays and the TRACE Analysis™ Software meet the needs of any research application that requires highly sensitive detection and quantification of the genome of one individual in the background of another individual or individuals. A genetic chimera is an organism with two or more genetically distinct cell populations, *i.e.*, cell populations with different genomes.

Chimerism can arise in humans through a variety of means, such as inheritance, maternal-fetal stem cell trafficking during gestation, blood vessel sharing in fraternal twin gestation, blood transfusions, bone marrow transplantation, cord blood transplantation, and solid organ transplantation. The presence of two distinct human genomes in a sample can also occur simply through the mixing of human cells from more than one individual, for example, when two cell lines are cross-contaminated, or in forensic tissue samples.

Product Overview

The QTRACE® and the DigitalTRACE™ Analysis System consist of DNA genotyping plates, individual quantification assays, one reference assay, and the TRACE Analysis™ Software.

The QTRACE® INDEL Assays are a set of 80 genetic markers that are able to differentiate, and then quantify, the contributors to a human-mixed DNA sample using quantitative polymerase chain reaction (qPCR). Each of the 80 assays is designed to a distinct bi-allelic insertion/deletion (INDEL) or copy number polymorphisms in the human genome. The QTRACE® INDEL assays are based upon self-quenched, hydrolysis probe chemistry. In a qPCR, a dye labeled oligonucleotide probe enables the detection of a specific PCR product as it accumulates during PCR cycling. qPCR has been shown in various studies to quantify target DNA samples over an 11-log dynamic range in optimized conditions (100 billion-fold differences in starting copy number; Nolan, *et al*, 2006). Comparison of cycle thresholds (CTs) from different samples can be used to determine the relative amounts of DNA in two different samples (Livak and Schmittgen, 2001). The high sensitivity is the result of the very large dynamic range of the real-time amplification method and is limited essentially by the input copy number of total genomic DNA that can be added to the PCR reaction. The RNase P assay serves as reference assay for the quantification.

The DigitalTRACE™ INDEL assays are a set of 70 digital polymerase chain reaction (dPCR) assays based upon self-quenched, hydrolysis probe chemistry. Each of the assays is designed to a distinct bi-allelic insertion/deletion (INDEL) or copy number polymorphisms in the human genome. In a dPCR reaction, a dye-labeled oligonucleotide probe enables the detection of a specific PCR product at the end PCR cycling. The high sensitivity is the result of the very large dynamic range of the real-time amplification method and is limited essentially by the input copy number of total DNA that can be added to the dPCR reaction. In the Monitoring test, the assays are formulated with the target in channel FAM and the reference assay (RNase P) in channel HEX.

The QTRACE® and DigitalTRACE™ HLA Assays are qPCR or dPCR assays based upon self-quenched, hydrolysis probe chemistry. These assays are designed to detect the loss of HLA heterozygosity after haploidentical HSCT.

The TRACE Analysis™ Software was designed specifically for the QTRACE® and DigitalTRACE™ INDEL Assay Set. The software provides a streamlined workflow for both the genotyping and quantification tests. The software guides the user through assay setup, performs data analysis, generates results reports and stores the data collected for samples over time.

The procedure for determining the level of a genome of interest in a sample consists of two parts: a genotyping test and a quantification test.

Genotyping Test

In the initial genotyping test, the DNAs that comprise a mixed DNA sample are analyzed using a QTRACE®, DigitalTRACE™ or MultiTRACE™ Genotyping Plate, to identify all of the informative assays for the samples. An informative assay is an assay for a marker allele that is present (positive) in one individual genome *and* absent (negative) in the other genome.

The QTRACE® Genotyping plate contains a duplicate set of assays: the 46 quantification assays and the reference (RNaseP) assay that serves as both a positive control and a No Template Control (NTC).

The QTRACE® Extended Panel Genotyping Plate, containing additional 34 markers, is available for rare cases where more markers may be needed, including patients receiving a second transplant or highly similar siblings.

The DigitalTRACE™ Genotyping Plate contains a set of 43 quantification assays and the reference (RNase P) assay that serves as both a positive control and a No Template Control (NTC).

The DigitalTRACE™ EP QIAcuity Genotyping Plate represents an extended panel of dPCR markers that can be used for dPCR monitoring in the case of need for additional markers allowing to distinguish between donor and recipient DNA. This plate contains a set of 27 quantification assays and the reference (RNase P) assay that serves as both a positive control and a No Template Control (NTC).

The MultiTRACE™ Genotyping Plate contains a set of 45 quantification assays and the reference (RNase P) assay that serves as both a positive control and a No Template Control (NTC).

Monitoring Test

In the monitoring (quantification) test, two or more of the informative assays identified in the genotyping test is used to quantify the DNA of interest in an unknown sample relative to a reference sample (calibrator). Any of the informative assays identified in the genotyping test can be used to perform a quantification test. The amount of the genome positive for the informative allele in the unknown sample is determined relative to the amount of that same genome in the reference sample, and the result is expressed as a percentage (ratio). For example, a result of 5% indicates that there is 5% of genome A in the unknown sample relative to the reference sample. The simplest case assumes that the reference sample has 100% of genome A.

The informativeness of a multi-locus genotyping panel is a measure of the probability of finding at least one informative assay between two individual genomes (or DNA samples). Informativeness is calculated from the population frequency estimates of the alleles used to

make up a multi-locus genotyping panel, and thus differs between ethnic populations. In addition, the informativeness of any panel of polymorphic loci is higher in unrelated individuals than in related individuals.

Note: Estimates are based upon both public and proprietary allele frequency data for the 80 assay panel.

The performance of the QTRACE® and DigitalTRACE™ Analysis System has been verified to a level of 0.1% minor component DNA in 150ngs total DNA.

QTRACE and DigitalTRACE System Workflow

The QTRACE® Analysis System



Enter Sample Information into QTRACE™ Software to Setup the Genotyping Experiment. Export a SetUp File for the qPCR machine and generate a lab protocol

Add 5X qPCR Master Mix + Sample DNA to the Genotyping plate per the Sample Layout generated by the software

Load genotyping plate into the qPCR machine. Open a QTRACE™ template, import the setup file and start the run.

Export data from qPCR Machine and import into QTRACE™ Software for informative marker identification

Enter Sample Information into QTRACE™ Software to Setup the Monitoring Experiment. Export a SetUp File for the qPCR machine and generate a lab protocol



Add 5X qPCR Master Mix + Sample DNA to the Plate per the Assay and Sample Layouts generated by the Software

Load Monitoring plate into the qPCR machine. Open a QTRACE™ template, import the setup file and start the run

Export data from qPCR Machine and import into QTRACE™ Software for analysis and reporting of chimeric mixture











The DigitalTRACE™ Analysis System














Materials

Key to Symbols









The following symbols appear within the labeling of the QTRACE® and DigitalTRACE™ System Products:

Symbol	Definition	Symbol	Definition
	Part Number		Store Below -15°C
	Lot Number		Store between 15°C – 30°C
	Expiration Date		Manufacturer
	See Operator's Manual		Contains Sufficient for <n> Tests

Materials Provided

REF	Name	Description	Storage Conditions	Unit
121045	QTRACE® Genotyping Plates, 4-Pack	Four ABI MicroAmp® Optical 96-well plates pre-arrayed with INDEL qPCR Assays; individually sealed		Box
121056	QTRACE® Genotyping Plates - FAST	Four ABI FAST MicroAmp® Optical 96-well plates pre-arrayed with INDEL qPCR Assays; individually sealed		Box
121066	QTRACE® Genotyping Plates – Roche 480	Four Roche 480 Optical 96-well plates pre-arrayed with INDEL qPCR Assays; individually sealed		Box
121129	MultiTRACE™ Genotyping Plate Pack – ABI 0.2mL v2	Four ABI MicroAmp® Optical 96-well plates pre-arrayed with INDEL qPCR Assays in multiplex; individually sealed		Box
121220	MultiTRACE™ Genotyping Plate Pack - ABI 0.1mL v2	Four ABI MicroAmp® Optical 96-well plates pre-arrayed with INDEL qPCR Assays in multiplex; individually sealed		Box
121221	MultiTRACE™ Genotyping Plate Pack - LC480 v2	Four Roche 480 Optical 96-well plates pre-arrayed with INDEL qPCR Assays in multiplex; individually sealed		Box
121226	MultiTRACE™ Genotyping Plate Pack – ABI 0.2mL v3	Four ABI MicroAmp® Optical 96-well plates pre-arrayed with INDEL qPCR Assays in multiplex; individually sealed		Box
121285	MultiTRACE™ Genotyping Plate Pack – ABI 0.1mL v3	Four ABI MicroAmp® Optical 96-well plates pre-arrayed with INDEL qPCR Assays in multiplex; individually sealed		Box
121228	MultiTRACE™ Genotyping Plate Pack – LC480 v3	Four Roche 480 Optical 96-well plates pre-arrayed with INDEL qPCR Assays in multiplex; individually sealed		Box
121139	QTRACE® Extended Panel Genotyping Plate - ABI 0.2mL	Four ABI MicroAmp® Optical 96-well plates pre-arrayed with INDEL qPCR Assays; individually sealed		Box
121045	DigitalTRACE™ QIAcuity Genotyping Plate	Two ABI MicroAmp® Optical 96-well plates pre-arrayed with INDEL dPCR Assays; individually sealed		Box

121056	DigitalTRACE™ EP QIAcuity Genotyping Plate	Two ABI MicroAmp® Optical 96-well plates pre-arrayed with INDEL dPCR Assays; individually sealed		Box
--------	--	--	---	-----

REF	Name	Description	Storage Conditions	Unit
311044	QTRACE® qPCR Master Mix	1.6ml Buffered solution of dNTPs, a DNA Polymerase, a passive reference dye and MgCl ₂		Tube
311043	QTRACE® RNaseP Assay	1.1ml buffered solution containing a mix of primers and probe for detecting the RNaseP gene		Tube
	QTRACE® INDEL Assays	375µl buffered solution containing a mix of primers and probe for detecting the variant of interest		Tube
	DigitalTRACE™ INDEL Assays	26µl buffered solution containing a mix of primers and probe for detecting the variant of interest in FAM and a reference assay in HEX		Tube
	QTRACE® HLA Assays	100µl buffered solution containing a mix of primers and probe for detecting the variant of interest		Tube
	DigitalTRACE™ HLA Assays	26µl buffered solution containing a mix of primers and probe for detecting the variant of interest in FAM and a reference assay in HEX		Tube
711295	QTRACE® Universal Positive Control	360µl buffered solution containing synthetic DNA serving as positive control		Tube
711294	DigitalTRACE™ Universal Positive Control	360µl buffered solution containing synthetic DNA serving as positive control		Tube

QTRACE® Reference Numbers:

REF	Name	REF	Name
211140	QTRACE® INDEL Assay 102	211030	QTRACE® INDEL Assay 748
211141	QTRACE® INDEL Assay 113	211053	QTRACE® INDEL Assay 755
211142	QTRACE® INDEL Assay 120	211031	QTRACE® INDEL Assay 768
211310	QTRACE® INDEL Assay 126	211167	QTRACE® INDEL Assay 777
211001	QTRACE® INDEL Assay 137	211032	QTRACE® INDEL Assay 784
211002	QTRACE® INDEL Assay 148	211168	QTRACE® INDEL Assay 795

211143	QTRACE® INDEL Assay 157	211033	QTRACE® INDEL Assay 803
211311	QTRACE® INDEL Assay 161	211052	QTRACE® INDEL Assay 819
211144	QTRACE® INDEL Assay 176	211034	QTRACE® INDEL Assay 824
211145	QTRACE® INDEL Assay 183	211035	QTRACE® INDEL Assay 832
211146	QTRACE® INDEL Assay 198	211036	QTRACE® INDEL Assay 840
211003	QTRACE® INDEL Assay 209	211037	QTRACE® INDEL Assay 854
211312	QTRACE® INDEL Assay 216	211317	QTRACE® INDEL Assay 861
211147	QTRACE® INDEL Assay 222	211169	QTRACE® INDEL Assay 874
211004	QTRACE® INDEL Assay 235	211170	QTRACE® INDEL Assay 884
211005	QTRACE® INDEL Assay 240	211038	QTRACE® INDEL Assay 888
211148	QTRACE® INDEL Assay 252	211171	QTRACE® INDEL Assay 892
211006	QTRACE® INDEL Assay 267	211039	QTRACE® INDEL Assay 907
211149	QTRACE® INDEL Assay 275	211040	QTRACE® INDEL Assay 916
211313	QTRACE® INDEL Assay 284	211172	QTRACE® INDEL Assay 923
211007	QTRACE® INDEL Assay 291	211173	QTRACE® INDEL Assay 936
211008	QTRACE® INDEL Assay 305	211041	QTRACE® INDEL Assay 948
211009	QTRACE® INDEL Assay 312	211042	QTRACE® INDEL Assay 954
211010	QTRACE® INDEL Assay 326	211174	QTRACE® INDEL Assay 962
211150	QTRACE® INDEL Assay 333	211175	QTRACE® INDEL Assay 971
211011	QTRACE® INDEL Assay 345	211176	QTRACE® INDEL Assay 987
211012	QTRACE® INDEL Assay 356	211177	QTRACE® INDEL Assay 990
211013	QTRACE® INDEL Assay 359	211078	QTRACE® HLA Assay H005
211014	QTRACE® INDEL Assay 361	211080	QTRACE® HLA Assay H007
211055	QTRACE® INDEL Assay 373	211081	QTRACE® HLA Assay H009
211064	QTRACE® INDEL Assay 386	211083	QTRACE® HLA Assay H017
211151	QTRACE® INDEL Assay 396	211130	QTRACE® HLA Assay H020
211015	QTRACE® INDEL Assay 408	211085	QTRACE® HLA Assay H022
211314	QTRACE® INDEL Assay 411	211087	QTRACE® HLA Assay H025
211016	QTRACE® INDEL Assay 425	211088	QTRACE® HLA Assay H028
211017	QTRACE® INDEL Assay 434	211131	QTRACE® HLA Assay H029
211152	QTRACE® INDEL Assay 441	211091	QTRACE® HLA Assay H036
211153	QTRACE® INDEL Assay 450	211092	QTRACE® HLA Assay H038
211315	QTRACE® INDEL Assay 457	211093	QTRACE® HLA Assay H039
211018	QTRACE® INDEL Assay 469	211094	QTRACE® HLA Assay H041
211154	QTRACE® INDEL Assay 472	211095	QTRACE® HLA Assay H043
211155	QTRACE® INDEL Assay 482	211096	QTRACE® HLA Assay H045
211156	QTRACE® INDEL Assay 493	211098	QTRACE® HLA Assay H051
211019	QTRACE® INDEL Assay 504	211099	QTRACE® HLA Assay H052
211054	QTRACE® INDEL Assay 519	211100	QTRACE® HLA Assay H053
211020	QTRACE® INDEL Assay 520	211101	QTRACE® HLA Assay H054
211021	QTRACE® INDEL Assay 531	211132	QTRACE® HLA Assay H101
211022	QTRACE® INDEL Assay 548	211133	QTRACE® HLA Assay H102
211157	QTRACE® INDEL Assay 555	211134	QTRACE® HLA Assay H103
211158	QTRACE® INDEL Assay 567	211279	QTRACE® HLA Assay H104
211159	QTRACE® INDEL Assay 574	211280	QTRACE® HLA Assay H105
211160	QTRACE® INDEL Assay 585	211281	QTRACE® HLA Assay H106
211161	QTRACE® INDEL Assay 597	211282	QTRACE® HLA Assay H107
211023	QTRACE® INDEL Assay 601	211283	QTRACE® HLA Assay H108
211024	QTRACE® INDEL Assay 615	211284	QTRACE® HLA Assay H109
211025	QTRACE® INDEL Assay 626	211290	QTRACE® HLA Assay H110
211026	QTRACE® INDEL Assay 634	211291	QTRACE® HLA Assay H111

211316	QTRACE® INDEL Assay 647	211292	QTRACE® HLA Assay H112
211027	QTRACE® INDEL Assay 650	211293	QTRACE® HLA Assay H113
211162	QTRACE® INDEL Assay 663	211296	QTRACE® HLA Assay H114
211164	QTRACE® INDEL Assay 678	211297	QTRACE® HLA Assay H115
211165	QTRACE® INDEL Assay 681	211298	QTRACE® HLA Assay H116
211166	QTRACE® INDEL Assay 694	211299	QTRACE® HLA Assay H117
211028	QTRACE® INDEL Assay 706	211300	QTRACE® HLA Assay H118
211065	QTRACE® INDEL Assay 710	211301	QTRACE® HLA Assay H119
211051	QTRACE® INDEL Assay 721	211319	QTRACE® HLA Assay H120
211029	QTRACE® INDEL Assay 736		

DigitalTRACE™ Reference Numbers:

REF	Name	REF	Name
811140	DigitalTRACE™ INDEL Assay 102	811051	DigitalTRACE™ INDEL Assay 721
811141	DigitalTRACE™ INDEL Assay 113	811029	DigitalTRACE™ INDEL Assay 736
811142	DigitalTRACE™ INDEL Assay 120	811030	DigitalTRACE™ INDEL Assay 748
811310	DigitalTRACE™ INDEL Assay 126	811053	DigitalTRACE™ INDEL Assay 755
811001	DigitalTRACE™ INDEL Assay 137	811031	DigitalTRACE™ INDEL Assay 768
811002	DigitalTRACE™ INDEL Assay 148	811167	DigitalTRACE™ INDEL Assay 777
811143	DigitalTRACE™ INDEL Assay 157	811032	DigitalTRACE™ INDEL Assay 784
811311	DigitalTRACE™ INDEL Assay 161	811168	DigitalTRACE™ INDEL Assay 795
811144	DigitalTRACE™ INDEL Assay 176	811033	DigitalTRACE™ INDEL Assay 803
811145	DigitalTRACE™ INDEL Assay 183	811052	DigitalTRACE™ INDEL Assay 819
811146	DigitalTRACE™ INDEL Assay 198	811034	DigitalTRACE™ INDEL Assay 824
811003	DigitalTRACE™ INDEL Assay 209	811035	DigitalTRACE™ INDEL Assay 832
811312	DigitalTRACE™ INDEL Assay 216	811036	DigitalTRACE™ INDEL Assay 840
811147	DigitalTRACE™ INDEL Assay 222	811037	DigitalTRACE™ INDEL Assay 854
811004	DigitalTRACE™ INDEL Assay 235	811317	DigitalTRACE™ INDEL Assay 861
811005	DigitalTRACE™ INDEL Assay 240	811169	DigitalTRACE™ INDEL Assay 874
811148	DigitalTRACE™ INDEL Assay 252	811170	DigitalTRACE™ INDEL Assay 884
811006	DigitalTRACE™ INDEL Assay 267	811038	DigitalTRACE™ INDEL Assay 888
811149	DigitalTRACE™ INDEL Assay 275	811171	DigitalTRACE™ INDEL Assay 892
811313	DigitalTRACE™ INDEL Assay 284	811039	DigitalTRACE™ INDEL Assay 907
811007	DigitalTRACE™ INDEL Assay 291	811040	DigitalTRACE™ INDEL Assay 916
811008	DigitalTRACE™ INDEL Assay 305	811172	DigitalTRACE™ INDEL Assay 923
811009	DigitalTRACE™ INDEL Assay 312	811173	DigitalTRACE™ INDEL Assay 936
811010	DigitalTRACE™ INDEL Assay 326	811041	DigitalTRACE™ INDEL Assay 948
811150	DigitalTRACE™ INDEL Assay 333	811042	DigitalTRACE™ INDEL Assay 954
811011	DigitalTRACE™ INDEL Assay 345	811174	DigitalTRACE™ INDEL Assay 962
811012	DigitalTRACE™ INDEL Assay 356	811175	DigitalTRACE™ INDEL Assay 971
811013	DigitalTRACE™ INDEL Assay 359	811176	DigitalTRACE™ INDEL Assay 987
811014	DigitalTRACE™ INDEL Assay 361	811177	DigitalTRACE™ INDEL Assay 990
811055	DigitalTRACE™ INDEL Assay 373	811078	DigitalTRACE™ HLA Assay H005
811064	DigitalTRACE™ INDEL Assay 386	811080	DigitalTRACE™ HLA Assay H007
811151	DigitalTRACE™ INDEL Assay 396	811083	DigitalTRACE™ HLA Assay H017
811015	DigitalTRACE™ INDEL Assay 408	811084	DigitalTRACE™ HLA Assay H020
811314	DigitalTRACE™ INDEL Assay 411	811085	DigitalTRACE™ HLA Assay H022
811016	DigitalTRACE™ INDEL Assay 425	811087	DigitalTRACE™ HLA Assay H025
811017	DigitalTRACE™ INDEL Assay 434	811088	DigitalTRACE™ HLA Assay H028
811152	DigitalTRACE™ INDEL Assay 441	811131	DigitalTRACE™ HLA Assay H029
811153	DigitalTRACE™ INDEL Assay 450	811091	DigitalTRACE™ HLA Assay H036
811315	DigitalTRACE™ INDEL Assay 457	811092	DigitalTRACE™ HLA Assay H038

811018	DigitalTRACE™ INDEL Assay 469	811093	DigitalTRACE™ HLA Assay H039
811154	DigitalTRACE™ INDEL Assay 472	811094	DigitalTRACE™ HLA Assay H041
811155	DigitalTRACE™ INDEL Assay 482	811095	DigitalTRACE™ HLA Assay H043
811156	DigitalTRACE™ INDEL Assay 493	811096	DigitalTRACE™ HLA Assay H045
811019	DigitalTRACE™ INDEL Assay 504	811098	DigitalTRACE™ HLA Assay H051
811054	DigitalTRACE™ INDEL Assay 519	811099	DigitalTRACE™ HLA Assay H052
811020	DigitalTRACE™ INDEL Assay 520	811100	DigitalTRACE™ HLA Assay H053
811021	DigitalTRACE™ INDEL Assay 531	811101	DigitalTRACE™ HLA Assay H054
811022	DigitalTRACE™ INDEL Assay 548	811133	DigitalTRACE™ HLA Assay H102
811157	DigitalTRACE™ INDEL Assay 555	811134	DigitalTRACE™ HLA Assay H103
811158	DigitalTRACE™ INDEL Assay 567	811279	DigitalTRACE™ HLA Assay H104
811159	DigitalTRACE™ INDEL Assay 574	811280	DigitalTRACE™ HLA Assay H105
811160	DigitalTRACE™ INDEL Assay 585	811281	DigitalTRACE™ HLA Assay H106
811161	DigitalTRACE™ INDEL Assay 597	811282	DigitalTRACE™ HLA Assay H107
811023	DigitalTRACE™ INDEL Assay 601	811283	DigitalTRACE™ HLA Assay H108
811024	DigitalTRACE™ INDEL Assay 615	811284	DigitalTRACE™ HLA Assay H109
811025	DigitalTRACE™ INDEL Assay 626	811290	DigitalTRACE™ HLA Assay H110
811026	DigitalTRACE™ INDEL Assay 634	811291	DigitalTRACE™ HLA Assay H111
811316	DigitalTRACE™ INDEL Assay 647	811292	DigitalTRACE™ HLA Assay H112
811027	DigitalTRACE™ INDEL Assay 650	811293	DigitalTRACE™ HLA Assay H113
811162	DigitalTRACE™ INDEL Assay 663	811296	DigitalTRACE™ HLA Assay H114
811163	DigitalTRACE™ INDEL Assay 670	811297	DigitalTRACE™ HLA Assay H115
811164	DigitalTRACE™ INDEL Assay 678	811298	DigitalTRACE™ HLA Assay H116
811165	DigitalTRACE™ INDEL Assay 681	811299	DigitalTRACE™ HLA Assay H117
811166	DigitalTRACE™ INDEL Assay 694	811300	DigitalTRACE™ HLA Assay H118
811028	DigitalTRACE™ INDEL Assay 706	811319	DigitalTRACE™ HLA Assay H120
811065	DigitalTRACE™ INDEL Assay 710		

REF	Name	Description
341048	TRACE Analysis™ Software	<p>Minimum System Requirements: Windows 7, 2 GB RAM, 250 MB free disk space, Network connection allowing TCP/IP traffic to and from port 3500, Microsoft .NET framework 4.5.</p> <p>Recommended System Requirements: Windows 7, 2 GB RAM, Internet connection for license validation and automatic updates, 500 MB free disk space, Microsoft .NET framework 4.5</p>

REF	Name	Description
331047	QTRACE® System Operator's Manual	Operator's Manual for QTRACE® Analysis System
331307	DigitalTRACE™ Operator's Manual	Operator's Manual for DigitalTRACE™ Analysis System

Materials Sold Separately

Additional Reagents Required but not Provided

Item name	Catalog number
Modified TE Buffer (10 mM Tris-HCl, pH 8.0, 0.1 mM EDTA) (also called TE 0.1 Buffer or TE-4 Buffer)	Not applicable

Molecular grade water (DNase and RNase free)	Not applicable
QIAcuity Probe PCR Kit (Qiagen)	250102
ddPCR Supermix for Probes (Biorad)	1863023

Additional Consumables Required but not Provided

Item name	Catalog number
Pipette Tips, disposable, sterile, aerosol-resistant, filtered, capable of dispensing up to 20, 200, and 1000 µL	Not applicable
1.5 mL microcentrifuge tubes	Not applicable
Lint-free tissue	Not applicable
Gloves, powder-free	Not applicable
96-well plates or strip tubes	Not applicable
Plate seals	Not applicable
QIAcuity Nanoplate 26k 24-well (Qiagen)	250001
QIAcuity Nanoplate 26k 8-well (Qiagen)	250031
QIAcuity Nanoplate 8.5k 96-well (Qiagen)	250021
QIAcuity Nanoplate 8.5k 24-well (Qiagen)	250011
DG8 Cartridges (Biorad)	1864008
DG8 Gaskets (Biorad)	1863009
ddPCR 96-well PCR Plates (Biorad)	12001925
PCR Plate Heat Seal, foil, pierceable (Biorad)	1814040
Droplet Generation Oil for Probes (Biorad)	1863005
ddPCR Droplet Reader Oil (Biorad)	1863004

Additional Equipment Required but not Provided

Item name	Catalog number
Adjustable single channel pipettes (0.5-1000 µL capacity)	Not applicable
Adjustable multi-channel, multi-dispensing pipettes (0.5-200 µL capacity)	Not applicable
Vortex mixer with flat rubber platform head	Not applicable
Centrifuge	Not applicable
Centrifuge with microtiter plate assembly	Not applicable
PC for the installation of TRACE Analysis™ Software	Not applicable

DNA Sample Requirements

Purified DNA should have an A_{260}/A_{280} ratio between 1.7 and 2.0.

We recommend using a fluorometric method to accurately quantify DNA.

If necessary, DNA should be diluted in 10mM Tris, pH 8.0; 0.1mM EDTA (TE) or nuclease-free H₂O before use.

Blood samples should be collected in ACD or EDTA anticoagulation tubes. Do NOT use heparinized samples, as this has an inhibitory effect on a PCR

The optimal amount of template DNA to use in genotyping is 5 ng per well for the QTRACE® Genotyping plates and 10 ng per well for the DigitalTRACE™ Genotyping plates. To streamline the process, validate your DNA purification procedure so you can use a set volume corresponding to 1-10 ng DNA.

The optimal amount of template DNA to use in chimeric mixture analysis depends on the desired sensitivity detection limit.

The following chart shows the relationship between input DNA amounts and sensitivity. These calculations assume at least 10 copies of the minor component DNA in the total amount of DNA.

Input DNA (ng)	# Cells	Sensitivity %
150	22 727	0.04%
100	15 151	0.07%
50	7 576	0.13%
25	3 788	0.26%
12.5	1 894	0.53%
6.25	947	1.06%

Sensitivity %	# Cells	Input DNA (ng)
0.05%	20 000	132
0.10%	10 000	66
1.00%	1 000	7
2.00%	500	3
5.00%	200	1

Warnings

For Research Use Only. Not for use in diagnostic procedures.

Use good laboratory practices for sample handling and tracking.

Use only recommended materials, procedures, and equipment.

Use sterile disposable pipettes and filtered pipette tips.

Wear appropriate personal protective equipment (*e.g.*, safety glasses, disposable gloves, and protective clothing) when handling samples and reagents.

Clean and disinfect all work surfaces with a 10% bleach (0.525% sodium hypochlorite) solution and follow with 70% ethanol, ensuring that all bleach residue is removed.

Assays should be run by individuals experienced in good laboratory practices and who have been previously trained to use the equipment by the original equipment manufacturer (OEM).

Operate, calibrate and maintain all instruments and equipment according to procedures provided by the manufacturers.

To reduce the risk of contamination, the area where amplified DNA is handled must be physically isolated from the work areas for sample preparation and qPCR setup.

Do not use components past their expiration date.

Do not dilute reagents.

Visually inspect wells or tubes after pipetting steps to detect operator errors with pipetting, sample transfer, etc.

To prevent repeated freeze/thaw cycling of reagents during frozen storage, do not store reagents within freezers that use an automatic defrost function (*i.e.*, frost-free).

Avoid microbial and nuclease contamination of reagents when removing aliquots from reagent tubes.

To prevent contamination, after aliquots are removed do not return the remaining volume to the original tube.

Comply with all local, state, or national laws and regulations related to chemical storage and disposal.

CHEMICAL HAZARD. Ethanol is a flammable liquid and vapor. Exposure can cause eye, skin, and upper respiratory tract irritation. Prolonged or repeated contact can dry the skin. Exposure can cause central nervous system depression and liver damage. Keep away from heat, sparks, and flame. Read the MSDS and follow the handling instructions. Wear appropriate protective eyewear, clothing and gloves.

CHEMICAL HAZARD. Bleach (sodium hypochlorite) is a corrosive liquid and vapor. Exposure can cause severe irritation or damage to eyes, skin and the respiratory system. Harmful if swallowed. Prolonged or repeated contact can lead to sensitization (*e.g.*, irritation) if skin damage occurs during exposure. Medical conditions that can be aggravated by exposure to high concentrations of vapor or mist include heart conditions or chronic respiratory problems such as asthma, emphysema, chronic bronchitis or obstructive lung disease. Read the MSDS and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

TRACE Analysis™ Software does not provide a mechanism to edit data files obtained from the qPCR and dPCR System or result files once they are created.

Shipping Conditions

The QTRACE®, DigitalTRACE™ and MultiTRACE™ Genotyping Plates are shipped at ambient temperature. The QTRACE® PCR Master Mix, QTRACE® INDEL Assays, QTRACE® HLA Assays, DigitalTRACE™ INDEL Assays and DigitalTRACE™ HLA Assays are shipped frozen.

Examine the shipment upon receipt and if the integrity of the products has been compromised during shipment, immediately contact your local customer support representative.

Storage and Handling Requirements

Upon receipt, the individual components should be stored according to the temperatures listed on the labels.

Note: When storage recommendations are observed, both unopened and opened/recapped tubes are stable until the expiration date indicated on the label. genotyping plates are stable when stored in the air-tight pouch containing a desiccant bag. Do not use any component after the expiration date.

Do not use any component that visibly shows signs of having been compromised (*e.g.*, particulate matter, presence of foreign debris, cloudy appearance, discoloration).

Technical Support

For technical assistance and more information:

Please contact your local distributor

France

bionobis

7 rue Nicolas Copernic
78280 Guyancourt - SQY
France

Tél : 01 30 43 01 55

Fax : 01 30 43 01 15

Email: info@bionobis.com

<https://www.bionobis.com>

Iberia

Werfen

Plaza Europa, 21-23
08908 L'Hospitalet de Llobregat
Barcelona, Spain

Tel: +34 93 401 01 01

Email: customerservice-es@werfen.com

<https://www.werfen.com/es/es>

Estonia

Quantum Eesti AS

Tiigrisilma 8
Räni alevik
Kambja vald
61708 Tartu, Estonia

Tel: +372 7 301 321

Email: quantum@quantum.ee

<https://quantum.ee/>

Turkiye

ATC Genomics

Galip Erdem Cad. 607. Sok. No:9
06550 Çankaya - Ankara, Türkiye

Tel: +90 312 496 43 14

Email: bilgi@atc.com.tr

<http://www.atc.com.tr>

or **JETA Molecular**

Krommewetering 101C

3543 AN, Utrecht

The Netherlands

Tel: +31 (0)6 54 13 66 97

Email: info@jetabv.com


<https://www.jetamolecular.com>

Genotyping Test


The QTRACE® Analysis System

The DNAs that comprise a mixed DNA sample are analyzed using a QTRACE® Genotyping Plate and QTRACE® qPCR Master Mix, to identify all of the informative assays for the samples. An informative assay is an assay for a marker allele that is present (positive) in one individual genome *and* absent (negative) in the other genome.


The genotyping plate contains a duplicate set of assays: the 46 quantification assays and the reference (RNaseP) assay that serves as both a positive control and a No Template Control (NTC), pre-dispensed and dried in an optical qPCR plate. The QTRACE® qPCR Master Mix comes supplied with dUTP and UNG for built-in carryover contamination control.



The multiplexed version of the genotyping plate (MultiTRACE™) contains a quadruplicate set of assays: the 45 quantification assays and the reference (RNaseP) assay that serves as both a positive control and a No Template Control (NTC), pre-dispensed and dried in an optical qPCR plate.




The Extended Panel genotyping plate can be used in rare cases including patients receiving a second transplant or highly similar siblings, where there may be need for additional markers to the standard QTRACE genotyping plate. It contains a duplicate set of 34 quantification assays and the reference (RNaseP) assay that serves as both a positive control and a No Template Control (NTC), pre-dispensed and dried in an optical qPCR plate.




The DigitalTRACE™ Analysis System

The DNAs that comprise a mixed DNA sample are analyzed using a DigitalTRACE™ or MultiTRACE™ Genotyping Plate, to identify all of the informative assays for the samples. An informative assay is an assay for a marker allele that is present (positive) in one individual genome and absent (negative) in the other genome.

The DigitalTRACE™ Genotyping Plate contains a set of 43 quantification assays and the reference (RNase P) assay that serves as both a positive control and a No Template Control (NTC). The assay sets are pre-arrayed in rows.



The DigitalTRACE™ EP QIAcuity Genotyping Plate represents an extended panel of dPCR markers that can be used for dPCR monitoring in the case of need for additional markers allowing to distinguish between donor and recipient DNA. This plate contains a set of 27 quantification assays and the reference (RNase P) assay that serves as both a positive control and a No Template Control (NTC). The assay sets are pre-arrayed in columns.



Alternatively, the MultiTRACE™ Genotyping Plates can be used for a dPCR genotyping. The MultiTRACE™ Genotyping Plates contain the quantification assays in FAM and the reference (RNase P) assay that serves as both a positive control and a No Template Control (NTC) in HEX.

Protocol

Before setting up an experiment, select the instrument as well as the type and version of the genotyping plate used.



Press the “Preferences” button and choose from the available instruments and genotyping plate types.


The QTRACE® Analysis System

To generate a new record in TRACE Analysis™ Software, in the Section labeled “Recipient” enter the Recipient Name, Recipient Identifier and a unique Sample Identifier. While “Date of Birth” is an optional field for all samples, a “Date of Transplant” must be entered, if you ultimately want to have the data stored and reported in a temporal manner.

Recipient	
Recipient Name	<input type="text"/>
Recipient ID	<input type="text"/>
Sample ID	<input type="text"/>
Concentration ng/ul	20
Date of Birth	<input type="text" value="XX-XX-XXXX"/> [13]
Date of Transplant	<input type="text" value="XX-XX-XXXX"/> [13]
Gender	<input checked="" type="radio"/> Male <input type="radio"/> Female
Comment	<input type="text"/>
Disease Type	<input type="text"/>
Donor	
Donor Name	<input type="text"/>
Donor ID	<input type="text"/>
Sample ID	<input type="text"/>
Concentration ng/ul	20
Gender	<input checked="" type="radio"/> Male <input type="radio"/> Female
Comment	<input type="text"/>


For a sample which should be genotyped against the Recipient sample, enter the Donor Name, Donor Identifier and unique Sample Identifier. You can comparatively genotype up to two samples on a single standard plate or up to four samples on a single multiplexed plate using TRACE Analysis™ Software. TRACE Analysis™ Software also allows you to virtually compare multiple samples, independent of when they were genotyped. Based on the data from verification studies, JETA Molecular recommends the use of 5ng DNA input per well for genotyping. (DNA inputs are customized in the software's [Preferences](#) for Concentrations).

Once all samples to be comparatively genotyped are entered into the Sample window, press the

“Add Typing Samples”  button to add the samples to the plate set up file. You will see the samples now added to the 96 well plate in the middle of the screen.

When using a multiplexed version of the QTRACE® Genotyping Plate, you may add the four DNA samples to be genotyped as two couples Recipient/Donor (added to the plate in two steps, the second couple after the first one). In this case, TRACE Analysis™ Software generates a separate genotyping report for each couple Recipient/Donor.

Alternatively, the four DNA samples may be added to the plate (in one single step) as one Recipient and three Donors, using the Add Donor button:



In this case, TRACE Analysis™ Software generates a single genotyping report containing all four DNA samples genotyped.


Sample Layout View for 2 Samples being genotyped in the standard QTRACE® genotyping plate (left) and in the multiplexed genotyping plate (right):

	1	2	3	4	5	6	7	8	9	10	11	12
A	James	James	James	James	James	James	James	James	James	James	James	James
B	James	James	James	James	James	James	James	James	James	James	James	James
C	James	James	James	James	James	James	James	James	James	James	James	James
D	James	James	James	James	James	James	James	James	James	James	James	James
E	James	James	James	James	James	James	James	James	James	James	James	James
F	James	James	James	James	James	James	James	James	James	James	James	James
G	James	James	James	James	James	James	James	James	James	James	James	James
H	James	James	James	James	James	James	James	James	James	James	James	James

	1	2	3	4	5	6	7	8	9	10	11	12
A	Brian	Brian	Brian	Only	Only	Only	Jim	Jim	Jim	Jim	Jim	Jim
B	Brian	Brian	Brian	Only	Only	Only	Jim	Jim	Jim	Jim	Jim	Jim
C	Brian	Brian	Brian	Only	Only	Only	Jim	Jim	Jim	Jim	Jim	Jim
D	Brian	Brian	Brian	Only	Only	Only	Jim	Jim	Jim	Jim	Jim	Jim
E	Brian	Brian	Brian	Only	Only	Only	Jim	Jim	Jim	Jim	Jim	Jim
F	Brian	Brian	Brian	Only	Only	Only	Jim	Jim	Jim	Jim	Jim	Jim
G	Brian	Brian	Brian	Only	Only	Only	Jim	Jim	Jim	Jim	Jim	Jim
H	Brian	Brian	Brian	Only	Only	Only	Jim	Jim	Jim	Jim	Jim	Jim

The colors in the small plate images (Sample View and Assay View) at the bottom of the window are enabled in the large plate image, by touching the colored plate image of interest. By touching the same image again, the large plate colors disappear. (For Genotyping, the assays are already in the wells and are not added by the operator).

Once the samples have been added to your plate, and the experiment name and operator-ID

have been entered, press the “Export Setup to PCR”  button.

Browse to the location where you want the file saved and name it as you wish. This file can then be imported into your qPCR machine's QTRACE® template to execute the qPCR analysis.

Once the file is saved, TRACE Analysis™ Software generates a protocol, based on the experimental inputs and the settings in the preferences menu.

The following protocol is an example output from TRACE Analysis™ Software for the Genotyping procedure using a [standard genotyping plate](#):

Set up all reactions in a pre-PCR lab, under ambient conditions without ice.

Open a QTRACE® Genotyping Plate Pack and remove the genotyping plate.

Label the genotyping plate with the genotyping test name.

Collect the two DNA samples to be screened, as well as QTRACE® qPCR Master Mix and de-ionized H₂O. Briefly vortex and centrifuge all tubes before opening

Label three, 1.5 mL microcentrifuge tubes:

Sample 1 Mix

Sample 2 Mix

NTC Mix

For each sample to be genotyped, prepare a qPCR Master Mix containing sample DNA, de-ionized H₂O and QTRACE® qPCR Master Mix as suggested by QTRACE® Software in Table 1. A No Template Control (NTC) Mix is prepared with de-ionized H₂O and QTRACE® qPCR Master Mix.

Table 1. qPCR Master Mix Composition

<u>Sample 1 qPCR Mix</u>	<u>1X</u>	<u>50X</u>
QTRACE® qPCR Master Mix	5.0ul	250.0ul
"recipient ID" DNA	0.8ul	40.0ul
H ₂ O	<u>19.2ul</u>	<u>960.0ul</u>
	25.0ul	1250.0ul

<u>Sample 2 qPCR Mix</u>	<u>1X</u>	<u>50X</u>
QTRACE® qPCR Master Mix	5.0ul	250.0ul
"Donor ID" DNA	0.1ul	5ul
H ₂ O	<u>19.9ul</u>	<u>955ul</u>
	25.0ul	1250.0ul

<u>NTC qPCR Mix</u>	<u>1X</u>	<u>3X</u>
QTRACE® qPCR Master Mix	5.0ul	15.0ul
H ₂ O	<u>20.0ul</u>	<u>60.0ul</u>
	25.0ul	75.0ul

Vortex each tube to thoroughly mix the contents and centrifuge briefly to collect the reaction mix at the bottom of the tube.

Remove the adhesive cover from the genotyping plate.

Refer to the DNA Sample Layout Plate Layout at the end of the protocol for pipetting the following:

Dispense 25 µl of the Sample 1 Mix into Wells A1-D6 and in Well D11 of the genotyping plate.

Dispense 25 µl of the Sample 2 Mix into Wells E1-H6 and in Well H11 of the genotyping plate. Dispense 25 µl of the 5X PCR Master Mix/NTC mixture to wells D12 and H12.

A repeat pipettor is recommended to minimize pipetting repetition and increase accuracy.

Visually inspect plate wells from the sides and bottom to confirm consistent volume. - *For use with qPCR machines which do not accept standard ABI MicroAmp Optical Plates:* Transfer the genotyping assay reactions to a 96-well plate for your instrument using a multichannel pipette.

Seal the plate completely with MicroAmp® Optical Adhesive Film using the MicroAmp® Adhesive Film Applicator.

IMPORTANT! Vortex the plate to mix the contents of each well. Centrifuge the plates briefly using a plate centrifuge to collect the contents at the bottom of the wells.

Load the plate into your qPCR machine. Open a pre-configured QTRACE® qPCR template and save the file with a unique name. (If you don't have a template, please see thermal cycling profiles below.

Import the Sample Setup sheet generated by TRACE Analysis™ Software.

Save the file and start the qPCR run.

The following protocol is an example output from TRACE Analysis™ Software for the Genotyping procedure using a multiplexed genotyping plate:

Set up all reactions in a pre-PCR lab, under ambient conditions without ice.

Open a MultiTRACE™ Genotyping Plate Pack and remove the genotyping plate.

Label the genotyping plate with the genotyping test name.

Collect the four DNA samples to be screened, as well as QTRACE® qPCR Master Mix and de-ionized H₂O. Briefly vortex and centrifuge all tubes before opening

Label five, 1.5 mL microcentrifuge tubes:

- Sample 1 Mix
- Sample 2 Mix
- Sample 3 Mix
- Sample 4 Mix
- NTC Mix

For each sample to be genotyped, prepare a qPCR Master Mix containing sample DNA, de-ionized H₂O and QTRACE® qPCR Master Mix as suggested by TRACE Analysis™ Software in Table 2. A No Template Control (NTC) Mix is prepared with de-ionized H₂O and QTRACE® qPCR Master Mix.

Table 2. qPCR Master Mix Composition – MultiTRACE™ plates

Sample 1 qPCR Mix	1X	25X
QTRACE® qPCR Master Mix	5.0ul	125.0ul
"Recipient 1 ID"		
DNA	0.3ul	7.5ul
H ₂ O	<u>19.7ul</u>	<u>492.5ul</u>
	25.0ul	625.0ul


Sample 2 qPCR Mix	1X	25X
QTRACE® qPCR Master Mix	5.0ul	125.0ul
"Donor 1 ID" DNA	0.2ul	5.4ul
H ₂ O	<u>19.8ul</u>	<u>494.6ul</u>
	25.0ul	625.0ul

Sample 3 qPCR Mix	1X	25X
QTRACE® qPCR Master Mix	5.0ul	125.0ul
"Recipient 2 ID"		
DNA	0.4ul	8.9ul
H ₂ O	<u>19.6ul</u>	<u>491.1ul</u>
	25.0ul	625.0ul


Sample 2 qPCR Mix	1X	25X
QTRACE® qPCR Master Mix	5.0ul	125.0ul
"Donor 2 ID" DNA	0.1ul	2.5ul
H ₂ O	<u>19.9ul</u>	<u>497.5ul</u>
	25.0ul	625.0ul

Thermal Cycling Protocol for QTRACE® System Products

The QTRACE® System will deliver optimal results when the following thermal profile is used in qPCR.



For labs which require use of molecular tests with built-in carryover contamination control, the QTRACE® qPCR Master Mix contains dUTP and Uracil N-Glycosylase (UNG) enzyme. In order to use this enhanced capability of the master mix, add a 2 minute hold at 50°C, prior to the enzyme activation/initial denaturation step.



The DigitalTRACE™ Analysis System


Genotyping Test Protocol - QIAcuity

Change the instrument type in the Preferences of the TRACE Analysis™ Software to QIAcuity and Plate type to QIAcuity, v1.


To generate a new record in TRACE Analysis™ Software, in the Section labeled “Recipient” enter the Recipient Name, Recipient Identifier and a unique Sample Identifier. While “Date of Birth” is an optional field for all samples, a “Date of Transplant” must be entered, if you ultimately want to have the data stored and reported in a temporal manner.

For a sample which should be genotyped against the Recipient sample, enter the Donor Name, Donor Identifier and unique Sample Identifier. You can genotype up to eight samples on a single plate using TRACE Analysis™ Software. TRACE Analysis™ Software also allows you to virtually compare multiple samples, independent of when they were genotyped. Based on the data from verification studies, JETA Molecular recommends the use of 10 ng DNA input per well for genotyping. (DNA inputs are customized in the software's [Preferences](#) for Concentrations). Enter the concentrations of your samples to let the software calculate with.

After all samples to be comparatively genotyped are entered into the Sample window, press the

“Add Typing Samples”  button to add the samples to the plate set up file. You will see the samples now added to the 96-well plate in the middle of the screen.

Sample Layout View for 8 Samples being genotyped in the DigitalTRACE™ QIAcuity plate:



	1	2	3	4	5	6	7	8	9	10	11	12
A	Tom 001	Tom 001	Tom 001	Tom 001	Tom 001	Tom 001	Tom 001	Tom 001	Tom 001	Tom 001	Tom 001	NTC
B	Alex 002	Alex 002	Alex 002	Alex 002	Alex 002	Alex 002	Alex 002	Alex 002	Alex 002	Alex 002	Alex 002	NTC
C	Julie 003	Julie 003	Julie 003	Julie 003	Julie 003	Julie 003	Julie 003	Julie 003	Julie 003	Julie 003	Julie 003	NTC
D	James 004	James 004	James 004	James 004	James 004	James 004	James 004	James 004	James 004	James 004	James 004	NTC
E	Peter 005	Peter 005	Peter 005	Peter 005	Peter 005	Peter 005	Peter 005	Peter 005	Peter 005	Peter 005	Peter 005	NTC
F	Barbara 006	Barbara 006	Barbara 006	Barbara 006	Barbara 006	Barbara 006	Barbara 006	Barbara 006	Barbara 006	Barbara 006	Barbara 006	NTC
G	John 007	John 007	John 007	John 007	John 007	John 007	John 007	John 007	John 007	John 007	John 007	NTC
H	Ashley 008	Ashley 008	Ashley 008	Ashley 008	Ashley 008	Ashley 008	Ashley 008	Ashley 008	Ashley 008	Ashley 008	Ashley 008	NTC

The colors in the small plate images (Sample View and Assay View) at the bottom of the window are enabled in the large plate image, by touching the colored plate image of interest. By touching the same image again, the large plate colors disappear. (For Genotyping, the assays are already in the wells and are not added by the operator).

Once the samples have been added virtually to your plate, and the experiment name have been

entered, press the “Export Setup to PCR”  button.

Browse to the location where you want the file saved and name it as you wish. This file can then be imported into your dPCR machine’s DigitalTRACE template to execute the dPCR analysis.

Once the file is saved, TRACE Analysis™ Software generates a protocol, based on the experimental inputs and the settings in the Preferences menu. Print out this protocol.

1. Set up all reactions in a pre-PCR lab, under ambient conditions without ice.
2. Collect all DNA samples to be screened, as well as QIAcuity typing plate together with 4x Probe PCR Master Mix and de-ionized H₂O.
3. Briefly vortex and centrifuge all tubes before opening.
4. For each sample to be genotyped, label a tube and a Mix containing sample DNA, de-ionized H₂O and 4x Probe PCR Master Mix as suggested by TRACE Analysis™ Software in Table 1. A No Template Control (NTC) Mix is prepared with de-ionized H₂O and 4x Probe PCR Master Mix.

* - Make a Ten-Fold Dilution (1:10) of Sample

Table 1. Master Mix Composition

Sample 1 dPCR Mix	1 x	13 x
4x Probe PCR Master Mix	3.0 µl	39.0 µl
001 DNA	0.5 µl	*6.5 µl
H ₂ O	8.5 µl	110.5 µl
	12.0 µl	156.0 µl

Sample 2 dPCR Mix	1 x	13 x
4x Probe PCR Master Mix	3.0 µl	39.0 µl
002 DNA	0.5 µl	*6.5 µl
H ₂ O	8.5 µl	110.5 µl
	12.0 µl	156.0 µl


NTC dPCR Mix	1 x	10 x
4x Probe PCR Master Mix	3.0 µl	30.0 µl
H ₂ O	9.0 µl	90.0 µl
	12.0 µl	120.0 µl

5. Vortex each tube to thoroughly mix the contents and centrifuge briefly to collect the reaction mix at the bottom of the tube.
6. Remove the adhesive cover from the genotyping plate.
7. Deliver 13.2 µl of each Sample Mix and NTC mix to the Typing plate as defined in TRACE Analysis™ Software’s Assay Layout view.
8. An automated multichannel pipette is recommended in this step to minimize pipetting repetition and increase accuracy.
9. Visually inspect plate wells from the sides and bottom to confirm consistent volume.
10. Seal the plate with an Adhesive Film.
11. **IMPORTANT! Vortex the plate to mix the contents of each well.** Centrifuge the plates briefly using a plate centrifuge to collect the contents at the bottom of the wells.
12. Remove the Adhesive cover very carefully.
13. Transfer 12 µl of each prepared reaction mix into a single column of a 96-well 8.SK Nanoplate. Seal the Nanoplate with the compatible plate sealer.
14. Load the Nanoplate into the QIAcuity digital PCR system.
15. Launch the QIAcuity Software Suite.
16. Open the DigitalTRACE typing template and import the Sample Setup sheet generated by TRACE Analysis™ Software.
17. Save the file and start the run.


Experiment Setup in QIAcuity Suite Software, v2.5

Create a new QIAcuity Plate by selecting "New Plate".

Load a template by selecting "Plate templates..." and import the appropriate DigitalTRACE Genotyping template:




Press "Save changes". The new plate will appear in the main window of the QIAcuity Software Suite:



Click on the plate name to open the plate configuration procedure. Type in a new plate name and save the changes.


Import the Sample Setup sheet (.csv) generated by TRACE Analysis™ Software by selecting "Plate layout" tab and "CSV import/export".

Select "Import from CSV".



Import the plate setup file by clicking on "Import".

Inspect if all your selected samples are represented on the Plate figure.



Select "Done" to finish the setup.
Start the run.

Thermal Cycling and Imaging Protocol for QIAcuity

The DigitalTRACE™ System will deliver optimal results when the following thermal profile is used in dPCR:

Number of cycles	Temperature	Time
1	95 °C	3 min
40	95 °C	15 s
	60 °C	30 s

Use the following imaging parameters:

Channel	Exposure duration	Gain
Green	500 ms	6
Yellow	500 ms	6
Orange	Off	Off
Red	300 ms	4
Crimson	400 ms	4

Genotyping Test Protocol - Biorad

Change the instrument type in the Preferences of the TRACE Analysis™ Software to Biorad QX-200 and Plate type to MultiTRACE, v3.

To generate a new record in TRACE Analysis™ Software, in the Section labeled “Recipient” enter the Recipient Name, Recipient Identifier and a unique Sample Identifier. While “Date of Birth” is an optional field for all samples, a “Date of Transplant” must be entered, if you ultimately want

to have the data stored and reported in a temporal manner.


Recipient	
Recipient First Name	
Recipient Last Name	
Recipient ID	
Sample ID	
Concentration ng/ul	100
Date of Birth	XX-XX-XXXX [15]
Date of Transplant	XX-XX-XXXX [15]
Gender	<input type="radio"/> Male <input type="radio"/> Female
Comment	
Disease Type	

Donor	
Donor First Name	
Donor Last Name	
Donor ID	
Sample ID	
Concentration ng/ul	100
Date of Birth	XX-XX-XXXX [15]
Gender	<input type="radio"/> Male <input type="radio"/> Female
Comment	

For a sample which should be genotyped against the Recipient sample, enter the Donor Name, Donor Identifier and unique Sample Identifier. You can comparatively genotype up to four samples on a single plate using TRACE Analysis™ Software. TRACE Analysis™ Software also allows you to virtually compare multiple samples, independent of when they were genotyped. Based on the data from verification studies, JETA Molecular recommends the use of 10 ng DNA input per well for genotyping. (DNA inputs are customized in the software's [Preferences](#) for Concentrations). Enter the concentrations of your samples to let the software calculate with. Once all samples to be comparatively genotyped are entered into the Sample window, press the “Screen” button to add the samples to the plate set up file. You will see the samples now added to the 96-well plate in the middle of the screen.

Sample Layout View for 4 Samples being genotyped in the multiplexed plate:

	1	2	3	4	5	6	7	8	9	10	11	12
A	Tim 1001 Mix028	Tim 1001 Mix036	Tim 1001 Mix044	Bert 1002 Mix028	Bert 1002 Mix036	Bert 1002 Mix044	Ken 1003 Mix028	Ken 1003 Mix036	Ken 1003 Mix044	Jenny 1004 Mix028	Jenny 1004 Mix036	Jenny 1004 Mix044
B	Tim 1001 Mix029	Tim 1001 Mix037	Tim 1001 Mix045	Bert 1002 Mix029	Bert 1002 Mix037	Bert 1002 Mix045	Ken 1003 Mix029	Ken 1003 Mix037	Ken 1003 Mix045	Jenny 1004 Mix029	Jenny 1004 Mix037	Jenny 1004 Mix045
C	Tim 1001 Mix030	Tim 1001 Mix038	Tim 1001 Mix046	Bert 1002 Mix030	Bert 1002 Mix038	Bert 1002 Mix046	Ken 1003 Mix030	Ken 1003 Mix038	Ken 1003 Mix046	Jenny 1004 Mix030	Jenny 1004 Mix038	Jenny 1004 Mix046
D	Tim 1001 Mix031	Tim 1001 Mix039	Tim 1001 Mix047	Bert 1002 Mix031	Bert 1002 Mix039	Bert 1002 Mix047	Ken 1003 Mix031	Ken 1003 Mix039	Ken 1003 Mix047	Jenny 1004 Mix031	Jenny 1004 Mix039	Jenny 1004 Mix047
E	Tim 1001 Mix032	Tim 1001 Mix040	Tim 1001 Mix048	Bert 1002 Mix032	Bert 1002 Mix040	Bert 1002 Mix048	Ken 1003 Mix032	Ken 1003 Mix040	Ken 1003 Mix048	Jenny 1004 Mix032	Jenny 1004 Mix040	Jenny 1004 Mix048
F	Tim 1001 Mix033	Tim 1001 Mix041	Tim 1001 Mix049	Bert 1002 Mix033	Bert 1002 Mix041	Bert 1002 Mix049	Ken 1003 Mix033	Ken 1003 Mix041	Ken 1003 Mix049	Jenny 1004 Mix033	Jenny 1004 Mix041	Jenny 1004 Mix049
G	Tim 1001 Mix034	Tim 1001 Mix042	Tim 1001 Mix050	Bert 1002 Mix034	Bert 1002 Mix042	Bert 1002 Mix050	Ken 1003 Mix034	Ken 1003 Mix042	Ken 1003 Mix050	Jenny 1004 Mix034	Jenny 1004 Mix042	Jenny 1004 Mix050
H	Tim 1001 Mix035	Tim 1001 Mix043	NTC RNaseP	Bert 1002 Mix035	Bert 1002 Mix043	NTC RNaseP	Ken 1003 Mix035	Ken 1003 Mix043	NTC RNaseP	Jenny 1004 Mix035	Jenny 1004 Mix043	NTC RNaseP

Once the samples have been added virtually to your plate, and the experiment name have been entered, press the “Export Setup to PCR”  button.

Browse to the location where you want the file saved and name it as you wish. This file can then be imported into your dPCR machine’s DigitalTRACE template to execute the dPCR analysis. Once the file is saved, TRACE Analysis™ Software generates a protocol, based on the experimental inputs and the settings in the preferences menu. Print out this protocol.

1. Set up all reactions in a pre-PCR lab, under ambient conditions without ice.
2. Open a MultiTRACE® Genotyping Plate Pack and remove the genotyping plate.
3. Label the genotyping plate with the genotyping test name.
4. Collect the four DNA samples to be screened, as well as 2x ddPCR Supermix Master Mix and de-ionized H₂O. Briefly vortex and centrifuge all tubes before opening.
5. Label five 1.5 ml microcentrifuge tubes:
 - a) Sample 1 Mix
 - b) Sample 2 Mix
 - c) Sample 3 Mix
 - d) Sample 4 Mix
 - e) NTC Mix
6. For each sample to be genotyped, prepare a ddPCR Mix containing sample DNA, de-ionized H₂O and 2x ddPCR Supermix as suggested by TRACE Analysis™ Software in Table 1. A No Template Control (NTC) Mix is prepared with de-ionized H₂O and 2x ddPCR Supermix.

Table 1. Master Mix Composition

Sample 1 dPCR Mix	1 x	27 x
2x ddPCR Supermix	11,0 µl	297,0 µl
1001 DNA	0,5 µl	13,5 µl
H ₂ O	10,5 µl	283,5 µl
	22,0 µl	594,0 µl
Sample 2 dPCR Mix	1 x	27 x
2x ddPCR Supermix	11,0 µl	297,0 µl
1002 DNA	0,5 µl	13,5 µl
H ₂ O	10,5 µl	283,5 µl
	22,0 µl	594,0 µl
Sample 3 dPCR Mix	1 x	27 x
2x ddPCR Supermix	11,0 µl	297,0 µl
1003 DNA	0,5 µl	13,5 µl
H ₂ O	10,5 µl	283,5 µl
	22,0 µl	594,0 µl
Sample 4 dPCR Mix	1 x	27 x
2x ddPCR Supermix	11,0 µl	297,0 µl
1004 DNA	0,5 µl	13,5 µl
H ₂ O	10,5 µl	283,5 µl
	22,0 µl	594,0 µl
NTC dPCR Mix	1 x	5 x

2x ddPCR Supermix	11,0 µl	55,0 µl
H ₂ O	11,0 µl	55,0 µl
	22,0 µl	110,0 µl

7. Vortex each tube to thoroughly mix the contents and centrifuge briefly to collect the reaction mix at the bottom of the tube.
8. Remove the adhesive cover from the genotyping plate.
9. Dispense 22 µl of the Sample 1 Mix into Wells A1-G3 by columns of the genotyping plate.
10. Dispense 22 µl of the Sample 2 Mix into Wells A4-G6 by columns of the genotyping plate.
11. Dispense 22 µl of the Sample 3 Mix into Wells A7-G9 by columns of the genotyping plate.
12. Dispense 22 µl of the Sample 4 Mix into Wells A10-G12 by columns of the genotyping plate.
13. Dispense 22 µl of the 5X PCR Master Mix/NTC mixture to wells H3, H6, H9, H12.
14. **A repeat pipettor is recommended to minimize pipetting repetition and increase accuracy.**
15. Refer to the DNA Sample Layout Plate Layout at the end of the protocol.
16. Visually inspect plate wells from the sides and bottom to confirm consistent volume.
17. Seal the plate completely with MicroAmp® Optical Adhesive Film using the MicroAmp® Adhesive Film Applicator.
18. **IMPORTANT! Vortex the plate to mix the contents of each well.** Centrifuge the plates briefly using a plate centrifuge to collect the contents at the bottom of the wells.
19. **Remove the Adhesive cover very carefully.**
20. **!!!The following steps from 21 to 26 are only for manual droplet generator users:**
21. Transfer 20 µl of each prepared sample to the sample wells (middle row) of the DG8 cartridge.
22. Add 70 µl of droplet generation oil to each oil well of the DG8 cartridge.
23. Hook the gasket over the cartridge holder using the holes on both sides.
24. Load the cartridge in the QX200 droplet generator.
25. When droplet generation is complete, remove the disposable gasket from the holder and discard it.
26. Pipet 40 µl of the contents of the droplets into a single column of a 96-well PCR plate.
27. Seal the PCR plate with foil plate seals that are compatible with the PX1 PCR plate sealer and the needles in the QX200 droplet reader.
28. Place the plate into the thermal cycler for PCR amplification.
29. Load the plate after amplification into QX200 droplet reader.
30. Import the Sample Setup sheet generated by TRACE Analysis™ Software.
31. Save the file and start the droplet reader.


Thermal Cycling Protocol for Biorad

The DigitalTRACE™ System will deliver optimal results when the following thermal profile is used in dPCR:



Number of cycles	Temperature	Time	Ramp Rate
1	95 °C	10 min	2 °C/s
40	94 °C	30 s	2 °C/s
	59 °C	60 s	2 °C/s
1	98 °C	10 min	2 °C/s

Biorad Droplet Reader Setup in QuantaSoft

To read the signal after the PCR cycling was completed, setup an experiment in QuantaSoft. Load a template by clicking on Load




Select QS Setup CSV Files as file format



Start the droplet reading run.

Genotyping Data Analysis and Report


The QTRACE® Analysis System

Import qPCR data by clicking the “Import qPCR Data”  button .

Browse to the location of your exported qPCR data file and select it. TRACE Analysis™ will import the data.

TRACE Analysis™ will perform a quality analysis on the data and will present the data in the plate view. The image below shows a genotyping result for a standard (left) and multiplexed (right) genotyping plate. There are three quality scores given to genotyping data: 1) positive, 2) negative and 3) atypical. These values are represented in the plate image by three different

colors:   




Dark green sections represent data scored as truly positive for a sample.

Light blue sections represent data scored as truly negative for a sample.

Red sections represent data scored as atypical for a sample. (Atypical assay results will exclude the assay from consideration as a potentially informative assay for all samples grouped in the analysis).

For more information on the Genotyping Algorithm used by TRACE Analysis™ Software, go to the [Data Analysis Algorithms](#) Section of this manual.

After inspecting the quality of the data, pressing the “Calculate”  button makes TRACE Analysis™ perform comparative genotyping analysis. It will determine and display markers which are informative for all samples in a group.

TRACE Analysis™ displays each assay which was informative for a single sample in the group being compared.

TRACE Analysis™ also displays the chromosomal location of the informative assays, as well as the positive or negative status of the assays for visual inspection.

The report generated from a QTRACE® or DigitalTRACE™ Genotyping experiment may be sorted to provide a custom view of the data.

ReportWindow

Experiment report

Scope:

☒ Entire experiment
 ☐ Transplantation

Samuel Haskell

Format:

☒ Full
 ☐ Summarized

Save as

Column

Order

Sort by:

Informative

Descending

Then by:

Marker

Locus

Informative for

InfoCq

Delta Cq

Used Assays


Assay	Locus	Informative for	InfoCq/ΔCq	Recipient 1	Cq	Donor 1	Cq	
824	22b	Samuel Haskell	27.51	0.2	Positive	27.5	Negative	40.0
401	22b	Samuel Haskell	28.48	0.75	Positive	28.4	Negative	40.0
504	17p	Samuel Haskell	29.18	1.44	Positive	29.1	Negative	40.0
434	3b	Samuel Haskell	29.96	2.26	Positive	30.0	Negative	40.0
514	22b	Mut. Haskell	24.89	-	Negative	40.0	Positive	24.9
			0.84					
148	17p	Mut. Haskell	27.29	-	Negative	40.0	Positive	27.3
			0.23					
854	12b	Mut. Haskell	27.89	0.17	Negative	40.0	Positive	27.7
519	8b	Mut. Haskell	27.98	0.40	Negative	40.0	Positive	28.0
489	7b	Mut. Haskell	28.09	0.17	Negative	40.0	Positive	28.1
749	18b	Mut. Haskell	28.27	0.75	Negative	40.0	Positive	28.2
173	18b	Mut. Haskell	28.31	0.79	Negative	40.0	Positive	28.3
312	18b	Mut. Haskell	28.42	0.9	Negative	40.0	Positive	28.4
415	2b	Mut. Haskell	28.48	0.96	Negative	40.0	Positive	28.5
425	6b	Mut. Haskell	28.56	1.01	Negative	40.0	Positive	28.6
235	2b	Mut. Haskell	28.55	1.49	Negative	40.0	Positive	28.5
137	2b	Mut. Haskell	28.58	1.53	Negative	40.0	Positive	28.2
240	8b	None	-	-	Negative	40.0	Negative	40.0
832	12b	None	-	-	Negative	40.0	Negative	40.0
148	3b	None	-	-	Positive	27.8	Positive	28.4
531	17p	None	-	-	Negative	40.0	Negative	40.0
395	3b	None	-	-	Negative	40.0	Negative	40.0
774	7b	None	-	-	Negative	40.0	Negative	40.0
545	7b	None	-	-	Negative	40.0	Negative	40.0

Type text to find...

The DigitalTRACE™ Analysis System

Genotyping Data Analysis and Report - QIAcuity

After the QIAcuity run has finished, check in the QIAcuity Software Suite if the automatic thresholds are correct. Adjust the thresholds manually if needed.




Select all wells on the plate and analyze per target (not per channel), export data by selecting Export to CSV.

List Signalmap Heatmap Histogram 1D Scatterplot 2D Scatterplot Concentration diagram

☐ Add to report ☐ Show mean values for replicates ☒ Export to CSV

Sample/NTC/Control	Reaction Mix	Target	IC	Control type	Concentration * copies/μL	CI (95%)	Partitions valid	positive	negative	Threshold
A1	c195 221010	d681	-	-	0.000	-	8191	0	8191	30.86
		d971	-	-	0.372	147.5%	8191	1	8190	53.81
		d113	-	-	0.000	-	8191	0	8191	20.66
		d597	-	-	149.4	9.9%	8191	392	7799	24.86
A2	c195 221010	d777	-	-	0.376	147.5%	8271	1	8270	21.93
		d396	-	-	0.000	-	8271	0	8271	34.68
		d892	-	-	0.000	-	8271	0	8271	21.93
		d333	-	-	180.2	9.1%	8271	466	7805	24.23


Import dPCR data to TRACE Analysis™ Software by clicking the “Import PCR Data”  button. Browse to the location of your exported dPCR data file and select it.

TRACE Analysis™ Software will perform a quality analysis on the data and will present the data in the plate view.

There are three quality scores given to genotyping data: 1) positive (green), 2) negative (blue) and 3) atypical (red).

These values are represented accordingly in the plate image by three different colors:





An atypical assay results will exclude the assay from consideration as a potentially informative assay for all samples grouped in the analysis.

After inspecting the quality of the data, pressing the “Calculate” **Calculate** button makes TRACE Analysis™ perform comparative genotyping analysis. It will determine and display markers which are informative for all samples in a group.

Press the "Report" **Report...** button to generate the Genotyping Report

TRACE Analysis™ Software displays each assay which was informative for a single sample in the group being compared, and it also displays the chromosomal location of the informative assays, as well as the positive or negative status of the assays for visual inspection.

The report generated from a TRACE Analysis™ Genotyping experiment may be sorted to provide a custom view of the data.

Experiment report

Scopes: ☒ Entire experiment
☐ Transplantation Tom
Format: ☒ Full
☐ Summarized

Save as: [icon] [icon] [icon] [icon] [icon]

Sort by: Locus Descending
Then by: Marker Ascending

Recipient 1: Tom
ID: 45468
Gender: Unknown
Disease type: -
Date of birth: -
Transplant date: -
Comments: -

Donor 1: Alex
ID: 654654
Gender: Unknown
Date of birth: -
Comments: -


Used Assays

Assay	Locus	Informative for	Recipient	CNV	Concen	Pos	Partition	Donor 1	CNV	Concen	Pos	Partitions
#221	Xq	None	Tom	Positive	0.9	161.45	424	#277	Positive	1.0	211.70	148
#333	3p	Tom	Tom	Positive	1.0	180.20	466	#275	Negative	0.0	0.00	0
#472	3p	Tom	Tom	Positive	0.8	141.90	369	#247	Negative	0.0	0.79	2
#113	8q	None	Tom	Negative	0.0	0.00	0	#261	Negative	0.0	0.00	0
#252	8q	None	Tom	Negative	0.0	0.00	0	#214	Negative	0.0	0.00	0
#336	8q	None	Tom	Negative	0.0	0.00	0	#224	Negative	0.0	0.00	0
#136	7q	Tom	Tom	Positive	1.1	187.80	482	#260	Negative	0.0	0.00	0
#990	7q	Tom	Tom	Positive	1.1	189.90	487	#239	Negative	0.0	0.00	0
#417	7p	None	Tom	Positive	1.1	184.90	500	#239	Positive	1.1	203.30	633
#892	6q	None	Tom	Negative	0.0	0.00	0	#271	Negative	0.0	0.00	0
#971	6q	None	Tom	Negative	0.0	0.17	1	#391	Positive	1.0	226.60	576
#222	5p	Tom	Tom	Positive	1.0	171.20	440	#242	Negative	0.0	0.00	0
#795	5q	None	Tom	Negative	0.0	0.10	1	#205	Negative	0.0	1.17	3
#176	3p	None	Tom	Positive	0.8	141.70	369	#277	Positive	0.9	205.90	510
#884	4q	None	Tom	Negative	0.0	0.00	0	#244	Positive	1.1	183.60	5393
#100	2q	Tom	Tom	Positive	1.2	190.10	376	#274	Negative	0.0	0.00	0
#174	2q	None	Tom	Negative	0.0	0.00	0	#238	Positive	1.0	213.50	174
#678	2q	None	Tom	Negative	0.0	0.00	0	#247	Positive	1.1	246.50	596
#885	2p	Tom	Tom	Positive	1.0	182.50	480	#204	Negative	0.0	0.00	0
#984	7q	None	Tom	Negative	0.0	0.10	1	#206	Negative	0.0	1.17	1


To create an anonymized report, go to the [Anonymized Reporting Section](#).

Genotyping Data Analysis and Report - Biorad

After the reading process has finished, check in the QuantaSoft Software if the automatic thresholds are correct. Adjust the thresholds manually if needed.




Select all wells on the plate and export data by selecting Export CSV.

Import dPCR data to TRACE Analysis™ Software by clicking the “Import PCR Data”  button. Browse to the location of your exported dPCR data file and select it.


TRACE Analysis™ Software will perform a quality analysis on the data and will present the data in the plate view.

There are three quality scores given to genotyping data: 1) positive (green), 2) negative (blue) and 3) atypical (red).

These values are represented accordingly in the plate image by three different colors:



An atypical assay results will exclude the assay from consideration as a potentially informative assay for all samples grouped in the analysis.

After inspecting the quality of the data, pressing the “Calculate”  button makes TRACE Analysis™ perform comparative genotyping analysis. It will determine and display markers which are informative for all samples in a group.

Press the "Report"  button to generate the Genotyping Report

TRACE Analysis™ Software displays each assay which was informative for a single sample in the group being compared, and it also displays the chromosomal location of the informative assays, as well as the positive or negative status of the assays for visual inspection.

The report generated from a TRACE Analysis™ Genotyping experiment may be sorted to provide a custom view of the data.

ReportWindow

Experiment report

Scope: ☒ Entire experiment
☐ Transplantation: **Papo Delgado**

Format: ☒ Full
☐ Summarized

Save as: [icon] [icon] [icon] [icon] [icon]

Column: **Marker** Order: **Descending**
 Sort by: **Marker**
 Then by: **Ascending**
 Locus
 Informative for
 Delta Cq

Recipient 1
 Name: **Papo Delgado**
 ID: **Papo001**
 Gender: **Unknown**
 Disease type: **-**
 Date of birth: **-**
 Transplant: **-**
 Date: **-**
 Comments: **-**

Donor 1
 Name: **Cope Delgado**
 ID: **Cope001**
 Gender: **Unknown**
 Date of birth: **-**
 Comments: **-**

Used Assays

Assay	Locus	Informative for	Recipient	CNV	Genotype	Peak	Area	Assayed/Donor 1	CNV	Genotype	Peak	Area	Assayed/Donor 1
P25	146	None	Positive	2.0	61.80	855	17094	Positive	2.0	68.00	2064	18796	
195	74	Cope Delgado	Negative	0.0	0.00	0	12887	Positive	1.1	36.40	491	16665	
987	16	None	Positive	1.0	30.00	304	12381	Positive	0.9	30.50	421	16382	
471	46	None	Negative	0.0	0.00	0	12178	Negative	0.0	0.00	0	17147	
862	176	None	Positive	1.0	67.30	643	17018	Positive	1.1	35.70	487	14264	
316	84	Papo Delgado	Positive	1.0	30.50	384	14811	Negative	0.0	0.00	0	16670	
173	174	Cope Delgado	Negative	0.0	0.00	0	12761	Positive	2.0	74.40	3315	16682	
816	208	Papo Delgado	Positive	1.0	30.40	391	13319	Negative	0.0	0.00	0	17253	
193	46	Cope Delgado	Negative	0.0	0.00	0	12884	Positive	1.0	35.20	524	17794	
884	26	None	Negative	0.0	0.07	1	17097	Negative	0.0	0.07	1	17476	
474	224	Papo Delgado	Positive	1.0	25.40	393	14843	Negative	0.0	0.00	0	17117	
824	276	Cope Delgado	Negative	0.0	0.00	0	13319	Positive	1.1	36.70	530	17253	
399	54	None	Negative	0.0	0.00	0	16650	Negative	0.0	0.04	0	16607	
177	176	None	Positive	1.0	64.80	465	12176	Positive	1.0	34.50	462	17258	
748	276	Papo Delgado	Positive	1.1	32.80	414	13280	Negative	0.0	0.00	0	16894	
721	94	Cope Delgado	Negative	0.0	0.00	0	13188	Positive	1.1	36.00	484	16379	
710	54	None	Positive	1.0	31.10	392	13080	Positive	1.0	30.40	486	16894	
494	46	None	Positive	1.1	18.50	464	16851	Positive	1.1	17.40	150	16683	
481	176	Papo Delgado	Positive	1.0	30.40	407	13880	Negative	0.0	0.00	0	17617	

Type text to find...

Multiple Donor Analysis Using QTRACE Plates

Enter Recipient and Donor 1 Sample Information, as you would normally. Note: the check boxes always indicate which samples will be added to the experiment.

Following the normal setup and export steps, TRACE Analysis™ will make a protocol for typing both samples on one plate.

GENOTYPING PROTOCOL FOR ULLA

- Set up all reactions in a pre-PCR lab, under ambient conditions without ice.
- Open a QTRACE® Genotyping Plate Pack and remove the genotyping plate.
- Label the genotyping plate with the genotyping test name.
- Collect the two DNA samples to be screened, as well as QTRACE® qPCR Master Mix and de-ionized H₂O. Briefly vortex and centrifuge all tubes before opening.
- Label three 1.5 ml microcentrifuge tubes:
 - Sample 1 Mix
 - Sample 2 Mix
 - NTC Mix
- For each sample to be genotyped, prepare a qPCR Master Mix containing sample DNA, de-ionized H₂O and QTRACE® qPCR Master Mix as suggested by QTRACE® Software in Table 1. A No Template Control (NTC) Mix is prepared with de-ionized H₂O and QTRACE® qPCR Master Mix.

Table 1. qPCR Master Mix Composition

Sample 1 qPCR Mix	1 µl	50 µl
QTRACE® qPCR Master Mix	5.0 µl	250.0 µl
Ulla Sample ID DNA	0.1 µl	5.0 µl
H ₂ O	19.9 µl	995.0 µl
	25.0 µl	1250.0 µl

Sample 2 qPCR Mix	1 µl	50 µl
QTRACE® qPCR Master Mix	5.0 µl	250.0 µl
Jonna Sample ID DNA	0.1 µl	5.0 µl
H ₂ O	19.9 µl	995.0 µl
	25.0 µl	1250.0 µl

NTC qPCR Mix	1 µl	50 µl
QTRACE® qPCR Master Mix	5.0 µl	250.0 µl
H ₂ O	20.0 µl	1000.0 µl
	25.0 µl	1250.0 µl

- Vortex each tube to thoroughly mix the contents and centrifuge briefly to collect the reaction mix at the bottom of the tube.

When the data is imported and approved,



TRACE Analysis™ generates a report for the first pair.

Chimerism Genotyping - Full Report

Experiment name: Recipient and Donor 1 Genotyping

Experiment date: 24 March 2017

Data folder: C:\Users\Doug\Desktop\QTRACE\Transplantation


Operator name: Doug

Recipient	Donor1
Name: Ulla	Name: Jonna
ID: Ulla UPN	ID: Jonna UPN
Gender: Female	Gender: Female
Disease type: -	Comments: -
Date of birth: 10 July 2017	
Transplant date: 10 July 2017	
Comments: -	

Used Assays

Assay	Locus	Informative for	InfoCq	BCq	Recipient	CoDonor 1	Cq	
548	1p	Ulla	26.64	0	Positive	26.6	Negative	40.0
548	17q	Ulla	26.89	0.25	Positive	26.8	Negative	40.0
605	14q	Ulla	27.07	0.44	Positive	27.3	Negative	40.0
504	17p	Ulla	27.14	0.51	Positive	27.3	Negative	40.0
406	2q	Ulla	27.19	0.54	Positive	27.3	Negative	40.0
948	5p	Ulla	27.29	0.66	Positive	27.3	Negative	40.0
425	6q	Ulla	27.32	0.69	Positive	27.3	Negative	40.0
634	13q	Ulla	27.43	0.78	Positive	27.4	Negative	40.0
784	10q	Ulla	27.67	0.94	Positive	27.4	Negative	40.0
240	8q	Ulla	27.62	0.99	Positive	27.4	Negative	40.0
768	18q	Ulla	28.19	1.54	Positive	28.2	Negative	40.0
346	18q	Ulla	29.15	2.52	Positive	29.3	Negative	40.0
519	8q	Jonna	25.09	-0.7	Negative	40.0	Positive	25.2
520	20q	Jonna	25.29	-0.51	Negative	40.0	Positive	25.3
326	6q	Jonna	25.77	-0.03	Negative	40.0	Positive	25.8
721	17q	Jonna	26.77	0.48	Negative	40.0	Positive	26.3
819	14q	Jonna	26.59	0.79	Negative	40.0	Positive	26.4
802	7p	Jonna	26.68	0.88	Negative	40.0	Positive	26.7
748	17q	Jonna	26.81	1.01	Negative	40.0	Positive	26.8
359	11p	Jonna	26.85	1.06	Negative	39.6	Positive	26.9
531	17p	Jonna	27.05	1.26	Negative	40.0	Positive	27.1
361	12q	Jonna	27.1	1.3	Negative	40.0	Positive	27.1

When the second donor is to be analyzed, select the recipient's record in the sample search window next to the home button. The samples which have already been genotyped will appear as they did when initially entered into TRACE Analysis™ Software. Notice that the check boxes by default are not checked, as these samples have already been genotyped by the software.

If you click the Add Donor  button, a new set of fields appear for a new sample record to be associated with the same recipient sample.

These Boxes are not checked by default, as their sample typing data is approved

Only this sample will be added to the genotyping experiment, but the data will be associated with this recipient

Enter the new sample information for the second donor, as well as the Date of the Second Transplant. When the Screen button is clicked, only the new sample is added to the experiment.

TRACE Analysis™ will now generate a genotyping protocol for only the second donor.

GENOTYPING PROTOCOL FOR ULLA

- Set up all reactions in a pre-PCR lab, under ambient conditions without ice.
- Open a QTRACE® Genotyping Plate Pack and remove the genotyping plate.
- Label the genotyping plate with the genotyping test name.
- Collect the DNA sample to be screened, as well as QTRACE® qPCR Master Mix and de-ionized H₂O. Briefly vortex and centrifuge all tubes before opening.
- Label two 1.5 ml microcentrifuge tubes:
 - Sample 1 Mix
 - NTC Mix
- For each sample to be genotyped, prepare a qPCR Master Mix containing sample DNA, de-ionized H₂O and QTRACE® qPCR Master Mix as suggested by QTRACE® Software in Table 1. A No Template Control (NTC) Mix is prepared with de-ionized H₂O and QTRACE® qPCR Master Mix.

Table 1. qPCR Master Mix Composition

Sample 1 qPCR Mix	1 µl	50 µl
QTRACE® qPCR Master Mix	5.0 µl	250.0 µl
Annika Sample ID DNA	0.1 µl	5.0 µl
H ₂ O	19.9 µl	995.0 µl
	25.0 µl	1250.0 µl

NTC qPCR Mix	1 µl	50 µl
QTRACE® qPCR Master Mix	5.0 µl	15.0 µl
H ₂ O	20.0 µl	60.0 µl
	25.0 µl	75.0 µl

- Vortex each tube to thoroughly mix the contents and centrifuge briefly to collect the reaction mix at the bottom of the tube.
- Remove the adhesive cover from the genotyping plate.
- Dispense 25 µl of the Sample 1 Mix into Wells A1-D11 of the genotyping plate.
- Dispense 25 µl of the 5X PCR Master Mix/NTC mixture to well D12.
- A repeat pipettor is recommended to minimize pipetting repetition and increase accuracy.
- Refer to the DNA Sample Layout Plate Layout at the end of the protocol.
- Visually inspect plate wells from the sides and bottom to confirm consistent volume.

When data is imported, one half of the plate image will be colored.



After data Approval, TRACE Analysis™ will display informative assays for the combinations of samples associated with the recipient sample.

Marker	Chr. location	Informative for	Ulla	Jonna	Annika
240	8q	Ulla	Positive	Negative	Negative
356	18q	Ulla	Positive	Negative	Negative
469	7p	Jonna	Negative	Positive	Negative
504	17p	Ulla	Positive	Negative	Negative
520	20q	Jonna	Negative	Positive	Negative
531	17p	Jonna	Negative	Positive	Negative
634	11q	Ulla	Positive	Negative	Negative
736	Xp	Annika	Negative	Negative	Positive
768	18q	Ulla	Positive	Negative	Negative
819	14q	Jonna	Negative	Positive	Negative
854	12q	Annika	Negative	Negative	Positive
954	6q	Annika	Negative	Negative	Positive

Clicking the data overview



button shows the genotyping results for all assays.

Recipient		Type			
Recipient Name	Chr. location	Informative for	Ulla	Jonna	Annika
Ulla	1p	None	Negative	Negative	Negative
Ulla UPN	17q	None	Positive	Negative	Positive
Ulla Sample ID	5q	None	Positive	Positive	Positive
Concentration ng/ul	2q	None	Positive	Positive	Positive
Date of Birth	8q	Ulla	Positive	Negative	Negative
Date of Transplant	4q	None	Negative	Positive	Positive
Gender	5q	None	Positive	Positive	Negative
Comment	20q	None	Positive	Positive	Positive
Disease Type	18q	None	Negative	Negative	Negative
Donor	6q	None	Negative	Positive	Positive
Donor Name	1p	None	Negative	Negative	Negative
Donor ID	18q	Ulla	Positive	Negative	Negative
Donor Sample ID	11p	None	Negative	Positive	Positive
Concentration ng/ul	12q	None	Negative	Positive	Positive
Gender	18q	None	Positive	Positive	Positive
Comment	2q	None	Positive	Negative	Positive
Disease Type	6q	None	Positive	Negative	Positive
Donor	3q	None	Positive	Positive	Positive
Donor Name	7p	Jonna	Negative	Positive	Negative
Donor ID	17p	Ulla	Positive	Negative	Negative
Donor Sample ID	8q	None	Negative	Positive	Positive
Concentration ng/ul	20q	Jonna	Negative	Positive	Negative
Gender	17p	Jonna	Negative	Positive	Negative
Comment	Xp	None	Positive	Negative	Positive

Custom Genotyping Panels

With TRACE Analysis™ Software version 1.08, users of TRACE Analysis™ Software can perform customized genotyping. They can choose which assays they want to use for genotyping experiments. They are able to create, store and use custom genotyping assay panels to perform genotyping using TRACE Analysis™ Software. This additional functionality requires special license keys to be issued from JETA Molecular. Please contact your local sales representative or JETA Molecular directly for access.

TRACE Analysis™ will generate unique protocols, based on the experimental set up and it will score the data, as it does for standard QTRACE Genotyping Plates manufactured by JETA. Users who perform off-line genotyping and then assign informative assays in the software, will no longer have to do this.


Genotyping Panels are created using entire rows or columns of a 96-well plate.

By default, each panel must contain one positive control and one negative control reaction.

The possible number of INDEL assays in a genotyping panel are 6, 10, 14, 22, 30, 34, 38 and 46.


Orientation	Available INDEL Assay Panel Sizes	Plate Throughput
1 Row	10 Assays	8 Samples
2 Rows	22 Assays	4 Samples
3 Rows	34 Assays	2 Samples
4 Rows	46 Assays	2 Samples
1 Column	6 Assays	12 Samples
2 Columns	14 Assays	6 Samples
3 Columns	22 Assays	4 Samples
4 Columns	30 Assays	3 Samples
5 Columns	38 Assays	2 Samples
6 Columns	46 Assays	2 Samples

Under File, there is an option called 'Plate editor'




New plate configurations are created by pressing the 'Add new' button. Give the plate a name.

Choose the desired plate orientation and pick from the available assays to create a custom panel.



Once a Plate has been saved using the Plate editor, the plate appears as an option in the Plate Configuration menu in the Preferences.



The following is an example of a single donor:patient pair being tested by a 14 member panel with a column configuration:

Sample View:

	1	2	3	4	5	6	7	8	9	10	11	12
A	Myriam Pre-Tx 148	Myriam Pre-Tx 819	Vanessa Donor 148	Vanessa Donor 819								
B	Myriam Pre-Tx 240	Myriam Pre-Tx 832	Vanessa Donor 240	Vanessa Donor 832								
C	Myriam Pre-Tx 267	Myriam Pre-Tx 840	Vanessa Donor 267	Vanessa Donor 840								
D	Myriam Pre-Tx 345	Myriam Pre-Tx 916	Vanessa Donor 345	Vanessa Donor 916								
E	Myriam Pre-Tx 408	Myriam Pre-Tx 948	Vanessa Donor 408	Vanessa Donor 948								
F	Myriam Pre-Tx 425	Myriam Pre-Tx 954	Vanessa Donor 425	Vanessa Donor 954								
G	Myriam Pre-Tx 531	Myriam Pre-Tx RNaseF	Vanessa Donor 531	Vanessa Donor RNaseF								
H	Myriam Pre-Tx 736	NTC RNaseF	Vanessa Donor 736	NTC RNaseF								

Assay View:

	1	2	3	4	5	6	7	8	9	10	11	12
A	Myriam Pre-Tx 148	Myriam Pre-Tx 819	Vanessa Donor 148	Vanessa Donor 819								
B	Myriam Pre-Tx 240	Myriam Pre-Tx 832	Vanessa Donor 240	Vanessa Donor 832								
C	Myriam Pre-Tx 267	Myriam Pre-Tx 840	Vanessa Donor 267	Vanessa Donor 840								
D	Myriam Pre-Tx 345	Myriam Pre-Tx 916	Vanessa Donor 345	Vanessa Donor 916								
E	Myriam Pre-Tx 408	Myriam Pre-Tx 948	Vanessa Donor 408	Vanessa Donor 948								
F	Myriam Pre-Tx 425	Myriam Pre-Tx 954	Vanessa Donor 425	Vanessa Donor 954								
G	Myriam Pre-Tx 531	Myriam Pre-Tx RNaseF	Vanessa Donor 531	Vanessa Donor RNaseF								
H	Myriam Pre-Tx 736	NTC RNaseF	Vanessa Donor 736	NTC RNaseF								

Table 1. qPCR Master Mix Composition			
Sample 1 qPCR Mix		1 x	15 x
QTRACE® qPCR Master Mix		5.0 µl	75.0 µl
Pre-Tx DNA		0.5 µl	*7.5 µl
H ₂ O		19.5 µl	292.5 µl
		25.0 µl	375.0 µl
Sample 2 qPCR Mix		1 x	15 x
QTRACE® qPCR Master Mix		5.0 µl	75.0 µl
Donor DNA		0.5 µl	*7.5 µl
H ₂ O		19.5 µl	292.5 µl
		25.0 µl	375.0 µl
NTC qPCR Mix		1 x	2 x
QTRACE® qPCR Master Mix		5.0 µl	10.0 µl
H ₂ O		20.0 µl	40.0 µl
		25.0 µl	50.0 µl
7. Vortex each tube to thoroughly mix the contents and centrifuge briefly to collect the reaction mix at the bottom of the tube.			
8. Remove the adhesive cover from the genotyping plate.			
9. Dispense 25 µl of the Sample 1 Mix into Wells A1-G3 by columns of the genotyping plate.			
10. Dispense 25 µl of the Sample 2 Mix into Wells A3-G4 by columns of the genotyping plate.			
11. Dispense 25 µl of the SX PCR Master Mix/NTC mixture to wells H2, H4.			
12. A repeat pipettor is recommended to minimize pipetting repetition and increase accuracy.			
13. Refer to the DNA Sample Layout Plate Layout at the end of the protocol.			

The software can also combine multiple templates, if needed.

The following is an example of using 3 different column templates in one experiment. The first 2 samples added use 2 columns, the next 2 use 3 column and the last pair use 1 column.

The protocol instructs the user to make 7 different mixes for typing. The mixes vary on size, depending on how many assays are in the panel

Table 1. qPCR Master Mix Composition

Sample 1 qPCR Mix	1 x	15 x	Sample 5 qPCR Mix	1 x	7 x
QTRACE® qPCR Master Mix	5.0 µl	75.0 µl	QTRACE® qPCR Master Mix	5.0 µl	35.0 µl
180425R6-3 DNA	0.5 µl	*7.5 µl	180425R6-5 DNA	0.5 µl	*3.5 µl
H ₂ O	19.5 µl	292.5 µl	H ₂ O	19.5 µl	136.5 µl
	25.0 µl	375.0 µl		25.0 µl	175.0 µl

Sample 2 qPCR Mix	1 x	15 x	Sample 6 qPCR Mix	1 x	7 x
QTRACE® qPCR Master Mix	5.0 µl	75.0 µl	QTRACE® qPCR Master Mix	5.0 µl	35.0 µl
180425D6-3 DNA	0.5 µl	*7.5 µl	180425D6-5 DNA	0.5 µl	*3.5 µl
H ₂ O	19.5 µl	292.5 µl	H ₂ O	19.5 µl	136.5 µl
	25.0 µl	375.0 µl		25.0 µl	175.0 µl

Sample 3 qPCR Mix	1 x	23 x	NTC qPCR Mix	1 x	6 x
QTRACE® qPCR Master Mix	5.0 µl	115.0 µl	QTRACE® qPCR Master Mix	5.0 µl	30.0 µl
180425R6-4 DNA	0.5 µl	*11.5 µl	H ₂ O	20.0 µl	120.0 µl
H ₂ O	19.5 µl	448.5 µl		25.0 µl	150.0 µl
	25.0 µl	575.0 µl			

and the pipetting instruction respond to the templates being used:

SECTION OF THE TABLE

- Remove the adhesive cover from the genotyping plate.
- Dispense 25 µl of the Sample 1 Mix into Wells A1-G2 by columns of the genotyping plate.
- Dispense 25 µl of the Sample 2 Mix into Wells A3-G4 by columns of the genotyping plate.
- Dispense 25 µl of the Sample 3 Mix into Wells A5-G7 by columns of the genotyping plate.
- Dispense 25 µl of the Sample 4 Mix into Wells A8-G10 by columns of the genotyping plate.
- Dispense 25 µl of the Sample 5 Mix into Wells A11-G11 by columns of the genotyping plate.
- Dispense 25 µl of the Sample 6 Mix into Wells A12-G12 by columns of the genotyping plate.
- Dispense 25 µl of the 5X PCR Master Mix/NTC mixture to wells H2, H4, H7, H10, H11, H12.
- A repeat pipettor is recommended to minimize pipetting repetition and increase accuracy.**
- Refer to the DNA Sample Layout Plate Layout at the end of the protocol.
- Visual inspection of plate wells from the sides and bottom to confirm consistent volume.

Here is the layout:

	1	2	3	4	5	6	7	8	9	10	11	12
A	180425R6-3 291	180425R6-3 519	180425D6-3 291	180425D6-3 519	180425R6-4 337	180425R6-4 408	180425R6-4 710	180425D6-4 337	180425D6-4 408	180425D6-4 710	180425R6-5 337	180425D6-5 337
B	180425R6-3 305	180425R6-3 634	180425D6-3 305	180425D6-3 634	180425R6-4 148	180425R6-4 519	180425R6-4 721	180425D6-4 148	180425D6-4 519	180425D6-4 721	180425R6-5 148	180425D6-5 148
C	180425R6-3 312	180425R6-3 824	180425D6-3 312	180425D6-3 824	180425R6-4 209	180425R6-4 520	180425R6-4 736	180425D6-4 209	180425D6-4 520	180425D6-4 736	180425R6-5 209	180425D6-5 209
D	180425R6-3 321	180425R6-3 832	180425D6-3 321	180425D6-3 832	180425R6-4 345	180425R6-4 615	180425R6-4 748	180425D6-4 345	180425D6-4 615	180425D6-4 748	180425R6-5 345	180425D6-5 345
E	180425R6-3 326	180425R6-3 832	180425D6-3 326	180425D6-3 832	180425R6-4 358	180425R6-4 626	180425R6-4 840	180425D6-4 358	180425D6-4 626	180425D6-4 840	180425R6-5 358	180425D6-5 358
F	180425R6-3 434	180425R6-3 907	180425D6-3 434	180425D6-3 907	180425R6-4 359	180425R6-4 626	180425R6-4 840	180425D6-4 359	180425D6-4 626	180425D6-4 840	180425R6-5 359	180425D6-5 359
G	180425R6-3 469	180425R6-3 RNaSP	180425D6-3 469	180425D6-3 RNaSP	180425R6-4 373	180425R6-4 634	180425R6-4 RNaSP	180425D6-4 373	180425D6-4 634	180425D6-4 RNaSP	180425R6-5 373	180425D6-5 RNaSP
H	180425R6-3 504	180425R6-3 RNaSP	180425D6-3 504	180425D6-3 RNaSP	180425R6-4 386	180425R6-4 706	180425R6-4 RNaSP	180425D6-4 386	180425D6-4 706	180425D6-4 RNaSP	180425R6-5 386	180425D6-5 RNaSP

Virtual Typing



Virtual Typing functionality was added in TRACE Analysis™ version 1.5. Virtual Typing allows the user to compare any DNA genotyped against any other using the same genotyping plate configuration. This functionality is helpful for validation of new lots of plates, as well as for research purposes. Reports can also be generated from any custom view created in virtual typing.

Validation of a new lot of plates may require verification of the typing patterns within a plate as well as with multiple known DNA samples. Inspecting 46 assays many times over against a spreadsheet or printed list of known genotypes is extremely tedious and prone to error. With Virtual Typing functionality it is easy to quickly see whether a new lot of genotyping plates gives the expected results.

Similarly, if users need to make control mixtures, the VirtualTyping function allows easy comparison of controls which can be made for specific assays of interest.

Assay	Chromosome location	Informative for ...	443 G7 DON 8889	536 D2 REF 0013
137	1p	-	Negative	Negative
148	17q	-	Positive	Positive
209	5q	G7 DON	Informative	Negative
235	2q	G7 DON	Informative	Negative
240	8q	G7 DON	Informative	Negative
267	4q	D2 REF	Negative	Informative
312	18q	-	Negative	Negative
326	6q	G7 DON	Informative	Negative
345	1p	D2 REF	Negative	Informative
356	18q	D2 REF	Negative	Informative
359	11p	-	Positive	Positive
361	12q	-	Negative	Negative
373	18q	-	Positive	Positive
386	5q	-	Positive	Positive
408	2q	D2 REF	Negative	Informative
425	6q	G7 DON	Informative	Negative
434	3q	D2 REF	Negative	Informative
469	7p	-	Negative	Negative
504	17p	D2 REF	Negative	Informative
519	8q	-	Negative	Negative
520	20q	D2 REF	Negative	Informative
531	17p	D2 REF	Negative	Informative
548	Xq	-	Negative	Negative
601	14q	-	Negative	Negative
615	1q	-	Positive	Positive

Virtual Typing Report

Sample 1

Name: G7 DON

ID: 443

Gender: Unknown

Experiment date: 08 August 2024

Experiment name: D2 G7 Genotyping

Comments: -

Sample 2

Name: D2 REF

ID: 536

Gender: Unknown

Experiment date: 08 August 2024

Experiment name: D2 G7 Genotyping

Comments: -

Assay	Chromosome location	Informative for	Sample 1 G7 DON 443	Sample 2 D2 REF 536
137	1p	-	Neg	Neg
148	17q	-	Pos	Pos
209	5q	G7 DON	Info	Neg
235	2q	G7 DON	Info	Neg
240	8q	G7 DON	Info	Neg
267	4q	D2 REF	Neg	Info
312	18q	-	Neg	Neg
326	6q	G7 DON	Info	Neg
345	1p	D2 REF	Neg	Info
356	18q	D2 REF	Neg	Info
359	11p	-	Pos	Pos
361	12q	-	Neg	Neg
373	18q	-	Pos	Pos
386	5q	-	Pos	Pos
408	2q	D2 REF	Neg	Info
425	6q	G7 DON	Info	Neg
434	3q	D2 REF	Neg	Info
469	7p	-	Neg	Neg
504	17p	D2 REF	Neg	Info
519	8q	-	Neg	Neg
520	20q	D2 REF	Neg	Info
531	17p	D2 REF	Neg	Info
548	Xq	-	Neg	Neg
601	14q	-	Neg	Neg
615	1q	-	Pos	Pos
624	11q	G7 DON	Info	Neg
625	7q	-	Neg	Neg

Monitoring Test

Once recipient specific markers have been found, quantification is performed to monitor engraftment. In the monitoring test, two or more of the informative assays is used to quantify the DNA of interest in an unknown sample relative to a reference sample. Any of the informative assays identified in the genotyping can be used to perform monitoring.

The QTRACE® Analysis System

The fraction of DNA positive for the informative marker in the unknown composition is determined relative to the pure reference sample DNA and the result is expressed as a percentage (ratio).

Monitoring Sample View

Monitoring Samples are arranged into groups byTRACE Analysis™ Software. The NTC reactions for each assay are in the lower right hand corner of the plate.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Recipe Recipe RNaseF	Recipe Recipe RNaseF	Recipe Recipe RNaseF	Recipe 1st Mo RNaseF	Recipe 1st Mo RNaseF	Recipe 1st Mo RNaseF						
B	Recipe Recipe 137	Recipe Recipe 137	Recipe Recipe 137	Recipe 1st Mo 137	Recipe 1st Mo 137	Recipe 1st Mo 137						
C	Recipe Recipe 148	Recipe Recipe 148	Recipe Recipe 148	Recipe 1st Mo 148	Recipe 1st Mo 148	Recipe 1st Mo 148						
D												
E												
F												
G												
H										NTC 148	NTC 137	NTC RNaseF

Monitoring Assay View

TRACE Analysis™ Software adds the assays to the plate in rows.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Recipe Recipe RNaseF	Recipe Recipe RNaseF	Recipe Recipe RNaseF	Recipe 1st Mo RNaseF	Recipe 1st Mo RNaseF	Recipe 1st Mo RNaseF						
B	Recipe Recipe 137	Recipe Recipe 137	Recipe Recipe 137	Recipe 1st Mo 137	Recipe 1st Mo 137	Recipe 1st Mo 137						
C	Recipe Recipe 148	Recipe Recipe 148	Recipe Recipe 148	Recipe 1st Mo 148	Recipe 1st Mo 148	Recipe 1st Mo 148						
D												
E												
F												
G												
H										NTC 148	NTC 137	NTC RNaseF

The DigitalTRACE™ Analysis System

The fraction of DNA positive for the informative marker in the unknown composition is determined relative to the reference gene (RNase P) signal and the result is expressed as a percentage (ratio).

Monitoring Sample View

Monitoring Samples are arranged into groups byTRACE Analysis™ Software. The NTC and UPC (Universal Positive Control) reactions for each assay are in the lower right hand corner of the plate.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Recipe Post 1 d137											
B	Recipe Post 1 d148											
C	Recipe Post 2 d137											
D	Recipe Post 2 d148											
E											UPC d137	
F											UPC d148	
G											NTC d137	
H											NTC d148	

Monitoring Assay View

TRACE Analysis™ Software adds the assays to the plate in columns.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Recipe Post 1 d137											
B	Recipe Post 1 d148											
C	Recipe Post 2 d137											
D	Recipe Post 2 d148											
E											UPC d137	
F											UPC d148	
G											NTC d137	
H											NTC d148	

See section [Drag and Drop Plate Layouts](#) to move a well from one location to another location.


Assigning Informative Markers

TRACE Analysis™ Software allows you to perform sample monitoring and leverage all the features of TRACE Analysis™, without the need to genotype the sample in advance. If you know that an assay is informative for your sample, you can designate its informative status and then use it as you normally would.


In order to use this feature of TRACE Analysis™:

Enter all the information necessary about the Recipient and the Donor samples and choose an appropriate Reference Sample.

Enter information related to the Post Sample(s) being tested.

In the Informative Markers section, click the "Add Marker"  button.


When this button is clicked, a window opens with all Assays available for assignment.




Assign the informative assays by clicking the circle and sample combination which is appropriate.

Press the 'X' to reset the choice for that assay.


In case of DigitalTRACE™ Assays, the assay name begins with the "d" symbol.



In the HLA tab, you can select for HLA TRACE™ Assays designed for HLA Loss of Heterozygosity monitoring.



In the Custom tab, custom assays can be added.



Once the marker(s) have been assigned, click ok.

Now the markers are selectable for the sample in the Informative Markers window.

TRACE Analysis™ Software will save this information, so it only has to be entered once for a given sample.

Protocol

In order to perform a monitoring experiment with TRACE Analysis™ Software, start typing the name of your Recipient Sample in the Sample Entry window. TRACE Analysis™ Software will start searching for the record as you are entering the name. The name can also be chosen from the drop down menu which appears when you click in the sample entry box.

The QTRACE® Analysis System

Once the Recipient Sample name is found and selected, the identifying information initially entered for the samples and the genotyping data results will appear.

(For samples wherein genotyping data was not determined using TRACE Analysis™ Software, please see section on [Assigning Informative Markers](#)).

The screenshot shows the Q-TRACE software interface. The main window is titled 'Q-TRACE™' and has a menu bar with 'File', 'View', and 'Help'. Below the menu bar is a home icon and a search bar. The interface is divided into several sections:

- Recipient**: Fields for Recipient Name (Recipient), Recipient ID (Recipient ID), Sample ID (Recipient), and Concentration ng/ul (75).
- Donor**: Fields for Donor Name (Donor), Donor ID (Donor ID), Sample ID (Donor Sample), and Concentration ng/ul (50). There is an 'Add Donor' button with a plus sign.
- Reference Samples**: A table with columns 'Sample ID' and 'Concentration ng/ul'. It contains rows for Recipient Sample (75), Donor Sample (50), Cell Line 1 (10), Cell Line 2 (10), and Cell Line 3 (10).
- Post Sample**: Fields for Sample ID, Concentration ng/ul (20), and Sample Date (XX-XX-XXXX).

The Reference Sample(s) and its/their concentration(s) for a group of samples in the experiment must be defined. This is typically a pure sample of the recipient's DNA.

In the "Post Sample" entry window, a unique Sample ID must be entered in the Sample ID field.

A sample Date must be entered.

Check the Sample type which is being tested - Sample Types can be defined by the User in the

[Preferences](#)  menu.

Click on the name(s) of the Assays you wish to use for the experiment.

Once all the information about the sample has been entered, click the “Add Quant Samples”

Quant

button and the sample will be added to the plate.

The reactions for a particular recipient sample are grouped together. In the example below, all reactions are performed in triplicate with two informative markers.

You can toggle between coloring for the Sample View and the Assay View by pressing the appropriate small plate at the bottom of the window.

One View (see "Sample Layout view") shows the placement of Pure (pre-transfer) and Mixed Samples (Post-Transfer). Pre samples go in wells A1-A3, B1-B3 and C1-C3. Post samples go in A4-A6, B4-B6 and C4-C6.

Sample Layout view

	1	2	3	4	5	6	7	8	9	10	11	12
A	Recipe Recipe RNasef	Recipe Recipe RNasef	Recipe Recipe RNasef	Recipe 1st Mo RNasef	Recipe 1st Mo RNasef	Recipe 1st Mo RNasef						
B	Recipe Recipe 137	Recipe Recipe 137	Recipe Recipe 137	Recipe 1st Mo 137	Recipe 1st Mo 137	Recipe 1st Mo 137						
C	Recipe Recipe 148	Recipe Recipe 148	Recipe Recipe 148	Recipe 1st Mo 148	Recipe 1st Mo 148	Recipe 1st Mo 148						
D												
E												
F												
G												
H										NTC 148	NTC 137	NTC RNasef


The other view shows the placement of the Assays. The Reference Gene Assay is used in A1-A6 and H12. Assay 014 is used in B1-B6 and H11. Assay 016 is used in C1-C6 and H10

Assay Layout view

	1	2	3	4	5	6	7	8	9	10	11	12
A	Recipe Recipe RNasef	Recipe Recipe RNasef	Recipe Recipe RNasef	Recipe 1st Mo RNasef	Recipe 1st Mo RNasef	Recipe 1st Mo RNasef						
B	Recipe Recipe 137	Recipe Recipe 137	Recipe Recipe 137	Recipe 1st Mo 137	Recipe 1st Mo 137	Recipe 1st Mo 137						
C	Recipe Recipe 148	Recipe Recipe 148	Recipe Recipe 148	Recipe 1st Mo 148	Recipe 1st Mo 148	Recipe 1st Mo 148						
D												
E												
F												
G												
H										NTC 148	NTC 137	NTC RNasef

The reactions in Row H constitute No Template Controls (NTCs) and have deionized water added in place of the sample.

Proceed to add more experiments to fill the plate as per you needs.

Once all samples have been added to your plate, press the “Export Setup to PCR”  button. Browse to the location where you want the file saved and name it as you wish. This file can then be imported into your qPCR machine’s QTRACE® template to execute the qPCR analysis.

Once the file is saved, TRACE Analysis™ generates a protocol, based on the experimental inputs and the settings in the preferences menu.

The DigitalTRACE™ Analysis System

For Post samples wherein genotyping data was either determined using qPCR or not determined using TRACE Analysis™ Software, please see section on [Assigning Informative Markers](#).

In order to perform a Monitoring experiment with TRACE Analysis™ Software using both the QIAcuity or Biorad dPCR instrument, select the name of your Recipient Sample in the Sample Entry window from the drop-down menu.

Once the Recipient Sample name is found and selected, choose from a pop-up window to use monitoring plate format which must be selected before in the Preferences. The selected sample's identifying information initially entered will appear.

Type in the Post sample information:

- A unique Sample ID must be entered in the Sample ID field.

- A sample Date must be entered.

- Check the Sample type which is being tested - Sample Types can be defined by the User in the Preferences menu.

- Click the plus sign next to the "Add Sample" tab

In the Informative Assays window, choose the loci to be tested by selecting assays from the list based on the screening test results.

Once all the information about the sample has been entered, click the "Add Quant Samples"




button and the sample will be added to the plate.

A negative control (NTC) is automatically added to the plate layout. A positive control (UPC) is added if it is enabled in the dPCR Data and Reports tab in the [Preferences](#). The use of UPC is optional.

You can toggle between coloring for the Sample View and the Assay View by pressing the appropriate small plate at the bottom of the window. One View shows the placement of monitoring and control samples. The other view shows the placement of the Assays.

Monitoring Samples are arranged into groups by TRACE Analysis™ Software, based on the preferences set for the monitoring machine. You can drag the wells to rearrange the samples within the plate. In this case, both controls were placed to the first column.


	1	2	3	4	5	6	7	8	9	10	11	12
A	Basti B Post20 d252											
B	Basti B Post20 d472											
C	Basti B Post20 d252											
D	Basti B Post20 d472											
E	UPC d472											
F	UPC d252											
G	NTC d472											
H	NTC d252											



Sample view


Experiment name:

Operator ID:



Assay view

Proceed to add more samples to fill the plate as per your needs.

Once all samples have been added to your plate, press the “Export Setup to PCR”  button. Browse to the location where you want the file saved and name it as you wish. This file can then be imported into your dPCR machine’s template file to execute the dPCR analysis.

Once the file is saved, TRACE Analysis™ Software generates a protocol, based on the experimental inputs and the settings in the preferences menu.

Monitoring Test Protocol - QIAcuity

The following protocol is an example output from TRACE Analysis™ Software for a Monitoring test using QIAcuity dPCR:

- Set up all reactions in a pre-PCR lab, under ambient conditions without ice.
- Briefly vortex and centrifuge all tubes before opening.
- Prepare DNA dilutions as specified in Table 1:

No.	Name	Sample ID	Concentration	Sample Volume	Water Volume	Total Volume	Amount per Reaction
1	Post Sample	Post_test	26 ng/μl	12,69 μl	48,91 μl	61,60 μl	150 ng
2	Post Sample	Post2_test	41 ng/μl	8,05 μl	53,55 μl	61,60 μl	150 ng
- Prepare for chosen informative dPCR Assay the following mixture:


Assay Mix	QIAcuity 4x Master Mix	dPCR 20x assay mix [μl]	Total Volume
d359	46,00 μl	9,20 μl	55,20 μl
d626	46,00 μl	9,20 μl	55,20 μl
- IMPORTANT! Vortex at least 5 seconds and spin briefly each prepared Master Mix + dPCR Assay mixture.**
- Deliver 13,2 μl of QIAcuity Probe Master Mix + dPCR Assay Mix to a 96-well plate as defined in TRACE Analysis™ Software's Assay Layout view.
- Add 30,8 μl of each Sample DNA dilution as indicated by TRACE Analysis™ Software's Layout View. Where NTC and UPC are indicated add 30,8 μl water for NTC and 30,8 μl from Universal Positive Control for UPC.
- IMPORTANT! Vortex at least 15 seconds the 96-well plate to mix the contents of each reaction. Centrifuge the 96-well plate briefly using a microcentrifuge.**
- Transfer 40 μl of each prepared reaction mix into a Nanoplate. Seal the Nanoplate with the compatible plate sealer.
- Load the Nanoplate into the QIAcuity digital PCR system.
- Launch the QIAcuity Software Suite.
- Open the DigitalTRACE template and import the Sample Setup sheet generated by TRACE Analysis™ Software.
- Save the file and start the run.

Experiment Setup in QIAcuity Suite Software, v2.5

Create a new QIAcuity Plate by selecting "New Plate".

Similarly as for the Genotyping, load a template by selecting "Plate templates..." and import the appropriate DigitalTRACE Monitoring template.


Press Save Plate. The new plate will appear in the main window of the QIAcuity Software Suite.



Click on the plate name to open the plate configuration procedure. Type in a new plate name and save the changes.

Import the Sample Setup sheet (.csv) generated by TRACE Analysis™ Software by selecting "Plate layout" tab and "CSV import/export".

Select "Import from CSV".



Inspect if all your selected samples are represented on the Plate figure. Select "Done" to finish the setup. Start the run.

Thermal Cycling and Imaging Protocol for QIAcuity

The DigitalTRACE™ System will deliver optimal results when the following thermal profile is used in dPCR:

Number of cycles	Temperature	Time
1	95 °C	3 min
40	95 °C	15 s
	60 °C	30 s

Channel	Exposure duration	Gain
Green	500 ms	6
Yellow	500 ms	6
Orange	Off	Off
Red	Off	Off
Crimson	Off	Off

After the QIAcuity run has finished, check in the QIAcuity Software Suite if the automatic thresholds are correct. Adjust the thresholds manually if needed.
Select all wells on the plate and analyze per target (not per channel), export data by selecting Export to CSV.

Monitoring Test Protocol - Biorad

The following protocol is an example output from TRACE Analysis™ Software for a Monitoring test using Biorad ddPCR:

1. Set up all reactions in a pre-PCR lab, under ambient conditions without ice.
2. Briefly vortex and centrifuge all tubes before opening.
3. Prepare DNA dilutions as specified in Table 1:

No.	Name	Sample ID	Concentration	Sample Volume	Water	Total Volume	Amount per Reaction
1	Basti Bob	Post201	20 ng/μl	17,25 μl	5,52 μl	22,77 μl	150 ng
2	Basti Bob	Post202	20 ng/μl	17,25 μl	5,52 μl	22,77 μl	150 ng

4. Prepare for chosen informative dPCR Assay the following mixture:

Assay	ddPCR supermix for probes (no dUTPs) [μl]	dPCR assay [μl]	Total Volume
d252	46,00 μl	4,60 μl	50,60 μl
d472	46,00 μl	4,60 μl	50,60 μl

5. Deliver 12.1 μl of ddPCR Supermix + dPCR Assay mixture to a 8 strip PCR tube as defined in TRACE Analysis™ Software's Assay Layout view.
6. Add 9.9 μl of each Sample DNA dilution and water for NTC wells as indicated by TRACE Analysis™ Software's Layout View.
7. **IMPORTANT!** Vortex the 8 strip tube to mix the contents of each reaction. Centrifuge the 8 strip tubes briefly using a microcentrifuge.
8. **!!!The following steps from 9 to 14 are only for manual droplet generator users:**
9. Transfer 20 μl of each prepared sample to the sample wells (middle row) of the DG8 cartridge.
10. Add 70 μl of droplet generation oil to each oil well of the DG8 cartridge.
11. Hook the gasket over the cartridge holder using the holes on both sides.
12. Load the cartridge in the QX200 droplet generator.
13. When droplet generation is complete, remove the disposable gasket from the holder and discard it.
14. Pipet 40 μl of the contents of the droplets into a single column of a 96-well PCR plate.
15. Seal the PCR plate with foil plate seals that are compatible with the PX1 PCR plate sealer and the needles in the QX200 droplet reader.
16. Place the plate into the thermal cycler for PCR amplification.
17. Load the plate after amplification into QX200 droplet reader.
18. Import the Sample Setup sheet generated by TRACE Analysis™ Software.
19. Save the file and start the droplet reader.


Use the same cycling parameters as for a Genotyping test:

Number of cycles	Temperature	Time	Ramp Rate
1	95 °C	10 min	2 °C/s
40	94 °C	30 s	2 °C/s
	59 °C	60 s	2 °C/s
1	98 °C	10 min	2 °C/s

The droplet reading process is the same as for the Genotyping test.
Start the droplet reading run.

After the reading process has finished, check in the QuantaSoft Software if the automatic thresholds are correct. Adjust the thresholds manually if needed.
Select all wells on the plate and export data by selecting Export CSV.

Monitoring Data Analysis and Report

Import qPCR data by clicking the “Import PCR Data”  button, and browse to the location of your exported qPCR data file and select it.

The QTRACE® Analysis System


TRACE Analysis™ Software uses a green/amber/red color coding for well highlighting.




Green wells always mean positive wells (Go). Red wells are data which has been rejected, manually or automatically (Stop). Amber wells highlight the need for data inspection (Caution).

Light blue wells  remain as negative wells.


Monitoring plates contain green and light blue colors, indicating true positives or negatives. Amber colored wells highlight replicates which violate user-defined replicate highlighting rules, based on %CV, SD or Cq Range.




Replicate wells in Monitoring which violate the replicate highlighting preferences set by the user will first be highlighted in Amber.




The user can then Accept/Reject individual wells in an attempt to conform to the acceptable variance allowed. Rejected wells turn Red and the remaining wells will either turn green, light blue or remain Amber - if the Rejected well does not cause the remaining replicates to meet acceptable variance.




Once all the data has been inspected, pressing the “Calculate”  button allows TRACE Analysis™ to perform the calculations necessary to determine the proportion of the mixed sample which contains the sample of interest.


The percentage of the total sample which contains the sample is displayed in the right hand Result window.



As more data is collected for a particular sample over time, TRACE Analysis™ Software provides a longitudinal view.


To view the composite set of data for an individual sample, press the “Overview”  button.

Pushing the Overview button takes you to a screen showing all the information input for a particular group of samples as well as all of their genotyping and monitoring data.



The DigitalTRACE™ Analysis System


TRACE Analysis™ Software uses a light green/dark green/red color coding for well highlighting.




Light green wells represent post sample wells. Dark green wells represent the positive control DNA, Red wells are data which has been rejected.

Light blue wells represent non-template control (NTC).

The user can then Accept/Reject individual wells in an attempt to conform to the acceptable variance allowed. Rejected wells turn Red.



When necessary, for example in case that no CNV information is available from the genotyping, TRACE Analysis™ Software allows the user to change the CNV information for each sample by selecting 1 or 2 in the indicated area under the Rawdata view **Show rawdata** in order to obtain an accurate quantification result:



Once all the data has been inspected, pressing the “Calculate” **Calculate** button allows TRACE Analysis™ to perform the calculations necessary to determine the proportion of the mixed sample which contains the sample of interest.

The percentage of the total sample which contains the sample is displayed in the right hand Result window.

Sample	Type	Date	Target	Chr.	CNV	Ratio	DNA (%)	DNA (%) + CNV
276056	DNA	01-05-2023	d267	4q	1.00	0.15600	31.20	31.20
279413	DNA	01-05-2023	d267	4q	1.00	0.16700	33.40	33.40
30p108	DNA	01-05-2023	d267	4q	1.00	0.13750	27.50	27.50


Click on **Select assays** and then on **Report...** to generate the Monitoring Report.

As more data is collected for a particular sample over time, TRACE Analysis™ Software provides a longitudinal view. For more details, see the [The QTRACE® Analysis System](#) part of this section.

Assay Filtering from Monitoring Reports

Assay Filtering is a feature which allows a user to selectively remove all data from a particular assay in the final report. This may be important, for example, if one pipette tip from a multi-channel pipette did not properly dispense the necessary reagents to a set of reactions. While the data may appear to have good precision - all negatives in this example - it may not be accurate.

In the following image, there is a failure of Assay 721 for all samples:



If the plate is approved, under the Results pane, there is 0% reported for Assay 721, while the other assays detect chimerism levels well above zero.

Result					
Date	Target	Chr.	ΔΔCq	DNA (%)	Reference
17-08-2018	434	3q	1.67	31.32	R Control
17-08-2018	721	Xq	0.00	0.00	R Control
17-08-2018	824	22q	1.70	30.76	R Control
17-08-2018	434	3q	10.65	0.06	R Control
17-08-2018	721	Xq	0.00	0.00	R Control
17-08-2018	824	22q	9.85	0.11	R Control
17-08-2018	434	3q	7.39	0.60	R Control
17-08-2018	721	Xq	0.00	0.00	R Control
17-08-2018	824	22q	7.34	0.62	R Control
17-08-2018	434	3q	3.80	7.16	R Control
17-08-2018	721	Xq	0.00	0.00	R Control
17-08-2018	824	22q	3.96	6.45	R Control
17-08-2018	434	3q	8.80	0.22	R Control
17-08-2018	721	Xq	0.00	0.00	R Control
17-08-2018	824	22q	9.09	0.18	R Control


If a report is made at this point, the data from Assay 721 will be included in the calculations and the result will underestimate the actual level of recipient in the same.

Cell fractions

Sample	Type	Date	Target	Chr.	ΔΔCq	DNA (%)	Reference	Informative for
R Sample 1	BM	17-Aug-2018	434	3q	1.67	31.32	R Control	R Control
R Sample 1	BM	17-Aug-2018	721	Xq	0	0	R Control	R Control
R Sample 1	BM	17-Aug-2018	824	22q	1.7	30.76	R Control	R Control
					Mean:	20.69		
R Sample 2	BM	17-Aug-2018	434	3q	10.65	0.06	R Control	R Control
R Sample 2	BM	17-Aug-2018	721	Xq	0	0	R Control	R Control
R Sample 2	BM	17-Aug-2018	824	22q	9.85	0.11	R Control	R Control
					Mean:	0.06		
R Sample 3	BM	17-Aug-2018	434	3q	7.39	0.6	R Control	R Control
R Sample 3	BM	17-Aug-2018	721	Xq	0	0	R Control	R Control
R Sample 3	BM	17-Aug-2018	824	22q	7.34	0.62	R Control	R Control
					Mean:	0.41		

In order to remove the data for Assay 721 from the Report, after approving the data, press the 'Select Assays' **Select assays** button.

The 'FilterAssaysWindow' will appear. In this window, users can choose to exclude the data from an entire assay from reporting.



In this example, once Assay 721 is de-selected and the changes saved, once the 'Report' button is pressed, the data is no longer present in the report and the average chimerism levels reflect only the included data.

Cell fractions

Sample	Type	Date	Target	Chr.	$\Delta\Delta C_q$	DNA (%)	Reference	Informative for
R Sample 1	BM	17-Aug-2018	434	3q	1.67	31.32	R Control	R Control
R Sample 1	BM	17-Aug-2018	824	22q	1.7	30.76	R Control	R Control
Mean:						31.04		
R Sample 2	BM	17-Aug-2018	434	3q	10.65	0.06	R Control	R Control
R Sample 2	BM	17-Aug-2018	824	22q	9.85	0.11	R Control	R Control
Mean:						0.09		
R Sample 3	BM	17-Aug-2018	434	3q	7.39	0.6	R Control	R Control
R Sample 3	BM	17-Aug-2018	824	22q	7.34	0.62	R Control	R Control
Mean:						0.61		

External Reference Functions

The External Reference sample within TRACE Analysis™ Software has been enhanced in version 2.0. In addition to being able to specify and use an external reference sample in a qPCR experiment, now the sample may be an impure sample (e.g., an earlier timepoint from the recipient) or a Universal Positive Control (Ref. number 711295).

Typically, the reference sample in qPCR is pure donor or pre-transplant, recipient DNA. In both cases, the reference sample is considered 100% and the calculation of % DNA, based on ddCq, assumes this.

A sample which is less than 100% recipient may be used as a reference, but the software needs to correct the % DNA calculation, based on the starting % of the sample.

TRACE Analysis™ Software now provides this flexibility and automatic correction in monitoring data.

Reference Samples		
Sample ID	Concentration ng/ul	% reference
	100	100
	100	100
10% Recipient	10	10

Here is an example of processing the same data wherein the reference sample is either a mixture or a pure sample.

Cell fractions

Sample	Type	Date	Target	Chr.	$\Delta\Delta Cq$	DNA (%)
15	DNA	21-Jan-2020	548	Xq	2.88	1.36
15	DNA	21-Jan-2020	361	12q	2.72	1.52
15	DNA	21-Jan-2020	916	10q	2.73	1.5
Mean:						1.46
15	DNA	21-Jan-2020	916	10q	2.73	15.04
15	DNA	21-Jan-2020	361	12q	2.72	15.16
15	DNA	21-Jan-2020	548	Xq	2.88	13.6
Mean:						14.6

10% Reference Sample

Pure Reference Sample

Simultaneous Genotyping and Monitoring

With the use of Custom Genotyping Panels, TRACE Analysis™ Software has the potential to allow concurrent genotyping and monitoring experiments on a single 96-well plate.

This example shows the use of a genotyping panel of 10 assays for a new recipient and also monitoring a different recipient on the same plate.

After setting up a custom panel based on one Row in the Plate editor, choose the name of the custom panel in the 'Plate configuration' drop down menu (Plate setup tab of Preferences).

Preferences

Plate setup Custom types Concentrations Data and Reports Data locations Material tracking Users Language

Machine for genotyping ABI 7500 v2.0.6

Machine for monitoring ABI 7500 v2.0.6

Layout type

Replicates qPCR Duplicates

Color mode Pastel

Plate for genotyping 1 Row

Plate for monitoring 96-low density

Ok Cancel Apply

Set up a genotyping test, followed by a monitoring test. When the 'Export' button is pressed, a unique, combined protocol is generated by TRACE Analysis™ Software. The combined protocol is consistent with the normal workflow for a monitoring experiment. That is, all reagent additions are performed manually, with master mix, assays and DNA dilutions being added to the plate.

Technical Support
 JETA Molecular
 Krommewetering 101C
 3543AN, Utrecht
 +31.654.136.697


GENOTYPING PROTOCOL FOR NEW TYPING R
QUANTIFICATION PROTOCOL FOR RETURNING PATIENT

1. Set up all reactions in a pre-PCR lab, under ambient conditions without ice.
2. Briefly vortex and centrifuge all tubes before opening.
3. Prepare DNA dilutions as specified in Table 1:
 * - Make a Ten-Fold Dilution (1:10) of Sample

Table 1. DNA dilutions

No.	Name	Sample ID	Concentration	Sample Volume	Water Volume	Total Volume	Amount per Reaction
1	New Typing R	10831 13 Aug	*10 ng/μl	6.05 μl	175.45 μl	181.50 μl	5 ng
2	New Typing D	12560 13 Aug	*10 ng/μl	6.05 μl	175.45 μl	181.50 μl	5 ng
3	Returning Patient	DNA 103 13 Aug	*10 ng/μl	13.20 μl	184.80 μl	198.00 μl	10 ng
4	Returning Patient	Post 1 13 Aug	100 ng/μl	19.80 μl	178.20 μl	198.00 μl	150 ng
5	Returning Patient	Post 2 13 Aug	100 ng/μl	19.80 μl	178.20 μl	198.00 μl	150 ng
6	Returning Patient	Post 3 13 Aug	100 ng/μl	19.80 μl	178.20 μl	198.00 μl	150 ng

After the results from the qPCR machine are imported into TRACE Analysis™ Software, the analysis algorithms appropriate for the samples are applied and the results displayed in the Result pane of the software.



TRACE Analysis™ Software also generates a report for both recipients being tested on the single plate.

Software Overview

Use of TRACE Analysis™ Software facilitates setting up qPCR and dPCR based tests for multiple platforms, analyzes data, calculates and displays analyzed data and stores sample-specific information for easy retrieval or exporting to a laboratory information management system. All data files generated by the software are stored in .xml format, for ease of data transfer.



Preferences

The Preferences Menu of TRACE Analysis™ Software allows for customization of many experimental parameters.

Through the Preferences window, you can choose default settings for common variables, as well as enter information which may be unique to their testing regime or laboratory.

There are four main tabs within the window which enable changes to the experimental plate setup, manual entry of sample types, disease states and reference samples, experimental protocol variables and data storage customization. These customized parameters are found on the Plate Setup, Custom Types, Concentrations and Data Locations tabs.

Plate Setup Preferences

Tab 1: Plate Setup Preferences

A screenshot of the 'Preferences' window in TRACE Analysis Software, specifically the 'Plate setup' tab. The window has a title bar with a close button (X). Below the title bar is a tabbed interface with the following tabs: 'Plate setup' (selected), 'Custom types', 'Concentrations', 'Data and Reports', 'Data locations', 'Material tracking', 'Users', and 'Language'. The 'Plate setup' tab contains several settings: 'Machine for genotyping' is a dropdown menu set to 'QIAcuity'; 'Machine for monitoring' is a dropdown menu set to 'QIAcuity'; 'Layout type' is a dropdown menu; 'Replicates' is a dropdown menu set to 'dPCR Singletons'; 'Color mode' is a dropdown menu set to 'Pastel'; 'Plate for genotyping' is a dropdown menu set to 'QIAcuity, v1, 96'; and 'Plate for monitoring' is a dropdown menu set to '96-low density'. At the bottom right of the window are three buttons: 'Ok', 'Cancel', and 'Apply'.

“Machine for genotyping”:

Define the default genotyping analysis platform by clicking the radio button appropriate for the machine in use. TRACE Analysis™ Software will generate the appropriate sample setup file for the machine, as well as be able to analyze the results exported from that machine.

“Machine for monitoring”:

Define the default monitoring analysis platform by clicking the radio button appropriate for the machine in use. TRACE Analysis™ Software will generate the appropriate sample setup file for the machine, as well as be able to analyze the results exported from that machine.

“Layout type”:

Define whether the samples are placed onto the plate in groups or in rows.

We strongly recommend organizing samples by groups, as this makes the pipetting scheme very easy to follow.

“Replicates”:

Define whether to perform quantitative analysis using singletons, duplicates or triplicates.

For a qPCR monitoring, JETA Molecular recommends the use of triplicates, as low concentration DNA samples may exhibit stochastic variance in Cq.

“Color mode”:

The color highlights which appear on the software interface and the protocols can be modified to use a scheme which is preferred by the user.

“Plate for genotyping”:

This preference selects the plate chemistry version in use for genotyping. Different plate versions have slightly different assay panel members. If you have created Custom Genotyping Panels, they will appear here.


“Plate for monitoring”:

This preference selects the plate version in use for monitoring. This allows you to select among different plate types used in quantitative analysis by digital PCR.

Custom Types Preferences


In the Custom types tab of the Preferences Menu, users may define the cell types being used, the diseases associated with the samples as well as any external reference materials which may be used.

Tab 2: Custom Types Preferences



Cell Types tab of the Custom types tab of the Preferences Menu:


Define the origin of the materials being tested. You can type a sample material in the empty box near the Add button, then press Add, and the sample type will now appear as an option in TRACE Analysis™. The small "x" on the same line as the sample type is used to remove that type from the software. It is best to add all known types from the earliest use of TRACE Analysis™.

Diseases tab of the Custom types tab of the Preferences Menu:

Define diseases associated with the samples to be tested.

Type the disease name in the empty box near the Add button, then press Add, and the disease will now appear as an option in TRACE Analysis™ Software. The small "x" on the same line as the disease name can be used to remove that disease from the software.

It is best to add all known diseases from the earliest use of TRACE Analysis™ Software.

External references tab of the Custom types tab of the Preferences Menu:

Define external reference DNAs to be used in the monitoring using qPCR.

Type the reference name in the empty box near the Add button, then press Add, and the reference will now appear as an option in TRACE Analysis™ Software. The small "x" on the same

line as the reference name can be used to remove that reference from the software. It is best to add all known references from the earliest use of TRACE Analysis™ Software.

Concentrations Preferences

Tab 3: Concentration Preferences

Setting	Value	Unit
Correction for excess master mix	15	%
Correction for excess DNA dilution	15	%
Sample Input for Genotyping	10	ng
Reference Sample Input for Monitoring (qPCR)	10	ng
PostTx Sample Input for Monitoring	150	ng
Default concentration of samples	50	ng/μL

Reset

Ok Cancel Apply

Concentrations tab of the Preferences Menu:

Define the variables used in protocol generation and experimental execution.

Define how much excess master mix and how much excess DNA dilution to use in the experimental protocol.

TRACE Analysis™ calculates the volumes needed, based on the experimental setup and then adds these additional factors to provide more than enough of each solution to execution the experiment.

Define the sample input for genotyping using different amplification methods. Based on the data from verification studies, JETA Molecular recommends the use of 5ng DNA input per well for genotyping using qPCR and 10ng for genotyping using dPCR.

Define the amount of reference sample used for monitoring using qPCR

Define the sensitivity to achieve in the experiment by entering the target DNA input per well for the monitoring samples.

While sensitivity may be thought of in terms of cell numbers or percentages or grams of DNA, the software uses nanograms of input to generate an appropriate protocol.

Define the default concentrations of samples.

This number will appear for all samples and can be altered if necessary.

Data and Reports Preferences

Tab 4: Data and Report Preferences

qPCR tab:

The screenshot shows the 'Preferences' dialog box with the 'Data and Reports' tab selected. Within this tab, the 'qPCR' sub-tab is active. The 'Replicate Highlighting Method' section has three radio buttons: '% CV' (set to 2), 'SD' (selected, set to 1), and 'Cq Range' (set to 0.5). The 'Report Highlighting' section has four checkboxes, all of which are unchecked: 'Replicate Highlighting', 'Reference Sample Cq Range' (with Low: 24, High: 28), 'Reference Sample ΔCq' (with Low: -1.5, High: 1.5), and 'Monitoring Sample RNaseP Cq Range' (with Low: 20, High: 24). The 'Laboratory Information' section contains text boxes for Department (Technical Support), Institution (JETA Molecular), Address (Krommewetering 101C), Postal code, City (3543 AN, Utrecht), Telephone (+316 5413 6697), and Comment. The 'Anonymous Reporting' section has an unchecked checkbox 'Allow Anonymous Reporting'. The 'HPRIM Reporting' section has an unchecked checkbox 'Allow HPRIM Reporting'. At the bottom are 'Ok', 'Cancel', and 'Apply' buttons.

In order to achieve more customized data analysis, TRACE Analysis™ has a tab called Data and Reports in the Preferences section of the software.

There is a Replicate Highlighting Section and a Report Highlighting Section.

In Replicate Highlighting, users can select between three methods to highlight variance in replicate samples. Replicates displaying variance in %CV, SD or Cq range may be user defined, based on the lab's practices. The chosen method is applied to all replicates in monitoring, whether it is the reference sample or the monitoring sample.

In Report Highlighting, users can select to have data highlighted on the experimental reports, if the data falls outside of input ranges. The user defines which ranges are critical and these will be readily apparent when data reports are reviewed

In Report Highlighting, users have the ability to choose to have portions of their report highlighted, if values generated fall outside an expected range. The user can choose up to four different data quality inspections to be performed on monitoring data, with values outside of the input ranges highlighted automatically on the report

Protocols and Reports can be customized by adding institutional information under the Laboratory Information section. Any information entered here will appear at the top of all protocols and reports generated by TRACE Analysis™ Software.

dPCR tab:

The image shows two side-by-side screenshots of the 'Preferences' dialog box, specifically the 'dPCR' tab. Both windows have the same layout with tabs for 'Plate setup', 'Custom types', 'Concentrations', 'Data and Reports', 'Data locations', 'Material tracking', 'Users', and 'Language'. The left window is for droplet-based digital PCR, showing 'Total valid droplets/partitions' with a low of 8000 and a high of 30000. The right window is for plate-based digital PCR, showing a low of 3000 and a high of 9000. Both windows include fields for 'Post transplant RNaseP concentration', 'NTC concentration', and 'UPC concentration'. There are also checkboxes for 'Anonymous Reporting', 'HPRIM Reporting', and 'Use UPC'. The 'Laboratory Information' section includes fields for 'Department', 'Institution', 'Address', 'Postal code, City', and 'Telephone'.

The number of valid droplets/partitions, post transplant RNase P concentration, UPC and NTC concentration can be defined here. If the actual droplet/partition number or RNaseP, NTC or UPC concentration falls outside of the defined range, data will be highlighted in the report.

In the report highlighting settings, the total valid droplet/partitions values would be different for plate-based digital PCR (example on the right) and for droplet-based digital PCR (example on the left).

The use of a Universal Positive Control in a dPCR monitoring can be enabled or disabled here.

Data Location Preferences

Tab 5: Data locations Preferences

The image shows a screenshot of the 'Preferences' dialog box, specifically the 'Data locations' tab. The tab is selected, and the text 'Choose the location where the software stores your data files' is displayed. Below this, there is a text field labeled 'Data store location' containing the path 'C:\Users\jeroen\Dropbox (JETA)\JETA Team Folder\JETA\02'. To the right of the text field is a 'Default' button. At the bottom of the dialog are 'Ok', 'Cancel', and 'Apply' buttons.

Data locations tab of the Preferences Menu:

Define where TRACE Analysis™ Software stores the data it generates.


For labs using multiple copies of TRACE Analysis™ Software, this location is likely best set as a shared location on a server. When all local copies of TRACE Analysis™ Software point to the same data storage location, all copies can read and write to the same data files, eliminating the need to transfer files between computers.

Reagent Tracking Preferences

This tab is where reagent lot numbers and expiration dates may be entered into TRACE Analysis™ Software. This information will populate protocols and reports, eliminating the need for the operator to write the information each time.

qPCR Reagent tracking tab


This tab enables to enter reagent lot numbers and expiration dates for qPCR reagents.



Name	Lot Number	Expiry Date	Comment
QTRACE INDEL Assay 634	P10700	31/07/2023	
QTRACE INDEL Assay 650	N00580	31/03/2022	
QTRACE INDEL Assay 663	Q10490	30/09/2024	
QTRACE INDEL Assay 670	P10450	28/02/2023	
QTRACE INDEL Assay 678	R11370	31/05/2025	
QTRACE INDEL Assay 681	R11380	31/05/2025	
QTRACE INDEL Assay 694	R11390	31/05/2025	
QTRACE INDEL Assay 706	N00060	31/12/2021	
QTRACE INDEL Assay 710	N01910	31/08/2022	
QTRACE INDEL Assay 721	N01920	31/08/2022	

dPCR Reagent tracking tab

This tab enables to enter reagent lot numbers and expiration dates for dPCR reagents.



Name	Lot Number	Expiry Date	Comment
d102	T10580	31/03/2026	
d113	S10670	30/06/2026	
d120	T10870	31/05/2027	
d137	S10680	30/06/2026	
d148	S10230	31/03/2026	
d157	T10210	31/03/2027	
d176	T10590	30/04/2027	
d183	S11300	30/11/2026	
d198	S11100	31/10/2026	
d209	S10110	31/01/2026	




Equipment & consumables tracking tab

This tab enables to enter serial number and location information for equipment and consumables.

Preferences

Plate setup Custom types Concentrations Data and Reports Data locations Material tracking **Users** Language

qPCR Reagent tracking dPCR Reagent tracking Equipment & consumables tracking

Name	Lot Number	Expiry Date	Comment	Print
ABI 7500	SN 27500XXX	21/10/2021 	Genetic Analysis Lab	<input checked="" type="checkbox"/>
P10	SN YXX	25/02/2021 	DNA Lab	<input checked="" type="checkbox"/>
P100	SN XXXX	24/03/2021 	PCR Setup Lab	<input checked="" type="checkbox"/>

Add Delete

Ok Cancel Apply

Users Preferences

This tab is where User profiles are managed by an individual with 'Supervisor' rights in TRACE Analysis™ Software.

Preferences

Plate setup Custom types Concentrations Data and Reports Data locations Material tracking **Users** Language

User Management


Name	Role
Technician	Analyst
Advanced User	AdvancedUser
Lab Director	Supervisor
Technician 2	Analyst

Change Password Add User Delete User

Ok Cancel Apply

Language Preferences

This tab allows the users to select among available language versions.




Software Access And User Profile Management

TRACE Analysis™ Software restricts access to the software and well as functions within the software.

When TRACE Analysis™ launches for the first time, User Profiles need to be established.

Once logged in with JETA's credentials, a Supervisor account needs to be created.

In the Preferences, under the new User tab, there is an option for creating New Users.



By clicking 'Add User' a dialog box opens wherein the details of the new user can be added.

The 'Add User' dialog box is shown with the following details:

- Name:** Lab Director
- Password:** [Redacted]
- Confirm password:** [Redacted]
- Role:** Supervisor (selected from a dropdown menu)
- Buttons:** Save, Cancel

After the Supervisor has created an account, multiple additional "Supervisor," "Advanced User," and "Analyst" accounts can be made.

The 'Preferences' dialog box, 'Users' tab, displays the following information:

Name	Role
Technician	Analyst
Advanced User	AdvancedUser
Lab Director	Supervisor
Technician 2	Analyst

Buttons below the table: Change Password, Add User, Delete User.

Bottom buttons: Ok, Cancel, Apply.

Analysts and Advanced Users are allowed to change their passwords, and to edit the Reagent tracking Preferences information. All other Preferences settings are controlled by the Supervisor. Advanced Users inspect and can validate a result.

	Edit all tabs in Preferences	Results Validation	Add or Delete Users/Passwords	Edit Reagent Tracking Data
Supervisor	Yes	Yes	Modify All Users and Passwords	Yes
Advanced User	No	Yes	Modify Personal Password	Yes
Analyst	No	No	Modify Personal Password	Yes


File Management

There are a few files that TRACE Analysis™ stores to enable tracking of data, logging of errors and supply a list of assays from JETA Molecular's kits.


Data storage and retrieval in TRACE Analysis™ Software version 2.0 is enhanced through the introduction of a relational database data storage format. Previous versions of software relied on generation of .xml files which were stored in various folders, making the information captured in the software difficult to access.

Information is now captured and stored in a database, making SQL queries of the data possible. This also facilitates enhanced performance of the software, allowing users more ability to modify or update experiments and sample records, as needed.

In order to query records in the database, under 'File-->Custom reports' leads to a search interface for TRACE Analysis™ Software.



Search queries can be built, executed and stored. The results of the queries can also be easily exported to .csv and .xls format files. The image below shows the results for the number of times an assay is informative for recipients in the database.



For the logging of errors, files are created in your temporary directory (C:\Users\Your Profile\AppData\Local\Temp), named QTRACE.log(n).

These file store per session in TRACE Analysis™ Software what actions the user performs. In the case of an error these log files can be send to support to detect the cause and solve the issue.


The list of assays from JETA's kits is stored by default at C:\ProgramData\QTRACE\QTRACEKit\JETAMarkers.txt. This file is created during installation and must at all times be present and unmodified for TRACE Analysis™ Software to work.


Experimental Files Saved in TRACE Analysis™ Software

TRACE Analysis™ Software automates the storage of all files associated with a recipient’s testing history by automatically creating folders which take the form: recipient Name_recipient ID


QTRACE > Data Folders > Gregor Mendel_7458		
Name	Date modified	Type
PCR Results Files	16/01/2025 09:14	File folder
PCR Run Files	16/01/2025 09:13	File folder
PCR Setup Files	11/11/2024 08:35	File folder
Protocols	11/11/2024 08:37	File folder
Reports	16/01/2025 09:14	File folder

Within each of these folders are all the files used to track the recipient. When PCR data is imported into TRACE Analysis™, the user also has the ability to identify the location of the original PCR machine file.

After clicking the "Import PCR Data"  button, TRACE Analysis™ first asks the user to identify the location of the exported results file. The imported file is automatically saved to the PCR Results File folder, while the data from the file is extracted and analyzed in TRACE Analysis™.




After selecting the results file, TRACE Analysis™ prompts the user to identify the location of the data collection file from the PCR machine. The imported file is automatically saved to the PCR Runs File folder.



The selection of the data file is optional, but if it is selected, TRACE Analysis™ will copy this file

and put it into the Run Files folder for each patient in the same experiment. With the reports that are generated from each experiment, they are differentially saved – only the report for the specific patient goes into that patient's file, while the other patient sample has a separate report, and it is saved in their separate folder.

TRACE Analysis™ Software also creates a folder called Exports. The data from every experiment is automatically converted into an xml file which is then easily available to integrate into a LIMS.



Data Exports - csv and xlsx

In addition to the previous .xml export files which were automatically created by TRACE Analysis™ Software, TRACE Analysis™ Software v1.08 now automatically generates .csv and .xlsx files containing the data from the experimental reports.

Name	Date modified	Type	Size
180614 L06300 UCLA801 T0.xlsx	08/04/18 4:33 PM	Microsoft Excel Worksheet	6 KB
180614 L06300 UCLA801 T0.csv	08/04/18 4:33 PM	Microsoft Excel Comma Separated Values File	3 KB
180614 L06300 UCLA801 T0.xml	08/04/18 4:33 PM	XML Document	14 KB

These new automated exports are stored in the Exports folder within the QTRACE folder.

The following is an example of the data found in the xlsx file from a Custom Genotyping experiment.

This file has recipient and donor IDs, informative assays, dCq values, assay status (positive, negative or atypical) and the mean Cq values for assay replicates.

	10831 21 Aug	2018-08-	10831 21 Aug						
1	12560 21 Aug	12560 21 Aug							
1	916 10q	12560 21 Aug	0.52	Negative	40	Positive	29.8002		
2	305 20q	12560 21 Aug	0.85	Negative	40	Positive	30.1296		
3	710 5q	12560 21 Aug	1.06	Negative	40	Positive	30.3468		
4	137 1p	12560 21 Aug	2.28	Negative	40	Positive	31.5592		
5	907 11q	10831 21 Aug	0.45	Positive	29.777	Negative	40		
6	520 20q	10831 21 Aug	0.78	Positive	30.1095	Negative	40		
7	209 5q	10831 21 Aug	2.39	Positive	31.7184	Negative	40		
8	361 12q	None	-	Negative	40	Negative	40		
9	706 16p	None	-	Positive	30.3095	Positive	30.2039		
10	755 11q	None	-	Negative	40	Negative	40		
11	NTC	14q	None	-	Negative	40	Negative	40	
12	POS	14q	None	-	Positive	29.3304	Positive	29.2834	

The following is an example of an xlsx export from a monitoring experiment.

Again, samples are identified, the mean Cq values and the %DNA value for each sample and

assay combination are exported, in addition to the Mean % DNA for all assays tested on a single sample.

RECIPIENT	DNA 103 21 Aug	DNA 103 21 Aug	2018-08- DNA 103 21 Aug	
DONOR		2 DNA 103 21 Aug	DNA 103 21 Aug	
SAMPLE		1 DNA 103 21 Aug	Reference	
MARKER	DNA 103 21 Aug	1	148	30.03
MARKER	DNA 103 21 Aug	2	386	30.99
MARKER	DNA 103 21 Aug	3	748	29.48
MARKER	DNA 103 21 Aug	4	RNaseP	29.85
SAMPLE		2 Post 1 21 Aug	Sample	
MARKER	Post 1 21 Aug	1	148	29.68
MARKER	Post 1 21 Aug	2	386	30.83
MARKER	Post 1 21 Aug	3	748	29.17
MARKER	Post 1 21 Aug	4	RNaseP	25.57
SAMPLE		3 Post 2 21 Aug	Sample	
MARKER	Post 2 21 Aug	1	148	31.82
MARKER	Post 2 21 Aug	2	386	33.04
MARKER	Post 2 21 Aug	3	748	31.48
MARKER	Post 2 21 Aug	4	RNaseP	24.41
SAMPLE		4 Post 3 21 Aug	Sample	
MARKER	Post 3 21 Aug	1	148	37.63
MARKER	Post 3 21 Aug	2	386	38.1
MARKER	Post 3 21 Aug	3	748	38.1
MARKER	Post 3 21 Aug	4	RNaseP	26.74
CELL_FRACTION		1 Post 1 21 Aug	148	6.58
CELL_FRACTION		2 Post 1 21 Aug	386	5.72
CELL_FRACTION		3 Post 1 21 Aug	748	6.36
CELL_FRACTION		4 Post 2 21 Aug	148	0.67
CELL_FRACTION		5 Post 2 21 Aug	386	0.56
CELL_FRACTION		6 Post 2 21 Aug	748	0.58
CELL_FRACTION		7 Post 3 21 Aug	148	0.06
CELL_FRACTION		8 Post 3 21 Aug	386	0.08
CELL_FRACTION		9 Post 3 21 Aug	748	0.03
CELL_FRACTION_MEAN		1 Post 1 21 Aug	6.22	
CELL_FRACTION_MEAN		2 Post 2 21 Aug	0.6	
CELL_FRACTION_MEAN		3 Post 3 21 Aug	0.06	

HPRIM Data Export

Fonction d'export HPRIM

- a) Activer cette fonction dans 'Preference\Data and Reports', en cochant 'Allow HPRIM Reporting'

The screenshot shows the 'Preferences' dialog box with the 'Data and Reports' tab selected. Under the 'Report Highlighting' section, the 'Allow HPRIM Reporting' checkbox is checked and circled in red. Other settings include 'Replicate Highlighting Method' set to 'SD' and 'Report Highlighting' set to 'Low'.

- b) Créer un nouveau couple donneur- receveur

Les nouveaux champs sont encadrés en rouge.

The screenshot shows the 'Recipient' and 'Donor' forms. The 'Recipient' form has fields for Recipient First Name, Recipient Last Name, Recipient ID, Sample ID, DNA ID, G-LIMS ID, Concentration ng/ul, Date of Birth, Date of Transplant, Gender, Comment, and Disease Type. The 'Donor' form has fields for Donor First Name, Donor Last Name, Donor ID, Sample ID, DNA ID, Concentration ng/ul, Date of Birth, Gender, and Comment. The new fields (DNA ID, G-LIMS ID, and Donor DNA ID) are highlighted with red boxes.

Remarque importante:

Le numéro 'DNA ID' est un champ ajouté pour faciliter la préparation de la plaque. Le personnel technique du laboratoire peut y indiquer le numéro de son choix, celui-ci est affiché uniquement sur le protocole et sur le plan de plaque. Si celui-ci n'est pas renseigné lors de la création du patient, c'est le 'Sample ID' qui sera affiché sur le protocole.

- c) Renseigner le nouveau champ '**Q-LIMS ID**'. Il apparait dans le bloc 'Recipient' après avoir sélectionné l'échantillon de référence pour le test de quantification.

The 'Recipient' form contains the following fields: Recipient First Name, Recipient Last Name, Recipient ID, Sample ID, DNA ID, G-LIMS ID (highlighted with a red rectangle), Concentration ng/ul (set to 100), Date of Birth (format XX-XX-XXXX), Date of Transplant (format XX-XX-XXXX), Gender (Male/Female radio buttons), Comment, and Disease Type (dropdown menu).

- d) Renseigner le nouveau bloc '**Prior Sample**'.

Celui-ci permet de quantifier l'échantillon précédent. La préparation de cet échantillon sur la plaque ainsi que le calcul du résultat sont effectués comme un échantillon 'Post Sample'.

Le % ADN est exporté dans le fichier HPRIM et il est imprimé sur le rapport.


Attention: Le résultat de l'échantillon '**Prior Sample**' n'est pas affiché sur le graphique de suivi.

The 'Reference Samples' table shows:

Sample ID	Concentration ng/ul	% reference
	100	100
	100	100
Pavel example	10	100
Ref Sample	10	100

The 'Prior Sample' form (highlighted with a red rectangle) contains: Prior Sample ID, DNA ID, Concentration ng/ul (set to 100), Prior Sample Date (format XX-XX-XXXX), and Prior Sample Type (radio buttons for BM, Blood, T Cells, B Cells, DNA). Below it is the 'Add Prior Sample' button. The 'Post Sample' form contains: Sample ID, DNA ID, Concentration ng/ul (set to 100), Sample Date (format XX-XX-XXXX), Sample Type (radio buttons for BM, Blood, T Cells, B Cells, DNA), and Donor Sample ID.

- e) Valider les résultats

- Après avoir importé les résultats de qPCR de votre patient et après avoir vérifié les données de votre plaque dans la fenêtre « Data analysis », cliquer sur le bouton 'Calculate' .

- Le logiciel affiche alors les résultats de votre génotypage ou de votre quantification dans la fenêtre 'Result'.


Exemple: résultat de génotypage

Result						
RECEVEUR 210 RECEVEUR 210						
Marker	Chr. location	InfoCq	ΔCq	Informative for	RECEVEUR 210 RECEVEUR 210	DONNEUR 210 DONNEUR 210
209	5q	31.46	2,6465	RECEVEUR 210 RECEVEUR 210	Positive	Negative
235	2q	31.03	2,2197	RECEVEUR 210 RECEVEUR 210	Positive	Negative
267	4q	30.7	1,6455	DONNEUR 210 DONNEUR 210	Negative	Positive

Exemple: résultat de quantification

Result							
4083384 FIRST 4083384 LAST							
Sample	Type	Date	Target	Chr.	ddCq	DNA (%)	Reference
4081082	Blood	05-05-2019	356	18q	4.45	4.59	4083384
4081082	Blood	05-05-2019	721	Xq	4.77	3.67	4083384
4099253	B Cells	09-05-2019	356	18q	6.52	1.09	4083384
4099253	B Cells	09-05-2019	721	Xq	3.97	6.37	4083384

- Cliquer ensuite sur le bouton 'Select assays' **Select assays**. Fonction autorisée ou non selon les droits de l'utilisateur.
- Une fenêtre s'ouvre pour sélectionner les marqueurs. Vous devez sélectionner les marqueurs de votre choix.




- Cliquer ensuite sur 'Save changes' pour valider vos résultats. Les boutons 'HPRIM Export' et 'Report' apparaissent alors.

Select assays **HPRIM Export** **Report...**

f) Exporter un fichier HPRIM

- Cliquer sur le bouton 'HPRIM Export'. La fenêtre ci-dessous s'ouvre, puis choisir le dossier dans lequel le fichier HPRIM sera enregistré:



Software Buttons



Home Screen

Return to the home screen of TRACE Analysis™ Software



Add New Sample

Add a new recipient record to the database. A record will always require a recipient name, unique recipient identifier and a sample identifier.



Add Typing Samples

After entering recipient and donor specific information, choosing Screen will start creating a genotyping experiment by placing your samples on the plate.



Add Quant Samples

After selecting a recipient to monitor and entering sample specific information, the Quant button adds your sample to a monitoring experiment



Export Setup to PCR

Export an experiment sample setup file for use with a PCR instrument.



Import PCR Data

Import the .txt or .csv results file from your PCR instrument to review the data collected.



Overview

Shows all available data for a transplantation: informative markers and quantitative analyses.

**Virtual Typing**

Enables side-by-side comparison of genotypes of selected samples.

**Browse Experiments**

Browse all experiments previously created in TRACE Analysis™ Software.

By pressing the open button, you can re-open the imported data files.

**Save**

Update information about the recipients or donors.

Do not update information about the recipient or donor if there is PCR data waiting to be analyzed for them.

**Print**

Print all typing and monitoring results from the currently selected recipient data.

**Preferences**

Set your preferred instrument, sample types, number of replicates and disease types.

**Help**

Review this manual directly via TRACE Analysis™ Software to search for useful tips, tricks and troubleshooting.

**About**

Technical information about TRACE Analysis™ Software, such as version, license, contact information.



Add Sample

Add another sample.



Remove

Remove a sample.



Reset

Reset the experiment plate completely.



Undo

Reset the last placement onto the experimental plate.




Redo

Reset the previous "Undo" action onto the experimental plate.


Drag and Drop Plate Layouts

TRACE Analysis™, Version 1.5 introduced movable wells for customized plate layout capability. Once monitoring samples are placed onto the plate, the user has the ability to move a well from one location to another location. Whole sample groups can be moved by left-clicking and holding the mouse button. This action treats all the wells in a particular sample group the same, moving the entire block of wells. Right-clicking and holding the mouse button allows a user to move one well at a time.


In a monitoring experiment, once the Quant button is pressed, samples and assays are added to the plate.



Entire Sample Blocks can be moved by holding the left-click mouse button while touching a well in the layout



Individual Sample Blocks can be moved by holding the right-click mouse button while touching a well in the layout



Anonymized Reporting

The ability to generate anonymized reports from TRACE Analysis™ Software is an option in the Preferences.

On the 'Data and Reports' tab of the Preferences, there is a check box which allows for reports to be generated without the names of the recipient and donor(s) appearing on them.

The screenshot shows the 'Preferences' dialog box with the 'Data and Reports' tab selected. Under 'Replicate Highlighting Method', 'SD' is selected with a value of 1. Under 'Report Highlighting', 'Replicate Highlighting' is checked. In the 'Laboratory Information' section, fields for Department, Institution, Address, Postal code, City, Telephone, and Comment are visible. In the 'Anonymous Reporting' section, the 'Allow Anonymous Reporting' checkbox is checked, and the 'HPRIM Reporting' checkbox is unchecked. At the bottom are 'Ok', 'Cancel', and 'Apply' buttons.

With the 'Allow Anonymous Reporting' option checked, Reports generated from TRACE Analysis™ Software will have the Recipient ID or the Donor ID replacing the name of the individuals.

The following is an example of a Transplantation Report (Recipient data overview) generated with the normal settings and the same report when 'Allow Anonymous Reporting' is checked. In this example, the recipient's name is "Receveur" and the donor's name is "Donneur." "EPT2017R" is the sample ID of the recipient pre-transplant sample.

Chim9-2	DNA	03-Nov-2017	137	1p	10.31	0.08	EPT2017R	Receveur
Chim9-2	DNA	03-Nov-2017	235	2q	10.43	0.07	EPT2017R	Receveur
Chim9-2	DNA	03-Nov-2017	434	3q	10.32	0.08	EPT2017R	Receveur
Mean:					0.08			

Chim10-2	DNA	03-Nov-2017	137	1p	9.35	0.15	EPT2017R	Receveur
Chim10-2	DNA	03-Nov-2017	235	2q	9.54	0.13	EPT2017R	Receveur
Chim10-2	DNA	03-Nov-2017	434	3q	9.58	0.13	EPT2017R	Receveur
Mean:					0.14			

Informative Assays

Assay	Locus	Informative for	InfoCq	ΔCq	Recipient	CqDonor 1	Cq
907	11q	Receveur	31.69	3.2	Positive	31.7	40.0
854	12q	Receveur	29.32	0.83	Positive	29.3	40.0
840	9q	Donneur	29.54	1.71	Negative	39.0	29.5
824	22q	Receveur	28.76	0.28	Positive	28.8	38.9
788	18q	Donneur	30.87	3.04	Negative	40.0	30.9
796	1p	Donneur	29.15	1.32	Negative	40.0	29.2
710	5q	Donneur	30.43	2.6	Negative	40.0	30.4
706	15p	Donneur	29.62	1.79	Negative	40.0	29.6
690	7q	Donneur	29.96	2.13	Negative	40.0	30.0
634	11q	Donneur	28.99	1.16	Negative	40.0	29.0
548	1q	Donneur	28.33	0.5	Negative	40.0	28.3
531	17p	Donneur	29.57	1.74	Negative	40.0	29.6
520	20q	Donneur	29.19	1.36	Negative	40.0	29.2
504	17p	Donneur	29.55	1.72	Negative	40.0	29.6

The following is the same report with anonymization of the individuals in the report. Where names existed before, now the sample ID of the individual is displayed.

Chim9-2	DNA	03-Nov-2017	137	1p	10.31	0.08	EPT2017R	[ID:EPT2017R]
Chim9-2	DNA	03-Nov-2017	235	2q	10.43	0.07	EPT2017R	[ID:EPT2017R]
Chim9-2	DNA	03-Nov-2017	434	3q	10.32	0.08	EPT2017R	[ID:EPT2017R]
Mean:						0.08		

Chim10-2	DNA	03-Nov-2017	137	1p	9.35	0.15	EPT2017R	[ID:EPT2017R]
Chim10-2	DNA	03-Nov-2017	235	2q	9.54	0.13	EPT2017R	[ID:EPT2017R]
Chim10-2	DNA	03-Nov-2017	434	3q	9.58	0.13	EPT2017R	[ID:EPT2017R]
Mean:						0.14		


Informative Assays

Assay	Locus	Informative for	InfoCq	ΔCq	Recipient	CqDonor 1	Cq
907	11q	[ID:EPT2017R]	31.69	3.2	Positive	31.7	Negative
854	12q	[ID:EPT2017R]	29.32	0.83	Positive	29.3	Negative
840	9q	[ID:EPT2017D]	29.54	1.71	Negative	39.0	Positive
824	22q	[ID:EPT2017R]	28.76	0.28	Positive	28.8	Negative
768	18q	[ID:EPT2017D]	30.87	3.04	Negative	40.0	Positive
736	Xp	[ID:EPT2017D]	29.15	1.32	Negative	40.0	Positive
710	5q	[ID:EPT2017D]	30.43	2.6	Negative	40.0	Positive
706	16p	[ID:EPT2017D]	29.62	1.79	Negative	40.0	Positive
650	7q	[ID:EPT2017D]	29.96	2.13	Negative	40.0	Positive
634	11q	[ID:EPT2017D]	28.99	1.16	Negative	40.0	Positive
548	Xq	[ID:EPT2017D]	28.33	0.5	Negative	40.0	Positive
531	17p	[ID:EPT2017D]	29.57	1.74	Negative	40.0	Positive
520	20q	[ID:EPT2017D]	29.19	1.36	Negative	40.0	Positive
504	17p	[ID:EPT2017D]	29.55	1.72	Negative	40.0	Positive


Modification of Existing Data

TRACE Analysis™ Software v2.0 allows for modification of some data associated with recipients, samples, transplantation dates and experiments.

To make changes to existing files, select the Data Management item from the File Menu



Selecting from the four different buttons allows users to edit some information related to these topics.



The 'D/R Info' button allows editing of Donor and Recipient information.

Sample Information may be edited by pushing the 'Samples' button.

Sample

SampleID

ExperimentID

3rd Party SampleID

Acquisition Date

Material

Concentration

Sample Percentage

DNA ID

Donor Sample ID

XX-XX-XXXX

15

Cancel

Save

After selecting sampleID select experimentID if sample is used in multiple experiments.

Change of sampleID is propagated to all samples with this sampleID.

The 'Transplantations' button allows edits to Date of Transplant.

Transplantation

Recipient ID

Date of Transplant

Date of Second Transplant

XX-XX-XXXX

15

XX-XX-XXXX

15

Export

Cancel

Save

☐ Anonymize export

The 'Experiments' button allows experiments to be edited and deleted.

Experiment

Name

Operator ID

Date of experiment

XX-XX-XXXX

15

Cancel

Save

Delete

Warning! This will delete the experiment with all its related data with no return.

Export

☐ Anonymize export

Data Analysis Algorithms

The algorithms implemented in the TRACE Analysis™ Software are designed to provide reliable analysis results based upon specific data criteria. Quality measurements, based upon the same data criteria, provide the operator with valuable information regarding data integrity.

The algorithms and quality measurements are invoked when the operator initiates data analysis through the TRACE Analysis™ Software. Separate algorithms and quality measurements exist for both the genotyping and quantitation tests. Warning messages generated from the quality measurements are provided when calculated results fall outside established data criteria ranges.

The QTRACE® Analysis System

The data criteria used in the algorithms and quality measurements are not operator configurable. Cq data will need to be exported from the qPCR Software and manipulated in a spreadsheet program if another method of data analysis is desired.

Genotyping Test Algorithm

The genotyping test algorithm determines the marker state for each assay/sample reaction followed by the informative state for each assay. Quality measurements are performed after the algorithm has been applied to the data. There are three possible marker states for each reaction in regards to amplification: positive, negative, or atypical.

Marker States

Positive

An Amplification Control reaction (positive control with QTRACE® RNaseP Assay; wells D11 and H11 in the standard QTRACE® Genotyping Plate) is positive for amplification with a Cq value less than 34.0.

A reaction with any QTRACE® INDEL Assay is positive for amplification if the ΔCq (the difference in Cq values between the reaction and the Amplification Control) is within an acceptable range. If the Amplification Control is positive, then each of the reactions on the plate must have a ΔCq of -2.0 to +3.5 to be considered positive.

Note: An NTC reaction (negative control with QTRACE® RNaseP Assay; wells D12 and H12) can only be negative or atypical.

Negative

An Amplification Control reaction is negative for amplification with a Cq value greater than or equal to 34.0.

All other reactions including the NTC reaction are negative with a Cq value greater than or equal to 38.0.

Atypical

An NTC reaction is classified as atypical if the Cq value is less than 38.0. A “+” sign for an NTC reaction represents the atypical result in the Analysis View. This is assigned to atypical NTC reactions when quality measurements are assessed.

A reaction with any QTRACE® INDEL Assay will have atypical amplification if the ΔC_q (the difference in C_q values between the reaction and the Amplification Control) is outside the acceptable ranges defined in the algorithm settings outlined above under the Marker States – Positive section. The exception to this rule is if the assay has already been classified as negative for amplification in which case the marker state will remain as negative. A large number of atypical reactions for a sample is usually indicative of low input copy number (DNA quantity) for that particular sample.

Note: An Amplification Control reaction can only be positive or negative.

Informative States

Once marker state has been determined, the genotyping test algorithm compares the marker states between two DNA samples on a plate for each assay to determine informativity. There are three possible informative states for any QTRACE® INDEL Assay that can be applied to each assay result:

informative, not informative, or undetermined. There are two possible informative states for the control assays using QTRACE® RNaseP Assay: pass or fail.

Informative

An assay is informative for a pair of DNA samples if one of the samples is positive for amplification and the other sample is negative.

Not Informative

An assay is not informative for a pair of DNA samples if both samples are positive for amplification or if both samples are negative for amplification.

Undetermined

An assay is undetermined for a pair of DNA samples if one of the two samples has an atypical reaction result. Assays with an undetermined result are not informative and should not be selected for use in quantitation for the associated pair of DNA samples.

Pass

The Amplification Control reactions in wells D11 and H11 will be labeled as pass if they are positive for amplification. The NTC reactions in wells D12 and H12 will be labeled as pass if they are negative for amplification.

Fail

The Amplification Control reactions in wells D11 and H11 will be labeled as fail if they are negative for amplification. The NTC reactions in wells D12 and H12 will be labeled as fail if they are atypical for amplification.

Genotyping Test Quality Measurements

There are two types of quality measurements provided in the genotyping test analysis: Amplification Control and NTC. Genotyping test quality measurements are not used during application of the algorithm to the data. They are presented in the Report View to inform the operator of data integrity. A warning message providing more detail accompanies each quality measurement that fails.

Amplification Control

The Amplification Control quality measurement serves as a positive PCR control for each sample.

The monomorphic RNase P locus for QTRACE® RNaseP Assay is present in all samples. The Amplification Control reactions are located in wells D11 and H11: one reaction for each sample. The result for this control is determined to be either pass or fail.

The Cq threshold for the Amplification Control is 34.0. If the Amplification Control has a Cq value of less than 34.0 then the Amplification Control is determined as pass. If the Amplification Control has a Cq value of greater than or equal to 34.0 then the Amplification Control is determined as fail.

If the Amplification Control for a particular sample fails, then all reactions for that sample that have a Cq value less than 38.0 will have an atypical marker state. This occurs because no reliable ΔCq can be calculated for any of the reactions run with this sample. It follows that the resulting informative state for each assay will be labeled as undetermined.

NTC

The NTC (No Template Control) quality measurement is determined as pass or fail per sample. The NTC reactions are located in wells D12 and H2: one reaction for each sample. In addition to the Report View.

The Cq threshold for the NTC is 38.0. If the NTC has a Cq value greater than or equal to 38.0 then the NTC is determined as pass. If the NTC has a Cq value less than 38.0 then the NTC is determined as fail.

Quantification Test Algorithm

The quantification test algorithm performs the steps necessary to generate percent determinations for the DNA of interest in a mixed sample.

Relative Quantification Background

Four possible combinations of samples and assays comprise a single quantification determination.

There will always be a reference sample and one unknown sample referred to generically as “sample”. The reference sample is most commonly a DNA specimen from the individual whose DNA is to be quantified in the mixed unknown sample. Relative quantification requires that a reference sample mean Cq value is compared to the sample mean Cq value. There will always be a reference assay (RNaseP Assay) and one informative quantification assay referred to generically as “assay”.

The reference assay is used to normalize for total DNA input amount among all reactions. To obtain the four possible combinations of samples and assays, the reference sample and unknown sample are each amplified with the reference assay and informative quantitative assay: Reference Sample/Reference Assay, Reference Sample/Assay, Sample/Reference Assay, and Sample/Assay. These reaction combinations all belong to a common Transplantation. A Transplantation contains a unique reference sample and all associated samples to be quantified in relation to the reference sample.

Quantification is performed using an established relative quantification method for real-time PCR also known as the $\Delta\Delta Cq$ method (Livak and Schmittgen, 2001). A sample calculation is outlined in the table below.

QTRACE® INDEL Assay 854 is the Assay and QTRACE® RNaseP Assay is the Reference Assay.

Well	Sample Name	Sample Type	Assay	Cq	Mean Cq	ΔCq	$\Delta\Delta Cq$	$2^{-[\Delta\Delta Cq]}$	Mean % DNA
B4	Sample 1	Post Sample	854	30.23	30.19		7.17	0.0069	0.69%
B5	Sample 1	Post Sample	854	30.12					
B6	Sample 1	Post Sample	854	30.23					
A4	Sample 1	Post Sample	RNaseP	21.72	21.74	8.46			
A5	Sample 1	Post Sample	RNaseP	21.63					
A6	Sample 1	Post Sample	RNaseP	21.86					
B1	Sample 2	100% Reference Sample	854	23.67	23.67				
B2	Sample 2	100% Reference Sample	854	23.72					
B3	Sample 2	100% Reference Sample	854	23.63					
A1	Sample 2	100% Reference Sample	RNaseP	22.34	22.39	1.29			
A2	Sample 2	100% Reference Sample	RNaseP	22.38					
A3	Sample 2	100% Reference Sample	RNaseP	22.44					

The $\Delta\Delta Cq$ is calculated with the following formula:

(Mean Cq, Sample/Assay – Mean Cq, Sample/ReferenceAssay) –

(Mean Cq, ReferenceSample/Assay – Mean Cq, ReferenceSample/ReferenceAssay)

The sample and reference sample are each normalized for total DNA input by comparing the mean Cq value of the assay reactions to the mean Cq value of the reference assay reactions.

Next, the normalized sample reaction value is compared to the normalized reference sample reaction value to generate a relative quantification value or $\Delta\Delta Cq$. This term is then converted into a linear term which represents the fold change in DNA input amount between the sample and reference sample with the following formula: $2^{-[\Delta\Delta Cq]}$. The fold change term is multiplied by the proportion of the reference sample that is purely reference sample (usually 100%) to generate the final percent determination for the DNA of interest in the unknown sample.

Algorithm Steps

The mean and standard deviation of the Cq values for replicate reactions are first calculated. The standard deviation calculation requires at least three replicate values. The Reference Sample/Reference Assay quality measurement is applied next. If the Reference Sample/Reference Assay quality measurement passes, the $\Delta\Delta Cq$, and DNA percents are calculated, displayed and stored. The Reference Sample/Reference Assay quality measurement requires that the mean Cq for the Reference Sample/Reference Assay must be ≤ 32.0 . Samples which have a Cq ≥ 39.0 are designated as negative.

The mean and standard deviation of the Cq values for replicate reactions are first calculated. The Reference Sample/Reference Assay quality measurement is applied next. If the Reference Sample/Reference Assay quality measurement passes, $\Delta\Delta Cq$, fold change, DNA percent, and mean DNA percent (if multiple quantification assays were used for the same unknown sample) are calculated.

The DigitalTRACE™ Analysis System

The correct placement of the threshold is crucial for an accurate data analysis. Positive droplets/partitions (above the threshold) are scored as "1", and negative droplets/partitions (below the threshold) are scored as "0". This digital signal is used to calculate the starting DNA concentration by statistically analyzing the numbers of positive and negative droplets in a sample. The users are strongly advised to check thoroughly all the thresholds after each dPCR run, before importing the run results into the TRACE Analysis™ Software. Incorrect thresholding will lead to wrong results.

The values gained from the Universal Positive Control (UPC) amplification are not used in the calculations.

Genotyping Test Algorithm

Based on the CNV value, a marker can be scored as positive, negative or atypical.

Note: An NTC reaction (negative control containing the RNaseP Assay) can only be negative or atypical.

Calculation of CNV

$CNV = (Concentration\ Marker / Concentration\ Reference)^2$

Marker States

Positive

For CNV values between 0.55 and 3.2, a marker is considered positive.

Negative

For CNV values lower than 0.02, a marker is considered negative.

Atypical

For CNV values between 0.02 and 0.55, between 1.4 and 1.6 and higher than 3.2, a marker is considered atypical.

Quantification Test Algorithm

The ratio (CNV Marker/CNV Reference) gained from the dPCR software (QuantaSoft or QIAcuity Software Suite) is used to calculate the DNA percentage. The ratio is based on the number of positive and negative partitions and is refined automatically by the use of Poisson statistics.

$DNA\ \% = (Concentration\ Marker / Concentration\ Reference) * 100 * 2$ for assays with CNV=1

$DNA\ \% = (Concentration\ Marker / Concentration\ Reference) * 100$ for assays with CNV=2

Number of Partitions

During the partitioning step, the sample is divided into thousands of individual partitions that serve both for the calculation of the sample concentration (positive partitions), but also for the accurate statistical refinement of the calculated result (empty partitions). Problems with the partitioning process can result in lower partition numbers and, therefore, inaccurate quantification results.

Copy Number Adjustment

TRACE Analysis™ Software allows the user to change the CNV information for each sample. See the [Monitoring Data Analysis and Report](#) section for details.

Using the correct CNV is necessary to obtain the accurate quantification result. An incorrect CNV will lead to an underestimated or overestimated percentage value.

Bibliography

Alizadeh, et al. Quantitative Assessment of Hematopoietic Chimerism after Bone Marrow Transplantation by Real-Time Quantitative Polymerase Chain Reaction. *Blood*, 2002; 99: 4618-4625.

George D, et al. Detection and Quantification of Chimerism by Droplet Digital PCR. *Chimerism*. 2013 Jul-Sep;4(3):102-8.

Jimenez-Velasco, et al. Reliable Quantification of Hematopoietic Chimerism after Allogeneic Transplantation for Acute Leukemia using Amplification by qPCR of Null Alleles and Insertion/Deletion Polymorphisms. *Leukemia*, 2005; 19: 336-423.

Koldehoff, et al. Quantitative Analysis of Chimerism after Allogeneic Stem Cell Transplantation by Real-Time Polymerase Chain Reaction with Single Nucleotide Polymorphisms, Standard Tandem Repeats, and Y-Chromosome-Specific Sequences. *Am J Hematol*, 2006; 81: 735-746.

Thiede, et al. Strategies and Clinical Implications of Chimerism Diagnostics after Allogeneic Hematopoietic Stem Cell Transplantation. *Acta Haematol*, 2004; 112: 16-23.

Verhoeven, et al. A Novel High-Throughput Droplet Digital PCR-Based Indel Quantification Method for the Detection of Circulating Donor-derived Cell-free DNA After Kidney Transplantation. *Transplantation*. 2022 Sep 1;106(9):1777-1786.

Willasch, et al. Sequence Polymorphism Systems for Quantitative Real-Time Polymerase Chain Reaction to Characterize Hematopoietic Chimerism-High Informativity and Sensitivity As Well As Excellent Reproducibility and Precision of Measurement. *Laboratory Hematology*, 2007; 13: 73-84.

Glossary

Bi-allelic	an allele which exists in two variant forms - a major and minor allele. Individuals may be homozygous for either variant or heterozygous
Calibrator	a sample used as the basis for comparative results
Chimerism	A biological state where two or more genetically distinct cell populations exist within a single individual
Comparative C_T ($\Delta\Delta C_T$) method	Method for determining relative target quantity in samples. The relative quantity of target in each sample is determined by comparing the C_q in each sample to the C_q in the reference sample.
C_q	Quantification Cycle; The fractional PCR cycle used for quantification; also referred to as the threshold cycle (C_T), or crossing point (C_p) and take-off point (TOP)
Experiment	a collection of genotyping and monitoring reactions that are carried out simultaneously (i.e, on a single plate or in a single set tubes for RGQ cyclers).
INDEL	Abbreviation for insertion/deletion polymorphism; a class of DNA mutation characterized by the loss or gain of genetic material at a specific locus
Informative assay	An assay capable of distinguishing between genetic material from two or more sources; An informative assay is an assay for a marker allele that is present (positive) in one individual genome and absent (negative) in the other genome(s).
Monitoring assay	an assay designed for use in quantification of a specific marker, allele, or analyte. The assay must demonstrate high specificity for accurate quantification and high sensitivity to achieve a desirable limit of detection
Reference assay	an assay designed specifically to detect the total amount of DNA in a reaction; used in quantification to normalize the amount of input DNA between reactions; also referred to as an endogenous control assay
Reference Gene	Gene used for normalization and relative quantification

Reference Sample	a control sample possessing the target DNA (typically 100% target DNA) to be quantified in the experimental sample by relative quantification; also referred to as a calibrator
Recipient Sample	a control sample possessing the target DNA (typically 100% target DNA) to be quantified in the experimental sample by relative quantification; also referred to as a calibrator
Relative Quantification	a method of quantification where quantity of an unknown sample is derived by comparison to a reference sample
Sample	A sample is a unique donor or recipient
Passive reference	A dye that produces fluorescence signal. Because the passive reference signal should be consistent across all wells, it is used to normalize the reporter dye signal to account for non-PCR related fluorescence fluctuations caused by minor well-to-well differences in concentrations or volume. Normalization to the passive reference signal allows for high data precision.
Normalized reporter (Rn)	Fluorescence signal from the reporter dye normalized to the fluorescence signal of the passive reference.
Baseline-corrected normalized reporter (ΔRn)	The magnitude of normalized fluorescence signal generated by the reporter. In experiments that contain data from qPCR, the magnitude of normalized fluorescence signal generated by the reporter at each cycle the PCR amplification. ΔRn (cycle) = Rn (cycle) – Rn (baseline), where Rn = normalized reporter
Threshold	The intensity of fluorescence that must be exceeded for each reaction to be seen as positive
Partition	A compartment where the PCR reaction takes place
Nanoplate	A compartment where the PCR reaction takes place
CNV	Copy Number Variation (CNV) refers to a type of genetic variation in which the number of copies of a particular segment of DNA differs between individuals
UPC	Universal Plasmid Control (UPC). A synthetic control sample possessing the targets for Assays in the DigitalTRACE™ panel (typically 100% target DNA)

End User Software License Agreement

JETA MOLECULAR BV END-USER SOFTWARE LICENSE AGREEMENT AND LIMITED PRODUCT WARRANTY

NOTICE TO USER: PLEASE READ THIS DOCUMENT CAREFULLY. THIS IS THE CONTRACT BETWEEN YOU AND JETA MOLECULAR BV (JETA) REGARDING THIS SOFTWARE PRODUCT. THIS AGREEMENT CONTAINS WARRANTY AND LIABILITY DISCLAIMERS AND LIMITATIONS. YOUR INSTALLATION AND USE OF THIS JETA SOFTWARE IS SUBJECT TO THE TERMS AND CONDITIONS CONTAINED IN THIS END-USER SOFTWARE LICENSE AGREEMENT, AND YOU WILL BE BOUND BY THESE TERMS AND CONDITIONS IF YOU INSTALL AND USE THE SOFTWARE. IF YOU DO NOT AGREE TO THESE TERMS AND CONDITIONS, YOU SHOULD SELECT THE DECLINE BUTTON BELOW, WHEREUPON INSTALLATION WILL STOP, AND PROMPTLY RETURN THIS SOFTWARE, TOGETHER WITH ALL PACKAGING, TO JETA AND YOUR PURCHASE PRICE WILL BE REFUNDED

This JETA End-User License Agreement accompanies a JETA software product ("Software") and related explanatory materials ("Documentation"). The term "Software" includes all executable files and applications included as part of the software product, and any upgrades, modified versions, updates, additions and copies of the Software licensed to you by JETA. The term "JETA" as used in this License, means JETA Corporation. The term "License" or "Agreement" means this End-User Software License Agreement. The term "you", "Customer" or "Licensee" means the purchaser of this license to the Software. A copy of this End-User Software License Agreement can be printed by copying the text and pasting it into word processing software that provides printing capabilities.

Third-party Products

This Software includes the following third-party software products (hereinafter collectively referred to as the "Third-Party Products"):

Microsoft - .NET Framework version 4.5

Moq

Copyright (c) 2007. Clarius Consulting, Manas Technology Solutions, InSTEDD

<http://www.moqthis.com/>

All rights reserved.

NUnit

Portions Copyright © 2002-2013 Charlie Poole or Copyright © 2002-2004 James W. Newkirk, Michael C. Two, Alexei A. Vorontsov or Copyright © 2000-2002 Philip A. Craig

PDFsharp

Copyright (c) 2005-2007 empira Software GmbH, Cologne (Germany)

SharpZipLib

WPF Toolkit

Title

Title, ownership rights and intellectual property rights in and to the Software and Documentation shall at all times remain with JETA Corporation and its subsidiaries. Title, ownership rights and intellectual property rights in and to the Third-Party Product shall at all times remain with their respective licensors. All rights not specifically granted by this License, including Federal and international copyrights, are reserved by JETA and its subsidiaries.

Copyright

The Software, including its structure, organization, code, user interface and associated Documentation, is a proprietary product of JETA Corporation or its suppliers, and is protected by international laws of copyright, as well as other intellectual property laws and treatises. The law provides for civil and criminal penalties for anyone in violation of the laws of copyright.

License

Use of the Software

1 Subject to the terms and conditions of this Agreement, JETA grants you a non-exclusive, nontransferable license only to install and use the Software on a single computer for the sole purpose of using the Software that you

have purchased. You must use the Software in accordance with all applicable laws and regulations. You may transfer the Software to another single computer (or network, if a network version), but the Software may never be installed on more than one computer (or more than one network, if a network version) at any one time.

2 If the Software uses registration codes, access to the number of licensed copies of Software is controlled by a registration code. For example, if you have a registration code that enables you to use three copies of Software simultaneously, you cannot install the Software on more than three separate computers.

3 YOU MAY MAKE ONLY ONE COPY OF THE SOFTWARE IN MACHINE-READABLE FORM SOLELY FOR BACKUP OR ARCHIVAL PURPOSES, PROVIDED YOU REPRODUCE ON ANY SUCH COPY ALL COPYRIGHT NOTICES AND ANY OTHER PROPRIETARY LEGENDS FOUND ON THE ORIGINAL. JETA SHALL NOT BE RESPONSIBLE FOR THE QUALITY, INTEGRITY, FUNCTIONALITY, OR PERFORMANCE OF ANY COPY OF THE SOFTWARE. YOU MAY NOT MAKE ANY OTHER COPIES OF THE SOFTWARE.

Restrictions

1 You may not copy, transfer, rent, modify, use, merge, or translate the Software, or the associated Documentation, in whole or in part, except as expressly permitted in this Agreement or in the license agreements covering the Third-Party Product.

2 You may not reverse assemble, decompile, or otherwise reverse engineer the Software.

3 You may not remove any proprietary, copyright, trade secret or warning legend from the Software or any Documentation.

4 You agree to comply fully with all export laws and restrictions and regulations of the United States and applicable foreign agencies or authorities, and you agree that you have the sole responsibility to obtain such licenses to export, re-export or import as may be required after delivery of the Software to you. You agree that you will not export or re-export, directly or indirectly, the Software into any country prohibited by the United States Export Administration Act and the regulations thereunder or other applicable United States law.

5 Subject to the respective license terms and conditions applicable to the Third-Party Product included in this Software, you may not modify, sell, rent, transfer (except temporarily in the event of a computer malfunction), resell for profit, or distribute this license or the Software, or create derivative works based upon the Software, or any part thereof or any interest therein.

Trial

If this license is granted on a trial basis, you are hereby notified that license management software may be included to automatically cause the Software to cease functioning at the end of the trial period.

Termination

You may terminate this Agreement by discontinuing use of the Software, removing all copies from your computers and storage media, and returning the Software and Documentation, and all copies thereof, to JETA. JETA may terminate this Agreement if you fail to comply with all of its terms, in which case you agree to discontinue using the Software, remove all copies from your computers and storage media, and return the Software and Documentation, and all copies thereof, to JETA.

U.S. Government End Users

The Software is a "commercial item," as that term is defined in 48 C.F.R. 2.101 (Oct. 1995), consisting of "commercial computer software" and "commercial computer software documentation," as such terms are used in 48 C.F.R. 12.212 (Sept. 1995). Consistent with 48 C.F.R. 12.212 and 48 C.F.R. 227.7202-1 through 227.7202-4 (June 1995), all U.S. Government end users acquire the Software with only those rights set forth herein.

European Community End Users

If this Software is used within a country of the European Community, nothing in this Agreement shall be construed as restricting any right available under the European Community Software Directive, O.J. Eur. Comm. (No. L. 122) 42 (1991).

LIMITED WARRANTY AND LIMITATION OF REMEDIES

Limited Warranty. JETA warrants that for a period of ninety days from the beginning of the applicable warranty period (as described below) the Software will function substantially in accordance with the functions and features described in the Documentation delivered with the Software when properly installed, and that for a period of ninety days from the beginning of the applicable warranty period (as described below) the tapes, CDs, diskettes or other media bearing the Software will be free of defects in materials and workmanship under normal use.

The above warranties do not apply to defects resulting from misuse, neglect, or accident, including without limitation: operation outside of the environmental or use specifications, or not in conformance with the instructions for any instrument system, Software, or accessories; improper or inadequate maintenance by the user; installation of software or interfacing, or use in combination with Software or products not supplied or authorized by JETA and modification or repair of the products not authorized by JETA.

Warranty Period Commencement Date. The applicable warranty period for the Software installed by JETA's personnel or authorized representative begins on the earlier of: the date of installation; or three (3) months from the date of shipment for the Software. For Software installed by you or anyone other than JETA, the warranty period begins on the date the Software is delivered to you. The applicable warranty period for media begins on the date the media is delivered to you.

JETA MAKES NO OTHER WARRANTIES OF ANY KIND WHATSOEVER, EXPRESS OR IMPLIED, WITH RESPECT TO THE SOFTWARE OR DOCUMENTATION, INCLUDING, BUT NOT LIMITED TO, WARRANTIES OF FITNESS FOR A PARTICULAR PURPOSE OR MERCHANTABILITY OR THAT THE SOFTWARE OR DOCUMENTATION IS NON-INFRINGEMENT. ALL OTHER WARRANTIES ARE EXPRESSLY DISCLAIMED. WITHOUT LIMITING THE GENERALITY OF THE FOREGOING, JETA MAKES NO WARRANTIES THAT THE SOFTWARE WILL MEET YOUR REQUIREMENTS, THAT OPERATION OF THE LICENSED SOFTWARE WILL BE UNINTERRUPTED OR ERROR FREE OR WILL CONFORM EXACTLY TO THE DOCUMENTATION, OR THAT JETA WILL CORRECT ALL PROGRAM ERRORS. JETA'S SOLE LIABILITY AND RESPONSIBILITY FOR BREACH OF WARRANTY RELATING TO THE SOFTWARE OR DOCUMENTATION SHALL BE LIMITED, AT JETA'S SOLE OPTION, TO (1) CORRECTION OF ANY ERROR IDENTIFIED TO JETA IN A WRITING FROM YOU IN A SUBSEQUENT RELEASE OF THE SOFTWARE, WHICH SHALL BE SUPPLIED TO YOU FREE OF CHARGE, (2) ACCEPTING A RETURN OF THE PRODUCT, AND REFUNDING THE PURCHASE PRICE UPON RETURN OF THE PRODUCT AND REMOVAL OF ALL COPIES OF THE SOFTWARE FROM YOUR COMPUTERS AND STORAGE DEVICES, (3) REPLACEMENT OF THE DEFECTIVE SOFTWARE WITH A FUNCTIONALLY EQUIVALENT PROGRAM AT NO CHARGE TO YOU, OR (4) PROVIDING A REASONABLE WORK AROUND WITHIN A REASONABLE TIME. JETA'S SOLE LIABILITY AND RESPONSIBILITY UNDER THIS AGREEMENT FOR BREACH OF WARRANTY RELATING TO MEDIA IS THE REPLACEMENT OF DEFECTIVE MEDIA RETURNED WITHIN 90 DAYS OF THE DELIVERY DATE. THESE ARE YOUR SOLE AND EXCLUSIVE REMEDIES FOR ANY BREACH OF WARRANTY. WARRANTY CLAIMS MUST BE MADE WITHIN THE APPLICABLE WARRANTY PERIOD.

LIMITATION OF LIABILITY

IN NO EVENT SHALL JETA OR ITS SUPPLIERS BE RESPONSIBLE OR LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY OR UNDER ANY STATUTE (INCLUDING WITHOUT LIMITATION ANY TRADE PRACTICE, UNFAIR COMPETITION OR OTHER STATUTE OF SIMILAR IMPORT) OR ON ANY OTHER BASIS FOR SPECIAL, INDIRECT, INCIDENTAL, MULTIPLE, PUNITIVE, OR CONSEQUENTIAL DAMAGES ARISING OUT OF THE POSSESSION OR USE OF, OR THE INABILITY TO USE, THE SOFTWARE OR DOCUMENTATION, EVEN IF JETA IS ADVISED IN ADVANCE OF THE POSSIBILITY OF SUCH DAMAGES, INCLUDING WITHOUT LIMITATION DAMAGES ARISING FROM OR RELATED TO LOSS OF USE, LOSS OF DATA, DOWNTIME, OR FOR LOSS OF REVENUE, PROFITS, GOODWILL OR BUSINESS OR OTHER FINANCIAL LOSS. IN ANY CASE, THE ENTIRE LIABILITY OF JETA AND ITS SUPPLIERS UNDER THIS LICENSE, OR ARISING OUT OF THE USE OF THE SOFTWARE, SHALL NOT EXCEED IN THE AGGREGATE THE PURCHASE PRICE OF THE SOFTWARE, IF PURCHASED INDEPENDENTLY, OR THE FAIR MARKET VALUE OF THE SOFTWARE, IF PURCHASED AS PART OF A JETA IN VITRO DIAGNOSTICS KIT.

SOME STATES, COUNTRIES OR JURISDICTIONS LIMIT THE SCOPE OF OR PRECLUDE LIMITATIONS OR EXCLUSION OF REMEDIES OR DAMAGES, OR OF LIABILITY, SUCH AS LIABILITY FOR GROSS NEGLIGENCE OR WILLFUL MISCONDUCT, AS OR TO THE EXTENT SET FORTH ABOVE, OR DO NOT ALLOW IMPLIED WARRANTIES TO BE EXCLUDED. IN SUCH STATES, COUNTRIES OR JURISDICTIONS, THE LIMITATION OR EXCLUSION OF WARRANTIES, REMEDIES,

DAMAGES OR LIABILITY SET FORTH ABOVE MAY NOT APPLY TO YOU. HOWEVER, ALTHOUGH THEY SHALL NOT APPLY TO THE EXTENT PROHIBITED BY LAW, THEY SHALL APPLY TO THE FULLEST EXTENT PERMITTED BY LAW. YOU MAY ALSO HAVE OTHER RIGHTS THAT VARY BY STATE, COUNTRY OR OTHER JURISDICTION.

CONFIDENTIALITY

Except as expressly allowed by this License, you will not use or disclose any Software, Documentation or related technology, idea, algorithm or information except to the extent you can document that it is generally available for use and disclosure by the public without charge or license. You recognize and agree that there is no adequate remedy at law for a breach of this provision, that such a breach would irreparably harm JETA and that JETA is entitled to equitable relief (without need to post a bond) with respect to any such breach or potential breach in addition to any other remedies.

INDEMNITY

You shall defend, indemnify and hold harmless JETA and its officers, directors, employees and agents (the "JETA Indemnitees"), from and against all suits, claims, demands, losses, liabilities, damages and expenses (including reasonable attorneys' fees and costs) that the JETA Indemnitees may suffer or incur in connection with: (i) any third-party claim arising from your breach of this License; and (ii) any infringement by you of JETA's or its licensors' intellectual property rights in the Software or the Third-Party Product.

BUSINESS OBJECTS SUBLICENSE GRANT, TERMS AND CONDITIONS

You agree not to alter, disassemble, decompile, translate, adapt or reverse-engineer the Business Objects Software or the report file (.RPT) format;

You agree not to distribute the Business Objects Software with any general-purpose report writing, data analysis or report delivery product or any other product that performs the same or similar functions as Business Objects' product offerings;

You agree not to use the Business Objects Software to create for distribution a product that is generally competitive with Business Objects' product offerings;

You agree not to use the Business Objects Software to create for distribution a product that converts the report file (.RPT) format to an alternative report file format used by any general-purpose report writing, data analysis or report delivery product that is not the property of Business Objects;

You agree not to use the Business Objects Software on a rental or timesharing basis or to operate a service bureau facility for the benefit of third-parties;

BUSINESS OBJECTS AND ITS SUPPLIERS DISCLAIM ALL WARRANTIES, EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, AND NONINFRINGEMENT OF THIRD PARTY RIGHTS. BUSINESS OBJECTS AND ITS SUPPLIERS SHALL HAVE NO LIABILITY WHATSOEVER UNDER THIS AGREEMENT OR IN CONNECTION WITH THE BUSINESS OBJECTS SOFTWARE.

GENERAL

This Agreement shall be governed by laws of the Netherlands, exclusive of its conflict of laws provisions. This Agreement shall not be governed by the United Nations Convention on Contracts for the International Sale of Goods. You hereby agree that the courts located in the Utrecht Netherlands, will constitute the sole and exclusive forum for the resolution of any and all disputes arising out of or in connection with this Agreement and you hereby irrevocably consent to the personal jurisdiction and venue of such courts and irrevocably waive any objections thereto. This Agreement contains the complete agreement between the parties with respect to the subject matter hereof, and supersedes all prior or contemporaneous agreements or understandings, whether oral or written. If any provision of this Agreement is held by a court of competent jurisdiction to be contrary to law, that provision will be enforced to the maximum extent permissible and the remaining provisions of this Agreement will remain in full force and effect. The controlling language of this Agreement, and any proceedings relating to this Agreement, shall be English. You agree to bear any and all costs of translation, if necessary. The headings to the sections of this Agreement are used for convenience only and shall have no substantive meaning. All questions concerning this Agreement shall be directed to: JETA Molecular BV, Krommewetering 101C, 3543AN,

Utrecht, Netherlands, Attention: Legal Department.

Unpublished rights reserved under the copyright laws of the United States.

JETA Molecular BV, Krommewetering 101C, 3543AN, Utrecht, Netherlands.

JETA and its logo design are registered trademarks of JETA Molecular BV in the U.S.

and/or certain other countries.

QTRACE® is a registered trademark of ElsworthMolecular Holding BV in the U.S. and/or certain other countries.

All other trademarks, copyrights, patents, service marks, logos and trade names are the sole property of their respective owners.