DigitalTRACE™ INDEL and HLA Assays

Intended use:

This kit is intended for the *in vitro*, non-automated, quantitative measurement of the levels of recipient-derived genomic DNA from white blood cells in peripheral blood and/or bone marrow.

Examination Method Principles:

The DigitalTRACE™ System consists of 95 individual quantification assays and the TRACE Analysis™ Software. The DigitalTRACE™ INDEL Assays are a set of 95 genetic markers that are able to differentiate, and then quantify, the contributors to a human-mixed DNA sample. Each of the 95 assays is designed to a distinct bi-allelic insertion/deletion (INDEL) or copy number polymorphisms or HLA alleles in the human genome. TRACE Analysis™ Software guides the user through reaction set-up for both genotyping and quantification and analyzes the collected data. The procedure for determining the level of a genome of interest in a sample consists of two parts: a genotyping test and a quantification test.

Materials provided:

DigitalTRACE™ INDEL Assays: 20X buffered solution sufficient for 24 dPCR reactions containing primers and probe for detecting the variant of interest and a reference gene.

Limitations:

The DigitalTRACE™ Analysis System kits may be used solely in accordance with the DigitalTRACE™ Kits manual and for use with components described in the DigitalTRACE™ Kits manual. JETA Molecular grants no license under any of its intellectual property to use or incorporate the enclosed components of this kit with any components not included within this kit except as described in the DigitalTRACE™ Kits manual and additional protocols available at www.jetamolecular.com.

Storage and Handling Requirements:

Examine the shipment upon receipt and if the integrity of the products has been compromised during shipment, immediately contact your local customer support representative. Upon receipt, the individual components should be stored according to the temperatures listed on the labels. **Note:** When storage recommendations are observed, both unopened and opened/recapped tubes are stable until the expiration date indicated on the label.

Do not use any component that visibly shows signs of having been compromised (e.g., particulate matter, presence of foreign debris, cloudy appearance, discoloration). **Do not use any component after its expiry date.**

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DigitalTRACE™ INDEL and HLA Assays

Intended use:

This kit is intended for the *in vitro*, non-automated, quantitative measurement of the levels of recipient-derived genomic DNA from white blood cells in peripheral blood and/or bone marrow.

Examination Method Principles:

The DigitalTRACE™ System consists of 95 individual quantification assays and the TRACE Analysis™ Software. The DigitalTRACE™ INDEL Assays are a set of 95 genetic markers that are able to differentiate, and then quantify, the contributors to a human-mixed DNA sample. Each of the 95 assays is designed to a distinct bi-allelic insertion/deletion (INDEL) or copy number polymorphisms in the human genome. TRACE Analysis™ Software guides the user through reaction set-up for both genotyping and quantification and analyzes the collected data. The procedure for determining the level of a genome of interest in a sample consists of two parts: a genotyping test and a quantification test.

Materials provided:

DigitalTRACE™ INDEL Assays: 20X buffered solution sufficient for 24 dPCR reactions containing primers and probe for detecting the variant of interest and a reference gene.

Limitations:

The DigitalTRACE™ Analysis System kits may be used solely in accordance with the DigitalTRACE™ Kits manual and for use with components described in the DigitalTRACE™ Kits manual. JETA Molecular grants no license under any of its intellectual property to use or incorporate the enclosed components of this kit with any components not included within this kit except as described in the DigitalTRACE™ Kits manual and additional protocols available at www.jetamolecular.com.

Storage and Handling Requirements:

Examine the shipment upon receipt and if the integrity of the products has been compromised during shipment, immediately contact your local customer support representative. Upon receipt, the individual components should be stored according to the temperatures listed on the labels. **Note:** When storage recommendations are observed, both unopened and opened/recapped tubes are stable until the expiration date indicated on the label.

Do not use any component that visibly shows signs of having been compromised (e.g., Particulate matter, presence of foreign debris, cloudy appearance, discoloration). **Do not use any component after its expiry date.**

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Analytical Performance characteristics

The DigitalTRACE™ Analysis System has been validated to reliably quantify the presence of 0.1% minor component DNA in the presence of 150ngs total DNA.

Additional Information about the QTRACE Analysis System is detailed in the Operator's Manual in the locations cited

	Location in Operator's Manual
Warnings	Warnings Section
Product Use Limitations	Product Use Limitations
Primary sample, collection and storage	DNA Sample Requirements
Materials Not Provided	Materials
Reagent Preparation	Genotyping Test and Monitoring Test Sections
Procedure	Genotyping Test and Monitoring Test Sections
Calculation of Examination Results	Data Analysis Algorithms Section, Genotyping and Monitoring Data Analysis and Report Sections
Interpretation of Results	Genotyping and Monitoring Data Analysis and Report Sections
Safe Waste Disposal	Waste Disposal
Literature References	Bibliography

NOTICE TO PURCHASER: LIMITED LICENSES: Not for resale.

DigitalTRACE™ and TRACE Analysis™ are registered trademarks of ElsworthMolecular Holding BV or its subsidiaries in the US and certain other countries. DigitalTRACE™ Products and DigitalTRACE™ Software are licensed exclusively from ElsworthMolecular Holding BV.

For Research Use Only. Not for use in diagnostic applications.

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Analytical Performance characteristics

The QTRACE Analysis System has been validated to reliably quantify the presence of 0.1% minor component DNA in the presence of 150ngs total DNA.

Additional Information about the QTRACE Analysis System is detailed in the Operator's Manual in the locations cited

	Location in Operator's Manual
Warnings	Warnings Section
Product Use Limitations	Product Use Limitations
Primary sample, collection and storage	DNA Sample Requirements
Materials Not Provided	Materials
Reagent Preparation	Genotyping Test and Monitoring Test Sections
Procedure	Genotyping Test and Monitoring Test Sections
Calculation of Examination Results	Data Analysis Algorithms Section, Genotyping and Monitoring Data Analysis and Report Sections
Interpretation of Results	Genotyping and Monitoring Data Analysis and Report Sections
Safe Waste Disposal	Waste Disposal
Literature References	Bibliography

NOTICE TO PURCHASER: LIMITED LICENSES: Not for resale.

DigitalTRACE™ and TRACE Analysis™ are registered trademarks of ElsworthMolecular Holding BV or its subsidiaries in the US and certain other countries. DigitalTRACE™ Products and DigitalTRACE™ Software are licensed exclusively from ElsworthMolecular Holding BV.

For Research Use Only. Not for use in diagnostic applications.

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